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Did you:

- select the appropriate format (e.g. 24W etc.)?
- include the necessary controls not already incorporated in the design?
- include the appropriate probes that target the appropriate species?
- check the quality and form (e.g. length, linkers, etc.) of any uploaded probes on this microarray?
- for your intended use, include only those uploaded probes for which you agree to take full responsibility?
- include the appropriate probe groups with the correct probes?
- fill all the unused features on the microarray you wanted filled?
- include the required technical replicates and/or experimental replicates?
- create a design that works with your equipment, kits, and protocols?

After you select the "Submit" button on the main, this microarray will be submitted to Agilent microarray design. To enable ordering an Agilent representative will contact you to confirm.

By selecting "Done" button below, you are confirming that you have reviewed your microarray to the checklist above. Independent of the checklist, you are responsible for your design's fit.

By selecting "Submit", you are also confirming acceptance for an Agilent representative to call. You will NOT be sent to you and invoiced until an Agilent representative has contacted you and you received.

If ordered, custom microarrays are warranted as articles of manufacture (workmanship) do not agree to the above checklist or do not wish to order and have an Agilent representative. "Cancel" button.

Done Cancel

Array Design is successfully submitted.

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Application Type: **microRNA** [Switch Application Type](#)

Search

Microarray Probe Group Probe

Design Wizards [View All](#)

Create a Microarray Design from Existing ProbeGroup(s) [Info](#)

Application Type: **CGH** [Switch Application Type](#)

Species [Info](#) [Search](#) [Upload](#) [Genomic Tiling](#) [Reannotate](#) [DNA Analytics Download](#)

Design ID:

Search Your Probes Interval Finder HD Probe Search

Standard HD Probe Search [Info](#) Advanced HD Probe Search [Info](#)

Job Information

Search Name [Info](#) *

Select HD Search By Genomic Intervals

Interval Options

Application Type: **Expression** [Switch Application Type](#)

Species

Build Number [Search](#) [Upload](#) [Simple Tiling](#) [GE Probe Design](#) [GE Probe Check](#) [Reannotate](#)

Probe Options

Filters [Info](#)

Prefer Catalog Probe

Use TM Filter [Info](#)

Similarity Filter [Info](#)

Base Composition Methodology [Info](#) Tm Matching Methodology [Info](#)

Probe Details

Design Job Name

Probe Length 60 bp

Probes per Target 1

Masking Vector [Info](#) Repeat [Info](#)

Probe Orientation Sense [Info](#) Antisense [Info](#)

Design Options Best Probe Methodology [Info](#) Best Distribution Methodology [Info](#)

Design with 3' Bias [Info](#)

Target File Details

Species H. sapiens

Upload in FASTA Format [Info](#)

Upload as Genbank Accessions [Info](#)

Transcriptome Details

Use Target File as Transcriptome [Info](#)

Select Agilent-provided Transcriptome (by Species) [Info](#)

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

<input checked="" type="radio"/> Base Composition Methodology Info		<input type="radio"/> Tm Matching Methodology Info	
Probe Details		Target File Details	
Design Job Name	<input type="text"/>	Species	<input type="text" value="H. sapiens"/>
Probe Length	<input type="text" value="60"/> bp	<input checked="" type="radio"/> Upload in FASTA Format Info	<input type="text"/>
Probes per Target Info	<input type="text" value="1"/>	<input type="radio"/> Upload as Genbank Accessions Info	<input type="text"/>
Masking	<input checked="" type="checkbox"/> Vector Info <input checked="" type="checkbox"/> Repeat Info	Transcriptome Details	
Probe Orientation	<input checked="" type="radio"/> Sense Info <input type="radio"/> Antisense Info	<input type="radio"/> Use Target File as Transcriptome Info	
Design Options	<input checked="" type="radio"/> Best Probe Methodology Info <input type="radio"/> Best Distribution Methodology Info	<input checked="" type="radio"/> Select Agilent-provided Transcriptome (by Species) Info	<input type="text" value="H. sapiens"/>
<input type="checkbox"/> Design with 3' Bias Info		<input type="radio"/> Upload Transcriptome File Info	<input type="text"/>
<input type="button" value="Submit"/> <input type="button" value="Cancel"/>			

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