

Big Blue Cell Lines

INSTRUCTION MANUAL

Catalog #726000 (Big Blue Rat Cell Line)

Revision A

For In Vitro Use Only

726000-12

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

BIG BLUE CELL LINES

CONTENTS

Materials Provided	1
Storage Conditions	1
Additional Materials Required	1
Introduction	2
Cell Culture Protocols	3
Starting the Cells in Culture	3
Passaging the Cells	4
Freezing the Cells for Long-Term Storage.....	4
Treating the Cells with a Test Compound	5
Troubleshooting	6
Preparation of Media and Reagents	6
References	7
Endnotes	7
MSDS Information	7
Quick-Reference Protocol	9

Big Blue Cell Lines

Catalog #726000

MATERIALS PROVIDED

Materials provided ^a	Quantity
Big Blue rat cell line	1-ml cryovial (~1 × 10 ⁶ cells)

^a Immediately on receipt, store the cryovial of cells in liquid nitrogen. Protective gloves, clothing, and a face mask or goggles should be worn when transferring the cryovial to and from liquid nitrogen.

STORAGE CONDITIONS

Store immediately in liquid nitrogen.

ADDITIONAL MATERIALS REQUIRED

G418 Sulfate [Stratagene Catalog #200399 (1 g)]

Liquid nitrogen storage tank

Growth medium[§]

Freezing medium[§]

Phosphate-buffered saline (PBS)[§]

Plasticware

Incubator (5–7% CO₂, 37°C)

Controlled freezing container [e.g., the StrataCooler Cryo preservation module (Stratagene Catalog #400005), or the equivalent]

Tabletop centrifuge

Rubber policeman

Trypsin–EDTA (0.05% trypsin and 0.53 mM tetrasodium ethylenediaminetetraacetic acid)

Big Blue Transgenic Rodent Mutagenesis Assay System Instruction Manual

RecoverEase DNA isolation kit [Stratagene Catalog #720202 (30 preparations), Catalog #720203 (10 preparations)]

[§]See Preparation of Media and Reagents.

Revision A

© Agilent Technologies, Inc. 2009.

INTRODUCTION

The Big Blue rat cell line is useful for *in vitro* mutagenesis studies, and is a part of the Big Blue transgenic rodent mutagenesis assay system, which is used to assess the genetic toxicity of test compounds.¹⁻⁴

The Big Blue rat cell line is derived from a Rat2 embryonic fibroblast cell line transfected with the Big Blue λ LIZ shuttle vector and the pSV2NEO plasmid.^{5,6} The Big Blue λ LIZ shuttle vector (Figure 1) is 45 kb in size and contains the *lacI* and *lacZ* genes. Each cell contains 50–70 copies of the λ LIZ shuttle vector integrated into the genome at two sites. The pSV2NEO plasmid provides resistance to antibiotic G418. The Big Blue rat cell line is a thymidine kinase mutant, and the chromosome number is polyploid.

The *lacI* gene in the Big Blue λ LIZ shuttle vector is the target for mutagenesis. After the cells are treated with a test compound, *in vitro* packaging rapidly and efficiently recovers the λ LIZ shuttle vector from the genomic DNA of the cell. Bacterial detection systems screen for mutations in the target gene. The target gene can be recovered from mutant plaques for sequence analysis. In addition to short-term *in vitro* mutagenesis studies, these cell lines are useful for various toxicological applications, studies of the effects of oncogenes and tumor suppressor genes on mutation frequencies, and studies of DNA repair.




Figure 1 The Big Blue λ LIZ shuttle vector.

CELL CULTURE PROTOCOLS

Starting the Cells in Culture

1. Prepare the growth medium.[§] Place 10 ml of growth medium in a 15-ml conical tube.
2. Thaw the frozen cryovial of cells within 40–60 seconds by rapid agitation in a 37°C water bath. Remove the cryovial from the water and immediately immerse the cryovial in 70% (v/v) ethanol at room temperature.

Note *All operations from this point should be carried out using sterile technique.*

3. Transfer the cell suspension to the conical tube containing the growth medium.
4. Centrifuge the cells at 200 × g for 5 minutes at room temperature. Remove the growth medium from the conical tube by aspiration.
5. Place 10 ml of growth medium in a 75-cm² tissue culture flask. Resuspend the cells in the conical tube in 5 ml of growth medium. Transfer the cell suspension to the tissue culture flask.

Note *Add 200 µg/ml of G418 Sulfate to the growth medium for the Big Blue rat cell line.*

6. Place the cells in a 37°C incubator with 5–7% CO₂.

[§] See *Preparation of Media and Reagents*.

Passaging the Cells

Note *Split the cells when the monolayer of cells is ~80–90% confluent or approximately every 3–4 days.*

1. Wash the cells with PBS. Remove the PBS by aspiration.
2. Add 1.5 ml of trypsin–EDTA to each 75-cm² tissue culture flask. Incubate the cells at room temperature for ~2–5 minutes. Tap the flask to release the adherent cells.

Note *Incubate the cells with trypsin–EDTA for the minimum time necessary for the adherent cells to release from the flask. Overexposure to trypsin can kill the cells.*

3. Dilute the cells with growth medium to a final volume of 5–10 ml. The serum in the medium inactivates the trypsin. Transfer the cells to the desired number of 75-cm² or 175-cm² tissue culture flasks. Depending on the desired use for the cells, use an inoculum of 1×10^3 to 1×10^6 cells. Add growth medium to the tissue culture flasks to a final volume of 15 ml for each 75-cm² flask or 30 ml for each 175-cm² flask.
4. Place the cells in a 37°C incubator with 5–7% CO₂.

Freezing the Cells for Long-Term Storage

Note *Maintain freezer stocks of early-passaged cells to ensure long-term stability of the cell line.*

1. Wash the cells with PBS. Remove the PBS by aspiration.
2. Add 1.5 ml of trypsin–EDTA to each 75-cm² tissue culture flask. Incubate the cells at room temperature for ~2–5 minutes. Tap the flask to release the adherent cells.

Note *Incubate the cells with trypsin–EDTA for the minimum time necessary for the adherent cells to release from the flask. Overexposure to trypsin can kill the cells.*

3. Dilute the cells with 8.5 ml of growth medium. The serum in the medium inactivates the trypsin.
4. Transfer the cell suspension to a 15-ml conical tube. Centrifuge the cells at $200 \times g$ for 5 minutes at room temperature. Remove the growth medium from the tube by aspiration.

5. Resuspend the cells in 1–2 ml of freezing medium.
6. Transfer the cells to cryovials. Place the cryovials in a chilled controlled freezing container. Incubate the cryovials at -80°C overnight.
7. The following day, transfer the cryovials to liquid nitrogen for long-term storage, taking the proper safety precautions for working with liquid nitrogen. If properly stored, the cells should remain stable for >1 year.

Treating the Cells with a Test Compound

Note *Theoretically, 150,000 cultured cells yields ~ 1 μg of genomic DNA. The Big Blue rat cell line achieves rescue efficiencies of $\sim 50,000$ – $100,000$ pfu/ μg of genomic DNA.*

1. Treat the cells with the test compound of interest when the monolayer of cells is ~ 30 – 40% confluent or is at the confluence stated in an established protocol for the test compound.^{5,6}
2. After treatment with the test compound and the desired expression time, isolate the DNA for mutagenesis studies. To isolate the DNA, remove the growth medium and wash the cells with PBS.
3. Harvest 1×10^7 to 2×10^7 cells in ~ 5 ml of PBS by scraping the surface of the tissue culture flasks with a rubber policeman.
4. Transfer the cell suspension to a 15-ml conical tube. Centrifuge the cells at $200 \times g$ for 5 minutes at room temperature. Remove the supernatant from the tube by aspiration.

Note *At this point, the pellet can be flash frozen in liquid nitrogen and stored at -80°C until DNA isolation is performed.*

5. Isolate the DNA from the cells. Proceed with step 3 of *Tissue Culture* in the *RecoverEase DNA Isolation Kit Instruction Manual*.

TROUBLESHOOTING

Observation	Suggestion(s)
The cells do not survive following resuspension in growth medium after storage in liquid nitrogen	Ensure that the cells are stored properly. Store the cryovials of cells in liquid nitrogen immediately on receipt.
	Ensure that the cells are thawed properly. Thaw the cells quickly and immediately dilute the cells in the growth medium.
The cells do not survive passaging	Ensure that the cells are passaged using the recommended conditions. When passaging, inoculate a fresh tissue culture flask with 1×10^3 to 1×10^6 cells.
	Split the cells when the monolayer of cells is at most 90% confluent to avoid starving the cells for growth medium.
	The cells may have been exposed to trypsin–EDTA for too long. Incubate the cells with trypsin–EDTA for the minimum time necessary for releasing the cells from the flask.
Low rescue efficiencies	Errors in preparing or packaging genomic DNA may cause low rescue efficiencies. Follow the instructions carefully in the <i>Big Blue Transgenic Rodent Mutagenesis Assay System Instruction Manual</i> .

PREPARATION OF MEDIA AND REAGENTS

<p>Freezing Medium</p> <p>9 ml of fetal bovine serum 1 ml of dimethylsulfoxide (DMSO) Filter sterilize</p>	<p>PBS</p> <p>137 mM NaCl 2.6 mM KCl 10 mM Na₂HPO₄ 1.8 mM KH₂PO₄ Adjust the pH to 7.4 with HCl</p>
<p>Growth Medium</p> <p>500 ml of Dulbecco’s Modified Eagle Medium (DMEM) (high glucose, with L-glutamine, without sodium pyruvate) 5 ml of 200 mM L-glutamine 5 ml of penicillin (5000 U/ml)–streptomycin (5000 µg/ml) mixture 50 ml of fetal bovine serum, heat inactivated</p>	

REFERENCES

1. Dyaico, M. J., Provost, G. S., Kretz, P. L., Ransom, S. L., Moores, J. C. *et al.* (1994) *Mutat Res* 307(2):461-78.
2. Dyaico, M. J., Stuart, G. R., Tobal, G. M., de Boer, J. G., Glickman, B. W. *et al.* (1996) *Carcinogenesis* 17(11):2347-56.
3. Kohler, S. W., Provost, G. S., Fieck, A., Kretz, P. L., Bullock, W. O. *et al.* (1991) *Proc Natl Acad Sci U S A* 88(18):7958-62.
4. Rogers, B. J., Provost, G. S., Young, R. R., Putman, D. L. and Short, J. M. (1995) *Mutat Res* 327(1-2):57-66.
5. Manjanatha, M. G., Chen, J. B., Harris, A. J., Shaddock, J. G., Jr. and Casciano, D. A. (1996) *Environ Mol Mutagen* 27:44.
6. Wyborski, D. L., Malkhosyan, S., Moores, J., Perucho, M. and Short, J. M. (1995) *Mutat Res* 334(2):161-5.

ENDNOTES

Triton® is a registered trademark of Union Carbide Chemicals and Plastics Co., Inc.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.

STRATAGENE

An Agilent Technologies Division

Big Blue Cell Lines

Catalog #726000

QUICK-REFERENCE PROTOCOL

Starting the Cells in Culture

- ◆ Place 10 ml of growth medium in a 15-ml conical tube
- ◆ Thaw the frozen cells quickly in a 37°C water bath
- ◆ Transfer the cell suspension to the conical tube containing the growth medium
- ◆ Centrifuge the cells at $200 \times g$ for 5 minutes at room temperature and remove the growth medium by aspiration
- ◆ Place 10 ml of growth medium in a 75-cm² tissue culture flask, resuspend the cells in the conical tube in 5 ml of growth medium, and transfer the cell suspension to the tissue culture flask

Note Add 200 $\mu\text{g/ml}$ of G418 Sulfate to the growth medium for the Big Blue rat cell line

- ◆ Place the cells in a 37°C incubator with 5–7% CO₂

Passaging the Cells

- ◆ Wash the cells with PBS
- ◆ Add 1.5 ml of trypsin–EDTA to each 75-cm² tissue culture flask, incubate the cells at room temperature for ~2–5 minutes, and tap the flask to release the adherent cells
- ◆ Dilute the cells with the growth medium and transfer the cells to the desired number of tissue culture flasks
- ◆ Place the cells in a 37°C incubator with 5–7% CO₂

Freezing the Cells for Long-Term Storage

- ♦ Wash the cells with PBS
- ♦ Add 1.5 ml of trypsin–EDTA mixture to each 75-cm² tissue culture flask containing cells, incubate the cells at room temperature for ~2–5 minutes, and tap the flask to release the adherent cells
- ♦ Dilute the cells with 8.5 ml of growth medium
- ♦ Transfer the cell suspension to a 15-ml conical tube, centrifuge at $200 \times g$ for 5 minutes at room temperature, and remove the growth medium by aspiration
- ♦ Resuspend the cells in 1–2 ml of freezing medium
- ♦ Transfer the cells to cryovials, place the cryovials in a chilled controlled freezing container, and incubate the cryovials at -80°C overnight
- ♦ The following day, transfer the cryovials to liquid nitrogen for long-term storage

Treating the Cells with a Test Compound

- ♦ Treat the cells with the test compound when the monolayer is ~30–40% confluent or is at the confluence stated in the protocol for the test compound
- ♦ To isolate DNA for mutagenesis studies after treatment with the test compound, remove the growth medium and wash the cells with PBS
- ♦ Harvest 1×10^7 to 2×10^7 cells in ~5 ml of PBS by scraping the surface of the tissue culture flasks with a rubber policeman
- ♦ Transfer the cell suspension to a 15-ml conical tube, centrifuge at $200 \times g$ for 5 minutes at room temperature, and remove the supernatant by aspiration
- ♦ If DNA isolation will be performed immediately, resuspend the cells in 5 ml of ice-cold lysis buffer and proceed with step 3 of *Tissue Culture* in the *RecoverEase DNA Isolation Kit Instruction Manual*
- ♦ If DNA isolation will not be performed immediately, flash freeze the pellet by immersing the 15-ml conical tube containing the cells in liquid nitrogen and store the frozen cells at -80°C until DNA isolation is performed