

# Universal Mouse Reference RNA

Catalog #740100



**Storage** Store the Universal Mouse Reference RNA at  $-80^{\circ}\text{C}$ . Store the RNase-free water at  $-20^{\circ}\text{C}$ .

## INTRODUCTION

Stratagene's Universal Mouse Reference RNA (UMRR) is composed of total RNA from 11 mouse cell lines. The reference RNA is designed to be used as a reference for microarray gene-profiling experiments. Since RNA species differ in abundance between cell lines, an ideal reference sample should represent these different RNAs. Equal quantities of DNase-treated total RNA from each cell line were pooled to make the Universal Mouse Reference RNA. This Universal Reference RNA is suitable for mouse microarray experiments.

## MATERIALS PROVIDED

Material Provided	Quantity
Reference RNA	2 tubes x 200 $\mu\text{g}$ each
RNase-free water	1.5 ml

Cell Line Derivations	
embryo	T-lymphocyte (thymus)
embryo, fibroblast	mammary gland
kidney	muscle myoblast
liver, hepatocyte	skin
lung, alveolar macrophage	testis
B-lymphocyte	

## ADDITIONAL MATERIALS REQUIRED

RNase-free 70% Ethanol

## PROTOCOL

The UMRR is provided in a solution of 70% ethanol and 0.1 M sodium acetate. Prepare the UMRR for use as follows:

1. Centrifuge the tube at  $12,000 \times g$  for 15 minutes at  $4^{\circ}\text{C}$ .
2. Carefully remove the supernatant.
3. Wash the pellet in 70% ethanol.
4. Centrifuge the tube at  $12,000 \times g$  for 15 minutes at  $4^{\circ}\text{C}$ .
5. Carefully remove the supernatant and air-dry the pellet at room temperature for 30 minutes to remove retained ethanol.
6. Resuspend the pellet in RNase-free water to the desired concentration.

Proceed with the preparation of labeled cDNA and interrogate the arrays according to the manufacturer's instructions.

## QUALITY CONTROL TESTING

The quality of the Universal Mouse Reference RNA is assessed by observing distinct 28S and 18S ribosomal bands on a  $1 \times$  MOPS agarose gel under denaturing conditions. The purity of the RNA is assessed by spectrophotometry ( $A_{260}/A_{280} \geq 1.8$ ). The RNA is then shown to be free of contaminating RNases by incubation in a suitable buffer at  $37^{\circ}\text{C}$  followed by gel analysis against known RNase-free controls. The RNA is further tested functionally by synthesizing labeled cDNA, which is then hybridized to a microarray to examine gene representation and coverage.

## LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.