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The screenshot shows the top navigation bar of the Agilent eArray website. The 'Application Type' is currently set to 'Expression', which is highlighted in a red box. A link 'Switch Application Type' is also visible next to it. Below the navigation bar, there is a search section with radio buttons for 'Microarray', 'Probe Group', and 'Simple Probe'. The 'Microarray' option is selected. There are input fields for 'Microarray Name', 'Species', and 'Design Number', each with a corresponding button ('Select and Add', 'Upload'). To the right, the 'Design Wizards' section offers three options for creating a microarray design, with a 'Next >>' button at the bottom.

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The screenshot shows a dialog box titled 'Switch Application Type' overlaid on the main application page. The dialog has a dropdown menu for 'Select View Type' currently set to 'Expression'. Below it is a checkbox for 'Set As Default View Type'. At the bottom of the dialog are 'Change' and 'Close' buttons. A red arrow points from the 'Switch Application Type' link in the top navigation bar to the dialog box.

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[Next >>](#)

Custom Design Guidance


Search Results: 3 matching results found

Job Type	Job Name	Status	Created Date	Job Position	Action
GeneIntervalSearch	chr21	Complete	21-Dec-2007		Delete

Search Results: 0 matching results found

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The screenshot shows the 'Create a Microarray Design from Existing ProbeGroup(s)' interface. The 'Define Design' step is selected in the left sidebar. The main form includes the following fields:

- Microarray Name:** A text input field, highlighted with a red box.
- Design Format:** A dropdown menu with 'Select Format' and a 'Show Details' link, also highlighted with a red box.
- Folder:** A dropdown menu currently set to 'Agilent'.
- Species:** A dropdown menu currently set to 'H. sapiens'.
- Description:** A text input field.
- Keywords:** A text input field.
- Comments:** A text input field.

The 'Linker Details' section is expanded and contains:

- Append linker to 3' end:** A checkbox with a link to 'info'.
- Linker length:** Two radio button options: 'Make probes of length' (set to 20) and 'Add linker of length' (set to 0).
- Linker sequence:** Three radio button options: 'Use Agilent linker sequence' (selected, showing 'ATAACGGAGGGG' and 'TAA'), 'Use Agilent linker sequence' (showing 'ATAACGGAGGGG'), and 'Use Customer linker sequence'.

At the bottom, there are 'Next >>' and 'Cancel' buttons.

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The screenshot shows the 'Create a Microarray Design from Existing ProbeGroup(s)' interface. The top panel is titled 'Select Probes and Layout Options' and contains a table with columns: Select, Probe Group name, Control Type, Replicate, and Evenly Distribute. Below this table is an 'Add' button, which is circled in red and labeled with a circled '1'. A blue arrow points from this button to the 'Add >' button in the 'Selected Probe Group' panel, which is circled in red and labeled with a circled '3'. The 'Selected Probe Group' panel also contains a 'Done' button, circled in red and labeled with a circled '4'. The 'Probe Group Category List' panel on the left shows a search for 'demo' with two results: 'Demo' (Yayoi, Incomplete) and 'Demonstration' (Agilent_Field, Incomplete). The 'Add >' button is circled in red and labeled with a circled '3'. The 'Done' button is circled in red and labeled with a circled '4'.

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1. Select Species
2. Define Design
3. Layout Probes
4. Create ArrayDesign

Select Probes and Layout Options

Click add button to select ProbeGroup from the system.

Probe Group Details:

Select	Probe Group name	Control Type	Replicate	Evenly Distribute
<input type="checkbox"/>	Demo		5	<input type="checkbox"/>
<input type="checkbox"/>	test21_catalogID		10	<input type="checkbox"/>

Remove Add

Enable MicroArray set:

Fill Microarrays:

Probe Group To Fill Microarray:

MicroArray Statistics:

Number of Microarrays:	8	Percentage Filled:	1.18 %
Total Number of Features:	15744	Number of Available Features:	15559
Number of Agilent Controls:	0	Number of User Controls:	0
Percentage filled using fill array:	1.18 %		

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Create a Microarray Design from Existing ProbeGroup(s)

1. Select Species
2. Define Design
3. Layout Probes
4. Create ArrayDesign

How do you want to save and create your Design?

<input checked="" type="radio"/>	Draft	Saves the design, allowing only you to make changes to it before completing.	More
<input type="radio"/>	Review	Saves the design, enabling users in your Workgroup to make versions of the design before completing.	More
<input type="radio"/>	Complete	Saves the design, this design cannot be edited, nor can it be ordered.	More
<input type="radio"/>	Submit	Saves the design, this design cannot be edited, and is available for ordering <input type="checkbox"/> Design check list	More

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Your Microarray Design has been created successfully.

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Microarray Name	Microarrayset name	Folder Name	Status	Design Number	Created Date	Actions
<input type="checkbox"/> Demo		Yayoi	Complete	0XXXXX	09-Jan-2008	Copy View Delete Submit Download Change Control Type

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Microarray Name: ABC

Status: **Draft** (Draft, Review, Complete)

Folder: [Dropdown]

Species: [Dropdown]

Created By: Yajoi Erika

Design Format: 4 X 44K [Show Details](#)

Control Grid: [IS-45220-4-V1_4x44K_GX_EQC_V20060608](#)

Description: [Text Area]

Keywords: [Text Area]

Attachment: +/-

Comments: [Text Area]

Feature Layout Info: Randomized Ordered for Agilent Scanner

Select	Probe Group name	Control Type	Replicate	Evenly Distribute
<input type="checkbox"/>	013791_M._mulatta (Rhesus) Oligo Microarray_5	[Dropdown]	5	<input type="checkbox"/>
<input type="checkbox"/>	Human 1A Microarray (V2)	[Dropdown]	2	<input type="checkbox"/>

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Download If you have difficulty downloading the desired file, hold down the <Ctrl> key until a File Download dialog box appears. This bypasses pop-up blocking software.

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Design Format: Status:

Search Results: 1 matching results found

<input type="checkbox"/>	Microarray Name	Microarrayset name	Folder Name	Status	Design Number	Created Date	Actions
<input type="checkbox"/>	Demo		Yayoi	Complete	0XXXXX	09-Jan-2008	Copy View Delete Submit Download Change Control Type

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Warning : You are about to submit this Microarray Design. Once saved no changes can be made

Please enter comments in the textarea:

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The screenshot shows a web interface for submitting a microarray design. On the left, a checklist titled "Did you:" contains ten items, each with an unchecked checkbox. A "Help" link is visible in the top right corner of the checklist area. Below the checklist, there are several paragraphs of text providing instructions and disclaimers. At the bottom of this section are "Done" and "Cancel" buttons. On the right, a warning dialog box is displayed with a blue header: "Warning : You are about to submit this Microarray Design. Once saved no changes can be made". Below the header is a text input field containing "ANY TXT|". Underneath the input field is the question "Are you sure you want to continue ?" followed by a checked checkbox and a "Design check list" link. At the bottom of the dialog are "Yes" and "No" buttons.

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 **Agilent Technologies**

You have submitted a design to manufacturing. However, an order needs to be placed for the synthesis to start. To place an order please use the request quote function on the eArray website. After doing so, you will receive a formal quote with further instructions to call in your PO number.

The details of the design are listed below:
 Array Set Name = Demo
 Design Number = 018940
 Slide Format = 8 X 15K
 Application = CGH
 Customer Design ID = NeAT_design
 Slide Layout = IS-15744-8-V1
 Submission Date = 09-Jan-2008 04:06

The details of the Customer are listed below.

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