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

The screenshot shows the 'Search Standard Probes' interface. At the top, there are navigation tabs: Home, Microarray, Probe Group, Probe (selected), My Functions, My Account, and Data. Below these are links for Expression and Switch Application Type. A secondary navigation bar contains links for Search, Upload, Simple Tiling, GE Probe Design, GE Probe Check, Reannotate, and GeneSpring Download. The main section has three radio buttons: Search Standard Probes (selected), Search Exon Probes, and GO Search Standard Probes. The search form includes: Search Type (dropdown menu set to 'All'), Exact Search checkbox, Search Term(s) input field, Species (dropdown menu set to 'Select'), Used in Array Designs input field with a 'Select and Add' link, Folder (dropdown menu set to 'Select'), subfolders checkbox, Used in Probe Groups (dropdown menu set to 'Select') with a 'Select and Add' link, and an 'Include' checkbox. At the bottom, there are 'Search' and 'Reset' buttons, with the 'Search' button highlighted by an orange box.

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The screenshot shows the 'Search Standard Probes' interface. At the top, there are navigation tabs: Home, Microarray, Probe Group, Probe (selected), My Functions, My Account, and Data. To the right, there are links for 'Expression' and 'Switch Application Type'. Below the tabs is a search bar with options: Search, Upload, Simple Tiling, GE Probe Design, GE Probe Check, Reannotate, and GeneSpring Download. The main search area has three radio buttons: 'Search Standard Probes' (selected), 'Search Exon Probes', and 'GO Search Standard Probes'. The search form includes: 'Search Type:' dropdown set to 'Accessions'; 'Search Term(s):' text input with 'AA213559|AA378382|AA4' and an 'Upload' link; 'Species' dropdown set to 'H. sapiens' with an asterisk; 'Used in Array Designs:' text input with a 'Select and Add' link; 'Exact Search' checkbox checked; 'Folder' dropdown set to 'AgilentCatalog' with an asterisk and an 'Include subfolders' checkbox; 'Used in Probe Groups' text input with a 'Select and Add' link. At the bottom, there are 'Search' and 'Reset' buttons.

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
The screenshot shows the Agilent eArray web interface. The 'Probe' tab is selected in the top navigation bar. Below it, the 'Search' option is highlighted in the sub-menu. The 'GO Search Standard Probes' radio button is selected. The search form includes the following fields and options:


- GO Term: * [Select and Add](#)
- Species [Info](#): *
- Folder [Info](#): * Include subfolders
- Used in Probe Groups [Info](#): [Select and Add](#)
- Exact Search: Exact Search

Buttons for 'Search' and 'Reset' are located at the bottom of the form.

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The screenshot shows the Agilent eArray web interface for searching standard probes. The 'Probe' tab is active. The search form includes the following fields and options:

- GO Term:** A text input field containing 'central nervous system de' with a 'Select and Add' link.
- Species:** A dropdown menu showing 'H. sapiens' with an 'info' link.
- Folder:** A dropdown menu showing 'AgilentCatalog' with an 'info' link and an 'Include subfolders' checkbox.
- Used in Probe Groups:** An empty text input field with a 'Select and Add' link.
- Exact Search:** A checkbox labeled 'Exact Search'.
- Buttons:** 'Search' and 'Reset' buttons.

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- Search Type:** A dropdown menu set to 'All'.
- Exact Search:** An unchecked checkbox.
- Search Term(s):** A text input field.
- Folder:** A dropdown menu set to 'AgilentCatalog'.
- Include subfolders:** An unchecked checkbox.
- Species:** A dropdown menu set to 'Select'.
- Used in Probe Groups:** A text input field with a 'Select and Add' link.
- Used in Array Designs:** A text input field with a 'Select and Add' link.

At the bottom of the form, there are two buttons: 'Search' and 'Reset'.

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The screenshot shows the 'Search Exon Probes' interface. At the top, there are navigation tabs: Home, Microarray, Probe Group, **Probe**, My Functions, My Account, and Data. Below these are links for Expression and Switch Application Type. A secondary navigation bar contains links for Search, Upload, Simple Tiling, GE Probe Design, GE Probe Check, Reannotate, and GeneSpring Download. The main search area has three radio buttons: Search Standard Probes, **Search Exon Probes**, and GO Search Standard Probes. The search form includes:

- Search Type:** A dropdown menu set to 'Accessions' with an 'Exact Search' checkbox checked.
- Search Term(s):** A text input field containing 'AA213559|AA378382|AA4' and an 'Upload' link.
- Species:** A dropdown menu set to 'H. sapiens' with an asterisk and an 'Info' link.
- Folder:** A dropdown menu set to 'AgilentCatalog' with an asterisk and an 'Include subfolders' checkbox.
- Used in Array Designs/Used in Probe Group(s):** Two empty text input fields with 'Select and Add' links.

At the bottom of the form are 'Search' and 'Reset' buttons.

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Used in Array Designs/Used in Probe Group(s):

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Home | Microarray | Probe Group | **Probe** | My Functions | My Account | Data | Expression | [Switch Application Type](#)

[Search](#) | [Upload](#) | [Simple Tiling](#) | [GE Probe Design](#) | [GE Probe Check](#) | [Reannotate](#) | [GeneSpring Download](#)

Search Standard Probes Search Exon Probes GO Search Standard Probes

Search Type: Accessions Exact Search Folder [Info](#): AgilentCatalog Include subfolders

Search Term(s): AA213559|AA378382|AA4 [Upload](#)

Species [Info](#): H. sapiens * Used in Probe Groups [Info](#): [Select and Add](#)

Used in Array Designs: [Select and Add](#)

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Search Results: 13 matching results found ②

Filter **Download** **Show Statistics** **Create Probe Group**
Array Calculator **Apply Default Filter**

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Create Probe Group

Probe Group Name	<input type="text"/>	Created Date	01/09/2008
Status	<input checked="" type="radio"/> Incomplete <input type="radio"/> Locked	Description	<input type="text"/>
Created By	Yayoi Fukuoka	Keyword Info	<input type="text"/>
Folder	Agilent		

Search Results: 17 matching results found

<input type="checkbox"/>	Probe ID	Accessions	Gene Name	Gene Symbol
<input type="checkbox"/>	A_16_P35089808	ref NM_005119	entg thyroid hormone receptor associated protein 3	entg THRAP3
<input type="checkbox"/>			entg membrane associated guanine	

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
Application Type: CGH [Switch Application Type](#)

Home Microarray **Probe Group** Probes My Account Site Maintenance

Search | Browse ProbeGroup

Browse WorkGroup By Category

— Category

Browse Catalog By Category

Category

- Applications
- Species

Browse By Folder

ROOT

- Workspace名
- AgilentCatalog

View Probe Group

Status: Select [Filter](#)

Search Results: 2 matching results found

[Compare](#) [Create Microarray](#) [Share](#) [Move](#)

<input type="checkbox"/>	Probe Group Name ▲	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	AllDemo_Rat	99328	true	LSCA_JapanSupportSpace	Incomplete	12-Feb-2008	Copy Edit View Delete Download
<input type="checkbox"/>	test	14680	true	LSCA_JapanSupportSpace	Incomplete	12-Feb-2008	Copy Edit View Delete Download

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