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The screenshot shows the Agilent eArray design tool interface. At the top, there are two radio buttons: "Base Composition Methodology" (selected) and "Tm Matching Methodology". Below this, the interface is divided into several sections:

- Probe Details:** Design Job Name (text input), Probe Length (60 bp), Probes per Target (1).
- Masking:** Vector (checked), Repeat (checked).
- Probe Orientation:** Sense (selected), Antisense.
- Design Options:** Best Probe Methodology (selected), Best Distribution Methodology, Design with 3' Bias (checked).
- Target File Details:** Species (H. sapiens), Upload in FASTA Format (Browse...), Upload as Genbank Accessions (Browse...).
- Transcriptome Details:** Use Target File as Transcriptome, Select Agilent-provided Transcriptome (by Species) (H. sapiens), Upload Transcriptome File (Browse...).

At the bottom, there are "Submit" and "Cancel" buttons.

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Probe Details	
Design Job Name	<input type="text"/>
Probe Length	<input type="text" value="60"/> bp
Probes per Target <a href="#">Info</a>	<input type="text" value="1"/>
Masking	<input checked="" type="checkbox"/> Vector <a href="#">Info</a> <input checked="" type="checkbox"/> Repeat <a href="#">Info</a>
Probe Orientation	<input type="radio"/> Sense <a href="#">Info</a> <input checked="" type="radio"/> Antisense <a href="#">Info</a>
Design Options	<input checked="" type="radio"/> Best Probe Methodology <a href="#">Info</a> <input type="radio"/> Best Distribution Methodology <a href="#">Info</a>
	<input checked="" type="checkbox"/> Design with 3' Bias <a href="#">Info</a>
<input type="button" value="Submit"/>	

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
Probe Details	
Design Job Name	<input type="text"/>
Probe Length	60 bp
Probes per Target <a href="#">Info</a>	1
Masking	<input checked="" type="checkbox"/> Vector <a href="#">Info</a> <input checked="" type="checkbox"/> Repeat <a href="#">Info</a>
Probe Orientation	<input type="radio"/> Sense <a href="#">Info</a> <input checked="" type="radio"/> Antisense <a href="#">Info</a>
Design Options	<input checked="" type="radio"/> Best Probe Methodology <a href="#">Info</a> <input type="radio"/> Best Distribution Methodology <a href="#">Info</a>
<input checked="" type="checkbox"/> Design with 3' Bias <a href="#">Info</a>	
<input type="button" value="Submit"/>	

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Probe Details	
Design Job Name	<input type="text"/>
Probe Length	<input type="text" value="60"/> bp
Probes per Target <a href="#">Info</a>	<input type="text" value="1"/>
Masking	<input checked="" type="checkbox"/> Vector <a href="#">Info</a> <input checked="" type="checkbox"/> Repeat <a href="#">Info</a>
Probe Orientation	<input checked="" type="radio"/> Sense <a href="#">Info</a> <input type="radio"/> Antisense <a href="#">Info</a>
Design Options	<input checked="" type="radio"/> Best Probe Methodology <a href="#">Info</a> <input type="radio"/> Best Distribution Methodology <a href="#">Info</a>
<input checked="" type="checkbox"/> Design with 3' Bias <a href="#">Info</a>	
<input type="checkbox"/> Allow Probes to be Trimmed <a href="#">Info</a>	
Preferred Probe Tm <a href="#">Info</a>	<input type="text" value="80.0"/> °C

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The screenshot shows a software interface with two main sections: 'Target File Details' and 'Transcriptome Details'. In the 'Target File Details' section, the 'Species' dropdown is set to 'H. sapiens'. The 'Upload in FASTA Format' radio button is selected, and the file path 'C:\Documents and Sett...' is entered in the text box, with a 'Browse...' button to its right. The 'Upload as Genbank Accessions' radio button is unselected. In the 'Transcriptome Details' section, the 'Select Agilent-provided Transcriptome (by Species)' radio button is selected, and the 'Species' dropdown is also set to 'H. sapiens'. There is an empty text box and a 'Browse...' button below this section. An orange arrow points from the 'Browse...' button in the 'Transcriptome Details' section to the 'Browse...' button in the 'Target File Details' section.

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Target File Details	
Species	H. sapiens
<input checked="" type="radio"/> Upload in FASTA Format <a href="#">Info</a>	C:\Documents and Sett... <a href="#">Browse...</a>
<input type="radio"/> Upload as Genbank Accessions <a href="#">Info</a>	<a href="#">Browse...</a>
Transcriptome Details	
<input type="radio"/> Use Target File as Transcriptome <a href="#">Info</a>	
<input checked="" type="radio"/> Select Agilent-provided Transcriptome (by Species) <a href="#">Info</a>	H. sapiens
<input type="radio"/> Upload Transcriptome File <a href="#">Info</a>	<a href="#">Browse...</a>

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<input checked="" type="radio"/> Base Composition Methodology <a href="#">Info</a>		<input type="radio"/> Tm Matching Methodology <a href="#">Info</a>	
<b>Probe Details</b>		<b>Target File Details</b>	
Design Job Name	<input type="text"/>	Species	H. sapiens <input type="text"/>
Probe Length	60 <input type="text"/> bp	<input checked="" type="radio"/> Upload in FASTA Format <a href="#">Info</a>	<input type="text"/> <input type="button" value="Browse..."/>
Probes per Target <a href="#">Info</a>	1 <input type="text"/>	<input type="radio"/> Upload as Genbank Accessions <a href="#">Info</a>	<input type="text"/> <input type="button" value="Browse..."/>
Masking	<input checked="" type="checkbox"/> Vector <a href="#">Info</a> <input checked="" type="checkbox"/> Repeat <a href="#">Info</a>	<b>Transcriptome Details</b>	
Probe Orientation	<input checked="" type="radio"/> Sense <a href="#">Info</a> <input type="radio"/> Antisense <a href="#">Info</a>	<input type="radio"/> Use Target File as Transcriptome <a href="#">Info</a>	
Design Options	<input checked="" type="radio"/> Best Probe Methodology <a href="#">Info</a> <input type="radio"/> Best Distribution Methodology <a href="#">Info</a>	<input checked="" type="radio"/> Select Agilent-provided Transcriptome (by Species) <a href="#">Info</a>	H. sapiens <input type="text"/>
<input checked="" type="checkbox"/> Design with 3' Bias <a href="#">Info</a>		<input type="radio"/> Upload Transcriptome File <a href="#">Info</a>	<input type="text"/> <input type="button" value="Browse..."/>
<input checked="" type="button" value="Submit"/> <input type="button" value="Cancel"/>			

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Design Job Submitted.
You will receive an e-mail when your design job is completed.
<input type="button" value="Close"/>

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
Search Results: 1 matching results found

Job ID	Job Name	Status	Design Type	Design Method	Creation Date	Position in Queue	Actions
PI377492179	Demo0207	Designing	Select Best Probes	Base Composition Method	06-Feb-2008	Not in Queue	

Refresh

Job ID	Job Name	Status	Design Type	Design Method	Creation Date	Position in Queue
PI377492179	Demo0207	Completed	Select Best Probes	Base Composition Method	06-Feb-2008	Not in Queue

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

Errors in file submitted for upload

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Agilent Technologies

Dear [redacted]

Your file submitted for upload on Mon Feb 11 19:53:10 MST 2008 could not go through because of errors. An error report is attached. To view the error messages within the help file, place your cursor over the red text.

Sincerely,  
Agilent Technologies eArray website

Line Number	ProbeID	Sequence	TargetID	Accessions	Gene
5	A.65_P100870	CGAGCAGAGAAGCAAGACTGACCAGACTCTCGGGCTGGAGCAGTGAAGCCTGCTAGTGTAT	ukuster4L019099	ref[NM_053252]gb AF305087 gb AF305088	ILOC1
8	A.65_P100474	AGTCGCTAATTGTACTGGTTTAAGTGTGCAATTAGTAGTGTCTCCCCAGGTAGATTGTAAG	ukuster4L002740	ref[NM_010923]ref[NM_180960]gb AB004048 gb AK003004	ILOC1

Line number 8: For Column : Accessions, Field : ref[NM\_010923ref[NM\_180960]gb|AB004048]gb|AK003004 is not valid as per File Format.

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
Transcriptome Details

- Use Target File as Transcriptome [Info](#)
- Select Agilent-provided Transcriptome (by Species) [Info](#)
- Upload Transcriptome File [Info](#)

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Refresh

Job ID	Job Name	Status	Design Type	Design Method	Creation Date	Position in Queue	Actions
PI377492179	Demo0207	Completed	Select Best Probes	Base Composition Method	06-Feb-2008	Not in Queue	<a href="#">Delete</a>   <a href="#">View Design</a>   <a href="#">Download</a>   <a href="#">Create Probe Group</a>

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**Probe Group Creation Job Submitted**

You will receive an email when your probe group has been created

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Dear [REDACTED]

Your probe design job Demo0207\_2 with ID PI377493211 has completed. Please log into eArray to view the results of the design. You can save the probes resulting from this design process into a probe group as part of the array creation process.

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
A screenshot of a search form. The navigation menu at the top includes Home, Microarray, Probe Group (highlighted), Probe, My Functions, and My Account. Below the menu is a blue bar with 'Search | Browse ProbeGroup'. The form has two input fields: 'Probe Group Name:' and 'Keyword [Info](#):'.

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The screenshot displays the 'View Probe Group' interface. On the left, the 'Browse By Folder' section is active, showing a tree view with 'ROOT' and a sub-folder 'Workgroup名' (highlighted with a red box). The main area shows a search filter for 'Status: Select' and a table of results. The table has columns for 'Probe Group Name', 'No. of Probes', 'High Density', 'Folder', 'Status', 'Created Date', and 'Actions'. One row is visible for 'Demo0207\_2' with 10 probes, low density, in the 'LSCA\_JapanSupportSpace' folder, with a status of 'Incomplete' and a creation date of '07-Feb-2008'. Action buttons like 'Compare', 'Create Microarray', 'Share', and 'Move' are present above and below the table.

<input type="checkbox"/>	Probe Group Name	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	Demo0207_2	10	false	LSCA_JapanSupportSpace	Incomplete	07-Feb-2008	<a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">View</a>   <a href="#">Delete</a>   <a href="#">Download</a>

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The screenshot displays the 'View Probe Group' interface. On the left, a sidebar shows a tree view under 'AgilentCatalog' with 'Workgroup名' selected. The main content area shows a search filter for 'Status' set to 'Select'. Below the filter, it indicates 'Search Results: 1 matching results found'. A table lists the probe group details:

<input type="checkbox"/>	Probe Group Name	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	Demo0207_2	10	false		Incomplete	07-Feb-2008	<a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">View</a>   <a href="#">Delete</a>   <a href="#">Download</a>

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**Edit Probe Group**

Probe Group Name Demo Created Date 07/27/2009

Status  Incomplete  Locked Description info

Created By  Keyword info

Folder

Save Probe Group Cancel

Search Results: 8 matching results found

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TDT  
 FASTA  
 COMPLETE  
 MINIMAL  
 BED

**DOWNLOAD** **CLOSE**

If you have difficulty downloading the desired file, hold the <Ctrl> key until a File Download dialog box appears bypasses pop-up blocking software.

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Attribute	TDT	FASTA	COMPLETE	MINIMAL	BED
ProbeID	●	●	●	●	●
Sequence	●	●	●	●	
TargetID	●		●		
Species	●				
GeneName	●				
GeneSymbol	●		●		
Description	●		●		
ControlType	●				
Accessions	●		●		
ProbeGroups	●				
Status	●				
ValidationMethod	●				
Chromosomal Location	●		●		●
CytoBand	●				
GoIDs	●				

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
The screenshot shows the 'Probe Details' form. It has a blue header with the text 'Probe Details'. Below the header, there are several fields and options:

- 'Design Job Name': A text input field.
- 'Species': A dropdown menu currently showing 'H. sapiens'.
- 'Upload Probe File': A radio button that is selected (indicated by a filled circle). To its right is a text input field and a 'Browse...' button.
- 'Select a Probe Group': A radio button that is unselected (indicated by an empty circle). To its right is a text input field.

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## Select Agilent-provided Transcriptome (by Species):

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

Job ID	Job Name	Status	Design Type	Creation Date ▼	Position in Queue	Actions
PI417509223	Demo	In Queue	Characterization	29-Jul-2009	Job 1 of 1	<a href="#">Delete</a>
PI417509223	Demo	Completed	Characterization	29-Jul-2009	NA	<a href="#">Delete</a>   <a href="#">View Design</a>   <a href="#">Download</a>   <a href="#">Create Probe Group</a>

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