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The screenshot shows the Agilent eArray web interface. At the top, there is a navigation bar with tabs for Home, Microarray, Probe Group, Probe, My Functions, My Account, and Data. The 'Expression' application type is selected and highlighted with a red box, with a link to 'Switch Application Type'. Below the navigation bar is a 'Releases & Announcements' section. The main content area is divided into two panels. The left panel is titled 'Search' and has three radio buttons: 'Microarray' (selected), 'Probe Group', and 'Probe'. Below the radio buttons is a form with three input fields: 'Microarray Name:', 'Species:', and 'Design ID:'. The 'Species:' field has a 'Select and Add' link, and the 'Design ID:' field has an 'Upload' link. At the bottom of the form are 'Search' and 'Reset' buttons. The right panel is titled 'Design Wizards' and has 'Refresh' and 'View All' links. It contains three radio buttons: 'Create a Microarray Design by Uploading Probes Info', 'Create a Microarray Design from Existing Probe Group(s) Info' (selected and highlighted with a red box), and 'Create a Microarray Design from Target Transcripts Info'. Below the radio buttons is a green 'Next >>' button and a green 'Custom Design Guidance' button. At the bottom of the right panel is a search results box containing the text 'Search Results: No matching results found'.

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Create a Microarray Design from Existing Probe Group(s)

1. Select Species
2. Define Design
3. Layout Probes
4. Create Microarray Design

Species is not required for Expression microarray designs.

Select Array Type : [Info](#) Standard Exon

Select Species :

You have selected to design an Exon Array. These arrays are not compatible with standard oligo-dT labeling kits. Agilent recommends the use of the Low-Input QuickAmp WT Labeling Kit (p/n 5190-2943, 5190-2944, 5190-2945 and 5190-3386). For sample processing, refer to the Agilent One Color Microarray-Based Exon Analysis Protocol or the Agilent Two Color Microarray-Based Exon Analysis Protocol versions 1.0 or later.

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Create a Microarray Design from Existing Probe Group(s)

1. Select Species	Design Details	
2. Define Design	Microarray Name: <input type="text"/>	
3. Layout Probes	Design Format: <input type="text" value="Select Format"/> Show Details	Control Grid: <input type="text"/>
4. Create Microarray Design	Folder: <input type="text" value="Agilent"/>	Description: <input type="text"/>
	Species: <input type="text" value="H. sapiens"/>	Keywords: <input type="text"/>
	Attachment: <input type="text" value="+/-"/>	Comments: <input type="text"/>

Linker Details:	
<input type="checkbox"/> Append linker to 3' end info	
Linker Length:	<input checked="" type="radio"/> Make probes of length <input type="text" value="60"/>
	<input type="radio"/> Add linker of length <input type="text" value="0"/>

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Linker Details:	
<input type="checkbox"/> Append linker to 3' end info	
Linker Length:	<input checked="" type="radio"/> Make probes of length <input type="text" value="60"/>
	<input type="radio"/> Add linker of length <input type="text" value="0"/>
Linker Sequence:	<input checked="" type="radio"/> Use Agilent linker sequence <input type="text" value="TATCCTACTATAC"/>
	<input type="radio"/> Use customer linker sequence <input type="text"/>

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Microarray Statistics:			
Number of Microarrays per Slide:	4	Percentage Filled by Selected Probe Groups:	2.75542%
Number of Slides:	1	Number of Available Features:	175896
Total Number of Features:	180880	Number of User Controls:	0
Number of Agilent Controls:	4974	Number of Features Occupied by Replicate Probes:	0
Total Percentage Filled:	2.75542 %		

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Create a Microarray Design from Existing Probe Group(s)	
1. Select Species	Create microarray design with this status:
2. Define Design	
3. Layout Probes	
4. Create Microarray Design	
	<input checked="" type="radio"/> Draft Allows only you to edit the design. Later, you can change the status to any of the others. More
	<input type="radio"/> Review Lets you and the other users in your workgroup make changes to the design and save new versions of it. More
	<input type="radio"/> Complete Prevents further edits to the design. You must submit the design to Agilent Manufacturing before you can request a quote. More
	<input type="radio"/> Submitted Submits the design to Agilent Manufacturing, and lets you request a quote for the design. The design cannot be edited. <input type="checkbox"/> Design Checklist More

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
Your Microarray Design has been created successfully.

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<input type="checkbox"/>	Microarray Name ▲	Microarray Set Name	Folder Name	Status	Design ID	GPL Identifie	Created Date	Actions
<input type="checkbox"/>	de			Draft			21-Dec-2010	Copy Edit View Delete
<input type="checkbox"/>	delete			Complete	OXXXXX		13-Dec-2010	Copy View Download

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Microarray Name: de Species: H. sapiens Design Format: 4 X 44K Show Details

Status: Draft Description: Comments: Keywords:

Folder: Attachment:

Control Grid: IS-45220-4-V1 4x44K GX EQC ERCC V20100113

Feature Layout: Randomized Customer Specified

Number of Available Features: 43333
Total Percentage Filled: 4.17293
Number of Features Occupied by Replicate Probes: 500

Select	Probe Group Name	Control Type	Replicate
<input type="checkbox"/>	demo demo	biological	1

Click Add to select probe group(s). Remove Add

Replicate Probe Group Details : Number of Replicate Probe Group(s) Selected :1

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
The screenshot shows the Agilent eArray web interface. At the top, there is a navigation menu with tabs for Home, Microarray, Probe Group, Probe, My Functions, My Account, and Data. To the right of the menu, there are links for Expression and Switch Application Type. Below the menu is a search bar and a list of links: Search, Create Microarray, Browse Microarray Designs, Browse Microarray Sets, Express Array, Notification, Array Calculator, and Request Quote. A blue warning banner reads: "Warning : You are about to save this microarray design with a status of Complete. Once you save it, no additional changes can be made." Below the warning, there is a text input field with the label "Please type comments here: :". To the right of the input field is a vertical scrollbar. Below the input field, the text "Are you sure you want to continue?" is displayed. At the bottom of the dialog, there are two buttons: "Yes" and "No".

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<input type="checkbox"/>	Category	File Type
<input type="checkbox"/>	BED	BED
<input type="checkbox"/>	EXTERNALFULLGEML	GEML_1.0
<input type="checkbox"/>	EXTERNALFULLGEML2	
<input type="checkbox"/>	FASTA	
<input type="checkbox"/>	GAL	
<input type="checkbox"/>	GENELIST	List
<input type="checkbox"/>	IMAGENE	Image
<input type="checkbox"/>	SEQUENCELIST	List
<input type="checkbox"/>	TDT	TDT File

Download If you have difficulty downloading the desired file, Select and hold the <Ctrl> key until you get File DownLoad Window. This will bypass pop-up blocking software.

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Design Format: Status: [Filter](#)

Search Results: 1 matching results found

[Share](#) [Move](#)

<input type="checkbox"/>	Microarray Name ▲	Microarrayset name	Folder Name	Status	Design Number	Created Date	Actions
<input type="checkbox"/>	Demo		Yayoi	Complete	0XXXXX	09-Jan-2008	Copy View Delete Submit Download Change Control Type

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Warning : You are about to submit this Microarray Design. Once saved no changes can be made

Please enter comments in the textarea :

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The screenshot shows a checklist titled "Did you:" with the following items:

- select the appropriate format (e.g. 244K etc.)?
- include the necessary controls not already incorporated in the design?
- include the appropriate probes that target the appropriate species?
- check the quality and form (e.g. length, linkers, etc.) of any uploaded probes on this microarray?
- for your intended use, include only those uploaded probes for which you agree to take full responsibility?
- include the appropriate probe groups with the correct probes?
- fill all the unused features on the microarray you wanted filled?
- include the required technical replicates and/or experimental replicates?
- create a design that works with your equipment, kits, and protocols?

Below the checklist is a warning dialog box with a blue header: "Warning : You are about to submit this Microarray Design. Once saved no changes can be made". It contains a text area with "ANY TXT", the question "Are you sure you want to continue?", a checked radio button for "Design check list", and "Yes" and "No" buttons.

At the bottom of the checklist area are "Done" and "Cancel" buttons.

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 **Agilent Technologies**

You have submitted a design to manufacturing. However, an order needs to be placed for the synthesis to start. To place an order please use the request quote function on the eArray website. After doing so, you will receive a formal quote with further instructions to call in your PO number.

The details of the design are listed below:
Array Set Name = Demo
Design Number = 018940
Slide Format = 8 X 15K
Application = CGH
Customer Design ID = NeAT_design
Slide Layout = IS-15744-8-V1
Submission Date = 09-Jan-2008 04:06

The details of the Customer are listed below:

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