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
The screenshot shows the Agilent eArray web interface. At the top, there is a navigation bar with tabs: Home, Microarray, Probe Group, **Probes**, My Account, and Site Maintenance. Below this is a secondary navigation bar with **Search**, Upload, and Genomic Tiling. The main content area has three radio buttons: Simple Probe Search, Chromosomal Location Search, and **HD Probe Search**. The HD Probe Search section contains several input fields: Search Name, Species (Select), Build Number, Filters (Average Spacing), and Genomic Intervals. There are also checkboxes for TM Filtered Probes and HM Filtered Probes, and buttons for Search and Reset.

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Simple Probe Search Chromosomal Location Search **HD Probe Search**

Search Name: Demo Select HD Search By: Genomic Intervals

Species: R. norvegicus Extended Interval Boundary Info: 5' Base Pairs: 0

Build Number: R. norvegicus, UCSC m4, HGSC Version 3.4, i 3' Base Pairs: 0

Filters Info: Total Probes: 100000 Genomic Intervals Info: chr1:1-500001|chr1:24494200 Upload

Include Regions: All Prefer Catalog Probes Info:

TM Filtered Probes HM Filtered Probes

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CGH

Your HD query has been submitted to the search queue.

Click on refresh button on next page to check whether the search has been completed.

When your query is complete, you will be able to view the results by selecting the 'View' hyperlink on the next page.

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Search Name	Search Type	Status	Job Position	Submitted Date	Actions
Demo_RatCGH	GeneIntervalSearch	Not Started	Job 1 of 1 Jobs	12-Feb-2008	Delete Search Criteria

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Status	Job Position	Submitted Date	Actions
Complete		12-Feb-2008	View Delete Search Criteria

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Refresh

Search Name	Search Type	Status	Job Position	Submitted Date	Actions
Demo	GeneIntervalSearch	Complete		12-Feb-2008	View Delete Search Criteria

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HD Summary Result categories

- **Length** – Total interval length considered for search result
- **RMLength** – Number of probes before removing redundancy.
- **Total Number of Probes** – Total number of resultant probes before filtering.
- **Average probes per 1000 bp** - The average number of probes retrieved per 1000 base pairs
- **Probes per Interval** - Number of probes per interval provided in search criteria.
- **Exonic Probes** - Number of Exonic type probes in final result.
- **Percent Exonic Probes** - Fraction of retrieved probes that are Exonic type.
- **Intragenic Probes** - Number of Intragenic type probes in final result.
- **Percent Intragenic** - Fraction of retrieved probes that are Intragenic type.
- **TM filtered probes** - Number of TM filtered probes in the final result
- **HM filtered probes** - Number of HM filtered probes in the final result

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HD Detailed Result categories

- **Interval Value** - Chromosomal location range or cytoband, or annotation entered as search criteria.
- **Chromosome** - Chromosome values of given interval.
- **Start** - Interval start position.
- **Length** - Total interval length considered for search result
- **RMLength** - Number of probes before removing redundancy.
- **Total Number of Probes** - Total number of resultant probes before filtering.
- **Intragenic Probes** - Number of Intragenic type probes in final result.
- **Exonic Probes** - Number of Exonic type probes in final result.
- **Average probes per 1000 bp** - The average number of probes retrieved per 1000 base pairs
- **Probes per Interval** - Number of probes per interval provided in search criteria.
- **TM filtered probes** - Number of TM filtered probes in the final result
- **HM filtered probes** - Number of HM filtered probes in the final result

Agilent Technologies

HD_Probesの検索

Step1_HD-CGH/HD-ChIPの検索閲覧

Summary Result | Detail Result | **Chromosome View**

Probe List

Search Results: 99328 matching results found

Pages [(1)First] 1 2 3 4 5 [Next>>] [Last (9933)]

Internal ID	Validated	Probe Type	Species	Gene Name	Gene Symbol	Accessions	Chromosomal Location	Cytoband	ProbeID	Exonic	Intragenic	IS TM Filtered	IS HM Filtered
A_86_P00250304	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100000061-100000121	rn q22	A_86_P00250304	0	3	0	0
A_86_P00250309	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100001091-100001147	rn q22	A_86_P00250309	0	3	0	0
A_86_P00250318	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100002957-100003017	rn q22	A_86_P00250318	0	3	0	0
A_86_P00250321	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100003813-100003866	rn q22	A_86_P00250321	0	3	0	0
A_86_P00250325	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100005266-100005324	rn q22	A_86_P00250325	0	3	0	0
A_86_P00250329	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100006881-100006921	rn q22	A_86_P00250329	0	3	0	0
A_86_P00250332	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100008110-100008170	rn q22	A_86_P00250332	0	3	0	0
A_86_P00250337	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100009638-100009696	rn q22	A_86_P00250337	0	3	0	0
A_86_P00250342	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100010726-100010786	rn q22	A_86_P00250342	0	3	0	0
A_86_P00250346	0	CGH	R. norvegicus	entg NEL-like 1 (chicken)	entg Nell1	ref NM_031069	chr1:100011860-100011914	rn q22	A_86_P00250346	0	3	0	0

Pages [(1)First] 1 2 3 4 5 [Next>>] [Last (9933)]

[Close](#)
[New HD Search](#)
[Create Probe Group](#)
[Download](#)

Note:- Your Search will be removed from the system after 48 Hrs!

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Refresh					
Search Name	Search Type	Status	Job Position	Submitted Date	Actions
Demo	GeneIntervalSearch	Complete		12-Feb-2008	View Delete Search Criteria

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Create a Probe Group			
Probe Group Name:	<input type="text"/>	Created Date:	02/12/2008
Type:	CGH	Status:	<input checked="" type="radio"/> Incomplete <input type="radio"/> Locked
Description:	<input type="text"/>	No. of Probes:	99328
Created By:	YAYOI FUKUOKA	Keywords:	<input type="text"/>
Folder:	LSCA_JapanSupportSpac		
<input type="button" value="Create Probe Group"/>		<input type="button" value="Cancel"/>	

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Probe Group Name: AllDemo_Rat

Your Create Probe Group Request has been submitted.

After your request has been processed, you will receive a notification about the status through email.

You will be able to search for the Probe Group only after you receive the email indicating successful creation of the Probe group.

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Application Type: CGH [Switch Application Type](#)

Home Microarray **Probe Group** Probes My Account Site Maintenance

Search | Browse ProbeGroup

Browse WorkGroup By Category

Category

Browse Catalog By Category

Category

- Applications
- Species

Browse By Folder

- ROOT
 - Workspace名**
 - AgilentCatalog

View Probe Group

Status:

Search Results: 2 matching results found

<input type="checkbox"/>	Probe Group Name	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	AllDemo_Rat	99328	true	LSCA_JapanSupportSpace	Incomplete	12-Feb-2008	Copy Edit View Delete Download
<input type="checkbox"/>	test	14680	true	LSCA_JapanSupportSpace	Incomplete	12-Feb-2008	Copy Edit View Delete Download

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Home Microarray **Probe Group** Probes My Account Site Maintenance Application Type: CGH [Switch Application Type](#)

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Compare Create Microarray Share Move

<input type="checkbox"/>	Probe Group Name ▲	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	AllDemo_Rat	99328	true	LSCA_JapanSupportSpace	Incomplete	12-Feb-2008	Copy Edit View Delete Download
<input type="checkbox"/>	test	14680	true	LSCA_JapanSupportSpace	Incomplete	12-Feb-2008	Copy Edit View Delete Download

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