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
A screenshot of a web form titled 'Probe Parameter Details'. The form has a blue header. Below the header, there are several input fields and options. The 'Species:' field is a dropdown menu with 'H. sapiens' selected. Below it is a checkbox labeled 'Remove replicate probes from upload' which is checked. At the bottom, there are three radio button options for 'Probe Precedence': 'Overwrite matching probes', 'Skip matching probes', and 'Cancel upload if any probes already exist'. The 'Overwrite matching probes' option is selected.

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
### Cancel upload if any probes already exist

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Upload Probe File Details	
Upload Type:	<input type="radio"/> Upload Probes Only <a href="#">Info</a> <input checked="" type="radio"/> Create New Probe Group <a href="#">Info</a> <input type="text" value="Demo"/>
Upload File:	<input type="text"/> <input type="button" value="Browse..."/>
File Format: <a href="#">Info</a>	COMPLETE ▾
File Type:	TDT ▾

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The screenshot shows two panels: 'Probe Parameter Details' and 'Upload Probe File Details'. In the 'Upload Probe File Details' panel, the 'Next' button is highlighted with an orange box.

Probe Parameter Details	Upload Probe File Details
Species: <a href="#">Info</a>   H. sapiens	Upload Type: <input type="radio"/> Upload Probes Only <a href="#">Info</a>
<input checked="" type="checkbox"/> Remove replicate probes from upload <a href="#">Info</a>	<input checked="" type="radio"/> Create New Probe Group <a href="#">Info</a>   Demo
Probe Precedence: <input checked="" type="radio"/> Overwrite matching probes <a href="#">Info</a>	Upload File: C:\Documents and Settings   Browse...
<input type="radio"/> Skip matching probes <a href="#">Info</a>	File Format: <a href="#">Info</a>   COMPLETE
<input type="radio"/> Cancel upload if any probes already exist	File Type: TDT
<input type="button" value="Next"/>	<input type="button" value="Cancel"/>

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The screenshot shows the 'Define Uploaded File Columns' dialog. A red circle highlights the 'My uploaded file contains' checkbox. Red arrows point from the text above to the 'Preview' table and the checkbox.

Define Uploaded File Columns	
Preview	
ProbeID	Sequence
A_52_P616047	TGTCCGAAGTGTCTCAGTGCCAACCAACGGAATGGCGAAGAACGGCAGTGAAGCAGACA
A_51_P288390	CCCTGAGTCCACATGCACACCTCAGCAGCCTATATTTATTTGTTTCTGGTTTAGCCTCAA
A_52_P276815	TCTCATCATTGTACGAGAAGTTTCTTCATAAGCACTATAGACAAAAACAATGGGGGAAGA
A_52_P363216	CAGATGATGTAGCAAGCCTTAGGATGGTTCATGATCAGTAGGAAGCAAATCTAAGACCT
A_52_P594049	GAAACAGCTGGAGGGAGAAGTGGGCAGTGCTGCTCTTTGCACTGTTTTGTTTTGTTTTA
A_52_P429944	TATCATTCTGAGCCAGCTCCACTTATTCTTGGTCCAAAGAGCAGTTTGTTCCATGTGG
Match columns from preview by selecting the appropriate columns from following dropdowns	
ProbeID	Sequence
<input type="checkbox"/> My uploaded file contains	ProbeID
	Sequence
	TargetID
	Accessions
	GeneSymbols
	Description
	ChromosomalLocation
	Ignore
<input type="button" value="Upload"/> <input type="button" value="Cancel"/>	

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File successfully submitted to the upload queue. You will be notified of file errors or upload success via email. You will be able to search for the uploaded probes/probe group only after you receive the email indicating upload success.

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The screenshot shows the Agilent eArray web interface. At the top, there are navigation tabs: Home, Microarray, Probe Group, Probes, My Account, and Site Maintenance. Below the tabs is a search section with radio buttons for Microarray (selected), Probe Group, and Simple Probe. There are input fields for Microarray Name, Species, and Design Number, along with buttons for Search and Reset. Below the search section is a 'Pending Jobs' section with a 'Refresh' and 'View All' link. A table displays search results for pending jobs. The first row is highlighted with a red border.

Job Type	Job Name	Status	Created Date	Job Position	Action
ProbeUpload	ID378389145	Upload Pending	13-Feb-2008	Job 1 of 1 Jobs	Delete
ProbeDesign	short2	Probe Design Completed	08-Feb-2008		Delete

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Errors in file submitted for upload

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**Agilent Technologies**

Dear [redacted]

Your file submitted for upload on Mon Feb 11 19:53:10 MST 2008 could not go through because of errors. error report is attached. To view the error messages within the help file, place your cursor over the red text.

Sincerely,  
Agilent Technologies eArray website

Line Number	ProbeID	Sequence	TargetID	Accessions	Gene
5	A_65_P100170	OSAGCAGAAAGCAAGACTGAACAGACTOTGGGOTGGAGCAGTGAGOOTGOTAGTGTTTAT	uc001ter4LJ019009	ref NM_053252 gb AF305087 gb AF305088	#LOC100287304
8	A_65_P100474	AGTOGOTAAATTGTACTGGTTTAAAGTGTGCATTAGTAGTGTOTCCDAGOTAGATTGTAAG	uc001ter4LJ027400	ref NM_010923 ref NM_180960 gb AB004048 gb AK003004	#Next

Line number 8: For Column : Accessions, Field : ref|NM\_010923|ref|NM\_180960|gb|AB004048|gb|AK003004 is not valid as per File Format.

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The screenshot shows the 'View Probe Group' interface. The left sidebar has a 'Workspace名' label highlighted. The main content area shows a table of probe groups:

<input type="checkbox"/>	Probe Group Name	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	Demo0213	10	false	LSCA_JapanSupportSpace	Incomplete	13-Feb-2008	<a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">View</a>   <a href="#">Delete</a>   <a href="#">Download</a>
<input type="checkbox"/>	catalogMouse	10	false	LSCA_JapanSupportSpace	Incomplete	11-Feb-2008	<a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">View</a>   <a href="#">Delete</a>   <a href="#">Download</a>

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Home Microarray **Probe Group** Probes My Account Site Maintenance

Application Type: **Expression** [Switch Application](#)

Search | Browse ProbeGroup

**Browse WorkGroup By Category**

—Category

**Browse Catalog By Category**

[-] Category

- [-] Applications
- [-] Species

**Browse By Folder**

- [-] ROOT
  - Workspace名**
  - AgilentCatalog

**View Probe Group**

Status:

Search Results: 4 matching results found

<input type="checkbox"/>	Probe Group Name	No. of Probes	High Density	Folder	Status	Created Date	Actions
<input type="checkbox"/>	Demo0213	10	false	LSCA_JapanSupportSpace	Incomplete	13-Feb-2008	<a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">View</a>   <a href="#">Delete</a>   <a href="#">Download</a>
<input type="checkbox"/>	catalogMouse	10	false	LSCA_JapanSupportSpace	Incomplete	11-Feb-2008	<a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">View</a>   <a href="#">Delete</a>   <a href="#">Download</a>

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Attribute	TDT	FASTA	COMPLETE	MINIMAL	BED
ProbeID	●	●	●	●	●
Sequence	●	●	●	●	
TargetID	●		●		
Species	●				
GeneName	●				
GeneSymbol	●		●		
Description	●		●		
ControlType	●				
Accessions	●		●		
ProbeGroups	●				
Status	●				
ValidationMethod	●				
Chromosomal Location	●		●		●
CytoBand	●				
GoIDs	●				

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