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


ChIP-on-chip to study DNA-Protein interaction

R. Young

DNA-protein interactions are studied within the context of the natural cell environment (in vivo).

Combining expression and binding of the transcription machinery data will help:

- Identify genomic targets and regulatory pathways
- Understand gene expression profiles of disease relevance
- Identify biomarker and toxicant signatures

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
Whole genome arrays
Promoter arrays

Scatter plot
(ChIP / reference)

Promoters bound by
Transcriptional regulators

- 1) Boyer LA, et al: Cell (2005) 122: 947-956
- 2) Lee TI, et al: Cell (2006) 125: 301-313
- 3) Boyer LA, et al: Nature (2006) 441: 349-353
- 4) Ji H, et al: Nucleic Acids Res (2006) 34: e146.
- 5) Tesar, PJ, et al: Nature (2007) 448: 196-199

- 6) Guenther MG et al : Cell (2007) 130, 77-88
- 7) Cole MF et al : Genes and Dev. (2008) 22: 746-755
- 8) Endo M et al : Development (2008) 135, 1513-1524
- 9) Jaenisch R et al : Cell (2008) 132: 567-582

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➤ Impact of DNA methylation:


- Protein recruitment leading to compact, silent chromatin
- Gene regulation, gene silencing, X-inactivation

Table 1. Overview of CpG distribution in the human genome

Subset	Length, Mb	GC content	Observed CpG fraction	Normalized CpG fraction
Whole genome	3.1*	0.38	0.009	0.25
1 kb upstream regions	15	0.53	0.042	0.60
1 kb downstream regions	15	0.45	0.013	0.26
Transcription units	930	0.42	0.011	0.26
Exons	45	0.50	0.028	0.45
Introns	880	0.41	0.010	0.24

Length refers to the total length of DNA examined.

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
Mapping of regulatory elements using ChIP-on-chip and/or ChIP-seq

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Genomic views of distant-acting enhancer:
“Transcribing the genome” Nature 435: 199-205 (10 September 2009)


ChIP-seq + ChIP-on-chip (Master of arrays and NGS)

Cell. 2010 Apr 30;141(3):432-45.

“c-Myc regulates transcriptional pause release.”

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


Research for “Genome and Transcribing the genome”..

Gene Expression + ChIP-on-chip (Time Course Exp.)

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

Promoter ChIP-on-chip Microarray Feature and Benefits

Feature	Benefit
• Higher Density	• Expanded coverage to detect more sites of significance
• Flexibility in Content and Format	• Option to balance coverage with cost-effectiveness and throughput
• Validated content	• Ability to detect DNA-protein binding sites with few probes and high confidence
• Updated content	• Comprehensive coverage based on current human Hg19 and mouse Mm9 databases

Probe coverage comparison

UCSC Genome Browser on Human Feb. 2009 (GRCh37/hg19) Assembly

Agilent Promoter ChIP-on-chip Microarrays Comparison

	Format	Design ID	# of Genes	Region Coverage	Resolution
Human	G3 1x1M	027811	~21,000	-9k to +2k	178bp (median)
	G3 2x400K	027954	~21,000	-3.5k to +1k	172bp (median)
	HD 1x244K	014706 014707	~17,000	-5.5k to +2.5K	~195bp
Mouse	G3 1x1M	028383	~19,000	-9k to +2k	180bp (median)
	G3 2x400K	028384	~19,000	-3.5k to +1k	193bp (median)
	HD 1x244K	014716 014717	~17,000	-5.5k to +2.5K	~195bp

SurePrint G3 Human Promoter ChIP-on-chip 1x1M Microarray Kit

- The same high sensitivity and accuracy provided by the earlier generation arrays but adds the benefits of extended coverage, updated content and increased cost savings. C+ version scanner is required.
- **Coverage includes extra miRNA and other small non-coding RNA regions.**

Agilent Part Number	G4873A
Design ID	27811
Sequences	~21,000 of the best-defined human genes represented as defined by RefSeq Probes 966,092 probes are represented, median spacing of 178bp.
Coverage	-9kb upstream to +2kb downstream of the transcriptional start sites, ~46 probes/gene
Composition	Enriched content annotated to UCSC hg19 (NCBI Build 37, February 2009)

SurePrint G3 Human Promoter ChIP-on-chip 2x400K Microarray Kit

- The same high sensitivity and accuracy provided by the earlier generation arrays but adds the benefits of extended coverage, updated content and increased cost savings. C+ version scanner is required.
- Each slide contains **two 400K** Microarrays

Agilent Part Number	G4874A
Design ID	27954
Sequences	~21,000 of the best-defined human genes represented as defined by RefSeq Probes 414,043 probes are represented, median spacing of 172bp.
Coverage:	-3.5kb upstream to +1kb downstream of the transcriptional start sites, ~20 probes/gene, <i>without</i> extra miRNA and other small non-coding RNA regions
Composition	Enriched content annotated to UCSC hg19 (NCBI Build 37, February 2009)

SurePrint HD Human Promoter ChIP-on-chip

1x244K Microarray, 2-design Set

- Specifically designed to identify human DNA binding proteins through traditional chromatin immunoprecipitation (ChIP) pull-down coupled with powerful DNA microarray location analysis.
- Earlier format (HD) optimized for B-version scanner

Agilent Part Number	G4495A
Design ID	014706, 014707
Sequences	~17,000 of the best-defined human genes represented as defined by RefSeq
Probes	-5.5kb upstream to +2.5kb downstream of the transcriptional start sites, ~25 probes/gene. Spacing ~195bp.
Composition	Enriched content sourced from UCSC hg18, NCBI build 36.1 (March 2006)

Human ENCODE Microarray

1x244K Microarray

- Designed to specifically identify human DNA binding proteins within the ENCODE regions.
- The Human ENCODE microarray is optimized using probes specifically validated for HD ChIP-on-chip technology.

Agilent Part Number	G4495A
Design ID	14792
Sequences	Specific ENCODE regions of chromosomes 1-22
Probes	~153 probes
Composition	ENCODE March 2006 human assemblies

SurePrint G3 Mouse Promoter ChIP-on-chip

Product Name	SurePrint G3 Mouse Promoter Microarray, 1x1M
Agilent Part Number	G4875A
Design ID	028383
Sequences	~19,000 of the best-defined mouse genes represented as defined by RefSeq
Probes	968,007 probes are represented, median spacing of 180bp.
Coverage:	-8kb upstream to +2kb downstream of the transcriptional start sites, ~49 probes/gene.

Product Name	SurePrint G3 Mouse Promoter Microarray, 2x400K
Agilent Part Number	G4876A
Design ID	28384
Sequences	~21,000 of the best-defined human genes represented as defined by RefSeq
Probes	415,814 probes are represented, median spacing of 193bp.
Coverage:	-3.5kb upstream to +1kb downstream of the transcriptional start sites, ~21.5 probes/gene.

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* Some other designs available through eArray Design Program

** Roughly calculated by coverage / probe numbers

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New


CH3

Updated pre-designed human CpG island-related microarray

Product Name	Human DNA Methylation Microarray
Agilent Part Number	G4495A
Design ID	023795
Format	244K: 1 x 244K
Arrays/Slide	244K: 1
Slide/Kit	1
Sequences	27,627 expanded CpG islands and 5081 UMR
Probes	237,227 biological probes median probe spacing of 97 bp. Probes annotated against UCSC HG19
Source	


Straussman et al., Nat Struct and Mol Biol. 2009 May; 564-71
Irizarry et al., Nat Genet. 2009 Feb; 178-86.

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Possibilities: These repressor elements may operate at a distance to silence adjacent promoters. They may also serve as promoters for antisense RNA transcripts.

Reference: Straussman et al., Nature SMB (2009)

Agilent Master Tiling Probe Database

To Utilize Only High Quality Probes

Genome
(3 billion base)

Tiling Methodology (Agilent Probe Selection)

- Tile **45~60-mers** across non-RM region, followed by selection
 - Uniqueness (homology)
 - Tm
 - Avoid secondary structure

eArray
design tool
for “custom”
array needs

High quality probes with high S/N


<100bp resolution)

Probe Sequence
Database

Any Genome. High Agilent Quality

Expanding portfolio of ChIP-on-chip and CpG Island Arrays

Organism	Description	Build	Database
Human	Feature-annotated, 23.8 million probes tiled at <100 bp	hg18	Agilent
Mouse	Feature-annotated, 27.4 million probes tiled at <100 bp	mm8	Agilent
Rat	Coordinate-annotated, 22.1 million probes tiled at <100 bp	rn4	Agilent
<i>C. elegans</i>	Coordinate-annotated, 718 thousand probes tiled at ~100 bp	Ce2	Agilent
Arabidopsis	Coordinate-annotated, 735 thousand probes tiled at ~100 bp	Ath1	Agilent
Drosophila	Coordinate-annotated, 447 thousand probes for whole genome with ~100bp spacing	Dm2	Agilent
Zebrafish	Coordinate-annotated, promoter probes only. Both expanded and proximal promoter areas represented	Zv4	Whitehead
<i>S. cerevisiae</i>	Feature-annotated. Whole genome with ~9 bp spacing	SacCer1	Agilent
<i>S. pombe</i>	Coordinate-annotated, whole genome	-	Agilent

Customized
Microarray Set

Access to Customization : It's eArray now!

Log-in


Select probe

- Probe database
 - Gene Expression
 - microRNA / CGH
 - ChIP / CH3
- Probe Design (GE)

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Current Available application :
Gene Expression, CGH, ChIP-CH3, microRNA and Target Enrichment

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

What Agilent Provides

Regulatory proteins bind to promoter DNA regions *in vivo*


c-Myc, H3K18ac,
H3K9ac etc. ...

Label enriched pool


- # Cross-link protein-DNA complexes
- # Lyse cells and sonicate DNA
- # IP chromatin to capture and purify bound DNA
- # Release and amplify DNA fragments

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

ChIP-on-chip実験フロー Workflow “Improvement” !

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Magna ChIP²™ kit (MILLIPORE) Chromatin IP DNA Microarray Universal kit

NEW

- The **first** ChIP kit specifically developed to couple with on-chip to **simplify** and **standardize** the genome-wide mapping studies.
- Designed to help address the challenges and make ChIP-on-chip technology easily accessible to a wider range of laboratories

The enhanced solution collaborated with MILLIPORE

Millipore Magna
ChIP² kit

Chromatin IP

LMPCR

Labeling

Hybridization
& Read

Agilent gDNA labeling kit
Agilent Hybridization kit
Agilent Microarray kit

Streamline Epigenetic Analysis

From ChIP kit to data analysis

- By providing optimized reagents, validated protocols, and expert support for the entire ChIP-on-chip workflow, Magna ChIP² kits help ensure **success, sensitivity** and **reproducibility** of ChIP-on-chip experiments using both Agilent and user-provided microarrays.

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
Agilent R&D data – CpG proof of concept

The HOX gene cluster in a lung cell carcinoma cell line characterized by higher CpG island promoter methylation


DNA methylation event detection

Not standardized yet..

Title, journal, date, last author	Detection algorithm	Normalization for input
“Developmental programming of CpG island methylation profiles in the human genome.”, Nat Struct Mol Biol. 2009 May, Cedar H	Tm adjusted Z-score	Global normalization in FE
“The presence of RNA polymerase II, active or stalled, predicts epigenetic fate of promoter CpG islands.”, Genome Res. 2009 Nov, Ushijima T	Me-value based on WI pXbar	No normalization performed
“Epigenetic profiling at mouse imprinted gene clusters reveals novel epigenetic and genetic features at differentially methylated regions.”, Genome Res. 2009 Aug, Beaudet AL	Predefined Peak shape detection (v2.0)	Unknown
“Human DNA methylomes at base resolution show widespread epigenomic differences”, Nature, 2009 Oct14, Ecker JR	MethylC-Seq	No consideration

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CGH module : Save aberration results in centralized data base

ChIP module : Algorithm “Pre-defined Peak Detection”

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*How do you begin to solve the challenges
of integrative multi-omics data analysis?*

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


GSE 11245 (Agilent Human Promoter ChIP-on-Chip Set 244K)


GSE 19789 (Affymetrix Human Human Gene 1.0 ST Array)

Diagnostic Application

Pharma and clinical:

- Focus on methylation patterns to characterize disease profiles/risk and create biomarkers (e.g. cancer)
- Understand role of methylation in X-inactivation diseases (e.g. Rett's) or imprinting disease (e.g. Prader-Willi, Lupus, Fragile X)
- The reversibility of epigenetic changes important for pharma manipulation

PNAS

Homeobox gene methylation in lung cancer studied by genome-wide analysis with a microarray-based methylated CpG island recovery assay

Tibor Rauch*, Zunde Wang*, Xinmin Zhang†, Xueyan Zhong*, Xiwei Wu‡, Sean K. Lau§, Kemp H. Kernstine¶, Arthur D. Riggs*, and Gerd P. Pfeifer*||

*Division of Biology, †Division of Information Sciences, §Division of Pathology, and ¶Division of Surgery, Beckman Research Institute of the City of Hope, Duarte, CA 91010; and *NimbleGen Systems, Inc., Madison, WI 53711

Contributed by Arthur D. Riggs, February 5, 2007 (sent for review December 11, 2006)

Homeobox genes are potentially useful as DNA methylation markers for early diagnosis of the disease.

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- Title: Chd1 regulates open chromatin and pluripotency of embryonic stem cells.
Journal: Nature. 2009 Aug 13;460(7257):863-8. Epub 2009 Jul 8.
Authors: Gaspar-Maia A, Alajem A, Polessko F, Sridharan R, Mason MJ, Heidersbach A, Ramalho-Santos J, McManus MT, Plath K, Meshorer E, Ramalho-Santos M.
Highly accessible chromatin is essential for the unique properties of stem cells
- Title: Lineage-specific DNA methylation in T cells correlates with histone methylation and enhancer activity.
Journal: Genome Res. 2009 Jul;19(7):1165-74. Epub 2009 Jun 3.
Authors: Schmidl C, Klug M, Boeld TJ, Andreesen R, Hoffmann P, Edinger M, Rehli M.
Histone and DNA methylation in cell type-specific genes in Tconv and Treg cells
- Title: DNA-binding specificity and in vivo targets of *Caenorhabditis elegans* nuclear factor I.
Journal: Proc Natl Acad Sci U S A. 2009 Jul 21;106(29):12049-54
Authors: Whittle CM, Lazakowitch E, Gronostajski RM, Lieb JD.
Unexpected small target set of NFI-1 for core cellular processes
- Title: The Level of the Transcription Factor Pax6 Is Essential for Controlling the Balance between Neural Stem Cell Self-Renewal and Neurogenesis.
Journal: PLoS Genet. 2009 Jun;5(6):e1000511
Authors: Sansom SN, Griffiths DS, Faedo A, Kleinjan DJ, Ruan Y, Smith J, van Heyningen V, Rubenstein JL, Livesey FJ
Core network regulating neural, neocortical stem cell
- Title: Biofilm Matrix Regulation by *Candida albicans* Zap1.
Journal: PLoS Biol. 2009 Jun;7(6):e1000133
Authors: Nobile CJ, Nett JE, Hernday AD, Homann OR, Deneault JS, Nantel A, Andes DR, Johnson AD, Mitchell AP.
Zinc-responsive regulatory protein controls matrix formation
- Title: Genomic distribution of CHD7 on chromatin tracks H3K4 methylation patterns.
Journal: Genome Res. 2009 Apr;19(4):590-601. Epub 2009 Feb 27
Authors: Schnetz MP, Bartels CF, Shastri K, Balasubramanian D, Zentner GE, Balaji R, Zhang X, Song L, Wang Z, Laframboise T, Crawford GE, Scacheri PC
Cell-specific binding of CHD7 correlates with H3K4me in tumor and ES cells

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Title: Epigenetic profiling at mouse imprinted gene clusters reveals novel epigenetic and genetic features at differentially methylated regions.

Journal: Genome Res. 2009 Aug;19(8):1374-83

Authors: Dindot SV, Person R, Strivens M, Garcia R, Beaudet AL

DNA and histone methylation at mouse imprinted geneclusters

Title: Methylation of polycomb target genes in intestinal cancer is mediated by inflammation.

Journal: Cancer Res. 2008 Dec 15;68(24):10280-9

Authors: Hahn MA, Hahn T, Lee DH, Esworthy RS, Kim BW, Riggs AD, Chu FF, Pfeifer GP

Methylation of polycombe target genes in intestical cancer

Title: Probe signal correction for differential methylation hybridization experiments.

Journal: BMC Bioinformatics. 2008 Oct 23;9:453

Authors: Potter DP, Yan P, Huang TH, Lin S

Data analysis of DNA methylation array

Title: Astrocyte-specific genes are generally demethylated in neural precursor cells prior to astrocytic differentiation.

Journal: PLoS ONE. 2008 Sep 11;3(9):e3189

Authors: Hatada I, Namihira M, Morita S, Kimura M, Horii T, Nakashima K

Demethylation in neural precursor cells before differentiation

Title: Epithelial Progeny of Estrogen-Exposed Breast Progenitor Cells Display a Cancer-like Methylome

Journal: Human Cancer Genetics Program, Department of Molecular Virology, Immunology and Medical Genetics, Comprehensive Cancer Center, The Ohio State University, Columbus, OH 43210, USA.

Authors: Cheng AS, Culhane AC, Chan MW, Venkataramu CR, Ehrich M, Nasir A, Rodriguez BA, Liu J, Yan PS, Quackenbush J, Nephew KP, Yeatman TJ, Huang TH

Estrogen-exposed breast progenitor cells show cancer-like methylome

Title: X-inactivation in female human embryonic stem cells is in a nonrandom pattern and prone to epigenetic alterations.

Journal: Proc Natl Acad Sci U S A. 2008 Mar 25;105(12):4709-14. Epub 2008 Mar 13

Authors: Shen Y, Matsuno Y, Fouse SD, Rao N, Root S, Xu R, Pellegrini M, Riggs AD, Fan G

X-inactivation in female embryonic stem cells

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ChIP-on-chip-related scientific publications (2010-2011)

Cell
PRESS

One more is coming in NAR. Stay tune!

Molecular Cell
Article

Letter
Re
G
of
Sh
How
of M

Genome-wide analysis of chromatin-associated protein kinase C-θ in human T lymphocytes
determines the role of microRNAs in gene expression

p300-
Function

Received for pub

Rahul K. Ve
and Devyan
From the [†]De

[‡]National Ce
Pashan, Pune 411021, India

Chromatin-Associated Protein Kinase C-θ Regulates an Inducible Gene Expression Program and MicroRNAs in Human T Lymphocytes

Elissa L. Sutcliffe,¹ Karen L. Bunting,^{3,8} Yi Qing He,^{1,8} Jasmine Li,¹ Chansavath Phetsouphanh,⁸ Nabila Seddiki,⁸ Anjum Zafar,¹ Elizabeth J. Hindmarsh,¹ Christopher R. Parish,¹ Anthony D. Kelleher,⁶ Russell L. McInnes,⁴ Toshiki Taya,⁵ Peter J. Milburn,² and Sudha Rao^{1,7,*}

¹Department of Immunology

²Biomolecular Resource Facility

John Curtin School of Medical Research, Australian National University, Canberra, ACT 2601, Australia

³Division of Hematology and Medical Oncology, Department of Medicine, Weill Cornell Medical College, Cornell University, New York, NY 10065, USA

⁴Agilent Technologies, 347 Burwood Highway, Forest Hill VIC 3131, Australia

⁵Agilent Technologies Japan, Hachioji Site, 9-1 Takakura-Cho, Hachioji-Shi, Tokyo 192-8510, Japan

⁶National Centre in HIV Epidemiology and Clinical Research, St. Vincent's Centre for Applied Medical Research, University of New South Wales, Sydney NSW 2010, Australia

⁷Discipline of Biomedical Sciences, Faculty of Applied Science, The University of Canberra, Canberra ACT 2601, Australia

IMMUNO

Integrated biochemical and computational approach identifies BCL6 direct target genes controlling multiple pathways in normal germinal center B cells

*Katia Basso,^{1,2} *Masumichi Saito,¹ *Pavel Sumazin,³ Adam A. Margolin,³ Kai Wang,³ Wei-Keat Lim,³ Yukiko Kitagawa,¹ Christof Schneider,¹ Mariano J. Alvarez,³ †Andrea Califano,³ and †Riccardo Dalla-Favera^{1,4,5}

¹Institute for Cancer Genetics and H. Irving Comprehensive Cancer Center, Columbia University, New York, NY; ²Department of Pediatrics, University of Padua, Padova, Italy; and ³Joint Centers for Systems Biology and Departments of ⁴Pathology and ⁵Genetics & Development, Columbia University, New York, NY

Transcriptional role of cyclin D1 in development revealed by a genetic-proteomic screen

Microarrays

- Human ChIP-on-chip promoter
- Mouse ChIP-on-chip promoter
- Drosophila ChIP-on-chip
- Yeast ChIP-on-chip

Agilent Technologies

論文紹介 Why It's Interesting

Agilentgenomics.jp/journalclub

... and many more to come

"Genes for embryo development are packaged in blocks of multivalent chromatin in zebrafish sperm." Genome Res. 2011

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
"Chromatin signatures of the Drosophila replication program Drosophila" Genome Res. 2010

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Methylation-related scientific publications (2010-2011)

Effects of genome arch
CpG islands to aberrant

Hideyuki Takeshima^a, Satoshi Y

Letter

Widespread and tissue-specific methylation changes

Shinji Maegawa,¹ George Hinkley,² Jixin Zha,³ and Jean-Philippe Giguere¹

RESEARCH

¹Department of Molecular and Cell Biology, The University of Texas at Austin, TX, USA

Genome-wide methylation analysis of sense-antisense pairs in

Yutaka Watanabe^a

Naoyuki Yamamoto^c, Kazufumi Watanabe^b, Hidemasa Kato^b, Kuniya Abe^{a,b}, Hidenori Kiyosawa^{a,b,c}

Cancer
Science

The official journal of the Japanese Cancer Association

JCA

Methylation silencing of angiopoietin-like 4 in rat and human mammary carcinomas

Naoko Hattori, Eriko Okochi-Takada, Mizuho Kikuyama, Mika Wakabayashi, Satoshi Yamashita and Toshikazu Ushijima¹

Division of Epigenomics, National Cancer Center Research Institute, Tokyo, Japan

OPEN ACCESS Freely available online

PLOS ONE

Broad Epigenetic Signature of Maternal Care in the Brain of Adult Rats

Patrick O. McGowan^{1,2,3*}, Matthew Suderman^{2,4,5*}, Aya Sasaki^{1,2,3*}, Tony C. T. Huang⁴, Michael Hallett⁵, Michael J. Meaney^{1,2,6,7}, Moshe Szyf^{2,4,7*}

¹Douglas Mental Health University Institute, Montreal, Quebec, Canada, ²Sackler Program for Epigenetics and Developmental Psychobiology at McGill University, McGill University, Montreal, Quebec, Canada, ³Centre for the Neurobiology of Stress, University of Toronto, Scarborough, Toronto, Ontario, Canada, ⁴Department of Pharmacology and Therapeutics, McGill University, Montreal, Québec, Canada, ⁵McGill Centre for Bioinformatics, McGill University, Montreal, Quebec, Canada,

⁶Singapore Institute for Clinical Sciences, Singapore, Republic of Singapore, ⁷Experience-Based Brain and Biological Development Program of the Canadian Institute for Advanced Research, Toronto, Ontario, Canada

A Genome-wide screen identifies frequently methylated genes in haematological and epithelial cancers

Thomas Dunwell¹, Luke Hesson¹, Tibor A Rauch², Lihui Wang³, Richard E Clark³, Ashraf Dalol¹, Dean Gentle¹, Daniel Catchpoole⁴, Eamonn R Maher¹, Gerd P Pfeifer², Farida Latif^{1*}, Naoyuki Yamamoto^c, Kazufumi Watanabe^b, Hidemasa Kato^b, Kuniya Abe^{a,b}, Hidenori Kiyosawa^{a,b,c}

Microarrays

- Human CpG Island
- Mouse CpG island
- Rat CpG island
- Custom Tiling

Agilent Technologies

論文紹介 Why It's Interesting

Agilentgenomics.jp/journalclub

... and many more to come

"Methylation silencing of angiopoietin-like 4 in rat and human mammary carcinomas" Can. Sci. 2011

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"Genome-wide analysis of expression modes and DNA methylation status at sense-antisense transcript loci in mouse." Genomics. 2010

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The Transcription Puzzle

Genomes don't just encode protein-coding RNAs

Different classes of ncRNAs identified to date

- Revealed novel sets of new small and large ncRNAs
- Key players for gene regulation, genome stability, and chromatin modification, etc.

The Biological Roles

- ncRNAs in translation
- ncRNAs in RNA splicing
- ncRNAs in gene regulation (*trans-* and *cis-*)
- ncRNAs and genome defense
- ncRNAs and chromosome structure
- Bifunctional RNA

ncRNA	Example
Short	t/r/sn/mi/piRNAs
Long	Xist/H19/Air/HOTAIR/...

Long noncoding RNAs in gene regulation


The silencing effect

Cell (2007)

“Functional Demarcation
of Active and Silent
Chromatin Domains in
Human HOX Loci by
Noncoding RNAs”

Cell (2010)

“A Large Intergenic
Noncoding RNA Induced
by p53 Mediates Global
Gene Repression in the
p53 Response”

Interdisciplinary Reviews, Focus Article, Published Online: Mar, 2011


“New” Long noncoding RNAs in gene regulation

The boosting effect

LincRNA: Large intervening (or Long intergenic) non-coding RNA

Cell (2010)

Ørom UA et.al., “Long Noncoding RNAs with Enhancer-like Function in Human Cells”


Nature (April, 2011)

Wang KC et.al., “A long noncoding RNA maintains active chromatin to coordinate homeotic gene expression”

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HOTTIP, ‘HOXA transcript at the distal tip’

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Enables **robust amplification for genome-wide analysis of DNA methylation** with PfuTurbo Cx Hotstart DNA Polymerase, a proofreading polymerase that not only knows what to do with deoxyuracil, but also has high proofreading activity !

SureCycler is “All-in-one” Thermal Cycler
It’s Turbo-Charged !

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

MULTIPLE TECHNIQUES FOR GENE EXPRESSION STUDIES

Best in Class just got better.....

Agilent Genomics Workflow Solution

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