

# Highly efficient multiplex PCR using novel reaction chemistry

## Application

Holger Engel  
Corinna Küppers  
Dirk Löffert

### **Abstract**

A new multiplex PCR kit has been introduced that allows amplification of multiple gene targets without the need to optimize PCR conditions. In this Application Note, the kit was used to amplify up to 19 targets in a single PCR step under standard conditions. Amount and purity of the PCR products has been confirmed using the Agilent 2100 bioanalyzer.

## Introduction

Multiplex PCR is a powerful technique enabling amplification of two or more products in a single reaction. Typically, it is used for genotyping applications where simultaneous analysis of multiple markers is required, such as typing of normal and genetically modified animals and plants, detection of pathogens or genetically modified organisms (GMOs), or for microsatellite analyses. Multiplex assays can be tedious and time-consuming to establish, requiring lengthy optimization procedures, such as adjusting primer concentrations,  $Mg^{2+}$  concentration, and the amount of enzyme. In many cases, the results are still disappointing and further extensive optimization may be required. However, due to novel developments in the reaction chemistry, multiplex assays are now simple and straightforward to establish. The combined use of a highly stringent hot start with a unique PCR buffer specially developed for multiplex reactions, makes the QIAGEN® Multiplex PCR Kit highly suited for multiplex PCR applications. The newly developed reaction buffer, containing a special multiplex PCR-enhancing synthetic factor, eliminates the need for optimization — even when using equimolar primer concentrations. In this study, we describe the effects of different reaction parameters and how the QIAGEN Multiplex PCR Kit minimizes the need for optimization of multiplex PCR assays. The results were visualized using the Agilent 2100 bioanalyzer.


## Material and methods

Genomic DNA was isolated from human K562 cells using the DNeasy® Tissue Kit and 20 ng were used as template for multiplex PCR. PCR was performed using either QIAGEN Multiplex PCR Master Mix (containing HotStarTaq® DNA Polymerase), or using standard Taq DNA polymerase in an otherwise identical reaction mixture, or using a Taq DNA polymerase with antibody-mediated hot start and KCl reaction buffer from Supplier I. Equal volumes of each multiplex PCR were analyzed on the Agilent 2100 bioanalyzer using the DNA 1000 LabChip® kit according to the supplied instructions.

## Results and discussion

### Effect of hot start on multiplex PCR specificity

Typically, a primer concentration of 0.2–0.5  $\mu M$  is used in conventional PCR. In contrast, the total primer concentration in multiplex PCR can be as high as 2–4  $\mu M$ , depending on the number of different primer pairs in the reaction. The large number of primers often results in the generation of non-specific PCR products and primer-dimers, reducing the specificity and sensitivity of the multiplex PCR. Using a stringent hot start to increase PCR specificity can prevent the generation of these non-specific products. Figure 1 shows a comparison of hot-start multi-




**Figure 1**  
**Effect of Hot Start on Multiplex PCR.** Multiplex PCR of 13 targets (99–955 bp) was carried out for 35 cycles in a 50  $\mu l$  reaction volume using the standard protocol in the QIAGEN Multiplex PCR handbook. Multiplex reactions were carried out using QIAGEN Multiplex PCR Master Mix (containing HotStarTaq DNA Polymerase) (QIAGEN Multiplex PCR Kit) or standard Taq DNA Polymerase (Taq). Equal volumes of each multiplex PCR were analyzed on the Agilent 2100 bioanalyzer using the DNA 1000 LabChip kit. A. Gel-like image of collected data. B. Electrophoretic trace of multiplex PCR.


plex PCR, using the QIAGEN Multiplex PCR Kit, and multiplex PCR, using a standard Taq DNA polymerase without a hot start. Analysis was carried out on the Agilent 2100 bioanalyzer using the DNA 1000 LabChip kit. Use of standard Taq DNA polymerase results in the generation of large amounts of primer-dimers and low efficiency of amplification. In contrast, the stringent hot start provided by HotStarTaq DNA Polymerase in the QIAGEN Multiplex PCR Kit leads to reproducible and efficient amplification of all 13 PCR products without generating nonspecific artifacts.

### Effect of novel multiplex PCR buffer chemistry on primer annealing

With multiplex PCR, the annealing efficiencies of the different primers in the reaction are usually dissimilar. Although it is possible to design primers that have a similar annealing temperature, the  $T_m$  of each primer does not always provide a good indication of its annealing efficiency. In addition, other factors such as the 3'-end sequence of primers may affect the efficiency of primer extension by Taq DNA polymerase. Efficient primer annealing and extension, irrespective of primer sequence, is achieved by a special synthetic factor in QIAGEN Multiplex PCR Master Mix. This synthetic factor increases the local concentration of primers at the template DNA and stabilizes specifically bound primers (figure 2). In addition, nonspecific primer binding is avoided by a balanced combination of salts and the specially optimized  $NH_4^+$  concentration in the QIAGEN Multiplex PCR Master Mix.  $NH_4^+$ , which exists



**Figure 2**  
Effect of Novel Synthetic Factor and  $NH_4^+$  Ions on Primer Annealing.  $Mg^{2+}$  and other salts bind to phosphate groups (P) on the DNA backbone, which stabilizes the annealing of the primers to the template.  $NH_4^+$ , which exists both as the ammonium ion and as ammonia under thermal-cycling conditions, can interact with the hydrogen bonds between the bases, destabilizing principally the weak hydrogen bonds at mismatched bases (B) of non-specifically bound primers. The synthetic factor (S) pushes the primers toward the template and stabilizes specifically bound primers, enabling efficient extension of all primers in the reaction.



**Figure 3**  
Effect of PCR Buffer on Multiplex PCR. Multiplex PCR of 19 targets (99–955 bp) was carried out for 35 cycles using standard conditions (Std) for the QIAGEN Multiplex PCR kit without optimization, or using the indicated  $Mg^{2+}$  concentrations with a hot-start enzyme and supplied KCl-based buffer from Supplier I. Equal volumes of each multiplex PCR were analyzed on the Agilent 2100 bioanalyzer using the DNA 7500 LabChip kit. A. Gel-like image of collected data. B. Electrophoretic trace of multiplex PCR.

predominantly as ammonia (NH<sub>3</sub>) under thermal-cycling conditions, interacts with the relatively weak hydrogen bonds formed when primers bind nonspecifically to the template DNA and destabilizes these nonspecifically bound primers. Figure 3 shows that in contrast to conventional PCR reagents, the formulation of the new QIAGEN Multiplex PCR Master Mix, ensures comparable efficiencies for annealing and extension of all primers in the reaction without further optimization.

## **Conclusion**

Establishing multiplex PCR assays is easy and fast using the QIAGEN Multiplex PCR Kit. Tedious optimization procedures are virtually eliminated. These properties make the kit well suited for multiplex applications, including genotyping of transgenic organisms, detection of pathogens or GMOs, and microsatellite genotyping (e.g., short tandem repeat (STR) and variable number tandem repeat (VNTR) analyses). The stringent hot start provided by HotStarTaq DNA Polymerase eliminates non-specific PCR products and primer-dimer formation in multiplex PCR, while a novel synthetic factor in the multiplex PCR buffer allows efficient annealing of multiple primers under identical cycling conditions. Visualization of multiplex PCR products can be performed efficiently using the Agilent 2100 bioanalyzer and DNA LabChip kits, because of the accurate sizing and quantitation that can be achieved.

*Holger Engel, Corinna Küppers and Dirk Löffert are scientists at the Technology Center Modification/Amplification Qiagen GmbH, Hilden, Germany*

**[www.agilent.com/chem/labonachip](http://www.agilent.com/chem/labonachip)  
[www.qiagen.com](http://www.qiagen.com)**



LabChip<sup>®</sup> and the LabChip Logo are registered trademarks of Caliper Technologies Corp. in the U.S. and other countries.

Patented or patent-pending technology and/or registered or registration-pending trademarks of the QIAGEN Group: QIAGEN<sup>®</sup>, DNeasy<sup>®</sup>, HotStarTaq<sup>®</sup>.

Purchase of QIAGEN products for PCR is accompanied by a limited license to use them in the Polymerase Chain Reaction (PCR) process for research and development activities in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Applied Biosystems or as purchased, i.e. an authorized thermal cycler. The PCR process is covered by U.S. Patents 4,683,195 and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG. Patents of third parties in certain countries may cover the process of multiplex PCR or of certain applications.

The information in this publication is subject to change without notice.

Copyright © 2003 Agilent Technologies  
All Rights Reserved. Reproduction, adaptation or translation without prior written permission is prohibited, except as allowed under the copyright laws.

Published April 1, 2003  
Publication Number 5988-9342EN



**Agilent Technologies**