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The screenshot illustrates the process of checking MSD logs. The main interface shows a schematic of the MSD components. A context menu is open over the schematic, with 'Report MS Faults' selected. A 'MSD Diagnosis' dialog box shows 'No Ms Fault is detected.' Below, two windows are shown: 'Maintenance Logbook Entry' with a list of maintenance messages and 'Ms Fault Report' showing a table of fault reports.

**Maintenance Logbook Entry**

Instrument: G6130B Serial No.: 5G12015101

Maintenance Message: [Dropdown menu]

Operator: Admin

- API-ES Nebulizer Needle adjusted
- API-ES Nebulizer Needle replaced
- APCI Nebulizer Needle adjusted
- APCI Nebulizer Needle replaced
- Corona Needle cleaned
- Corona Needle replaced
- Capillary cleaned
- Capillary replaced
- EM Horn replaced
- Rough Pump Fluid Level checked
- Rough Pump Fluid replaced
- Gas Conditioner replaced
- Calibrant Levels checked
- Calibrant Bottle (A) replaced
- Calibrant Bottle (B) replaced
- CDS flushed
- MS Selection Valve replaced

**Ms Fault Report**

Tune: reported Ms Fault during tune (no entry found)  
Logbk.: last reported Ms Fault in error logbook (no entry found)  
Actual: actual Ms Fault state (6/21/2017 3:54:53 PM)

Description	Tune	Logbk.	Actual
5. Difficulty with calibrant delivery system. Calibration delivery system leak occurred.			X