

A Label-Free Screening Approach to Epigenetic Drug Discovery: Using Mass Spectrometry to Monitor Protein and Nucleic Acid Modification Events

William A. LaMarr, Peter T. Rye, Lauren E. Frick, Agilent Technologies, Inc. Wakefield, MA 01880

SLAS 2012
Poster #MP116



Introduction

Epigenetics, the study of changes in gene expression caused by mechanisms other than deviations in DNA sequence, is rapidly emerging as a field with tremendous drug development potential across a broad range of therapeutic areas. Epigenetic discovery via traditional screening technologies is complicated by the ability of these enzymes to impart multiple modification events onto singular substrates resulting in complicated and dynamic enzyme kinetic environments. In this study, we investigated the use of mass spectrometry as an almost universal platform for the detection and quantification of multiple protein and nucleic acid modification events (i.e. histone acetylation status, histone methylation status, DNA methylation status). The use of mass spectrometry facilitates direct detection of native enzyme substrate and product pairs, and allows for rapid assay development of a wide array of chemical analytes (peptides, oligonucleotides, whole proteins, etc...). The label-free methodology mass spectrometry provides an enable targeted screening across a large collection of epigenetic processes and represents a ubiquitous approach to this complex research area.

Experimental

All assays were investigated via solid phase extraction electrospray mass spectrometry using the Agilent RapidFire High-throughput MS System (RapidFire/MS) which facilitates cycle times of < 10 seconds per sample.

Assay Specific Conditions




Analyte Class	RF Cartridge	Wash Buffer	Elute Buffer	MS System
Peptides	C4	50% MeCN 0.09% formic acid 0.01% TFA	100% MeCN 0.09% formic acid 0.01% TFA	QQQ (6460)
Co-factors	Graphitic Carbon	dH ₂ O	25% MeCN 25% isoctone 5 mM NH ₄ acetate	QQQ (6460)
Whole Proteins	C4	dH ₂ O 0.09% formic acid 0.01% TFA	100% MeCN 0.09% formic acid 0.01% TFA	Q-TOF (6530)
DNA Oligos	Cyano	dH ₂ O 6 mM octylammonium acetate	80% MeCN 6 mM octylammonium acetate	QQQ (6460)

Mass spectrometry parameters were independently developed for each specific analyte of interest on the appropriate MS platform.

Histone Protein Acetylation Status



SIRT-1 Protein Deacetylation

The sirtuin enzymes, a class of NAD⁺-dependent histone deacetylases, are a focal point of epigenetic research because of their roles in regulating gene expression and cellular differentiation by deacetylating histones and a host of transcription factors, including p53. SPE-MS/MS technology offers two independent label-free screening methodologies to study sirtuin activity: 1) peptide specific quantitation and 2) 0-acetyl-ADP-ribose co-product quantitation.



p300 Protein Acetylation



Histone acetyltransferases (HATs) catalyze the transfer of an acetyl group from an acetyl-coenzyme A donor molecule to specific lysine residues within proteins. The acetylation state of proteins, particularly histones, is known to modulate their intermolecular binding properties and control various cellular processes, most notably transcriptional activation. In addition, deregulation of HAT activity has been linked to the development of a number of cancers; therefore, compounds that affect these enzymes have strong potential as therapeutic agents. SPE/MS/MS technology offers three distinct label-free HAT screening approaches: 1) AcCoA/CoA pair of cofactors, 2) direct peptide measurement and 3) whole protein measurement.



Histone Protein Methylation Status


LSD-1 Histone Demethylation

Histone methylation status can be a major determinant in gene expression by altering histone protein structure thereby regulating the access to specific DNA sequences. The ability to regulate the transcription rate of a multitude of genes has been linked to a variety of disease states, making enzymes that control histone methylation attractive drug discovery targets. SPE/MS/MS technology presents the unique ability to uniquely identify and quantitate independent methylation events on single lysine residues within histone proteins.



Nucleic Acid Methylation Status

Over the past couple years, the application of the RapidFire system to studying epigenetic modifications has grown significantly. Primarily, however, the target substrates have been peptide-based. To help establish the applicability of the platform for investigating epigenetic changes to DNA, preliminary experiments using DNA standards were conducted.



Conclusions

We have demonstrated the utility of solid phase extraction electrospray mass spectrometry in multiple epigenetic screening applications. In addition to methylation and acetylation changes, histone tails can be modified by a diverse range of mechanisms including phosphorylation, deimination, ubiquitination, sumoylation, ADP-ribosylation, etc...; all of which may be monitored by mass spectrometry. With cycle times of < 10 seconds per sample, the RapidFire High-throughput MS System could be used to establish a range of peptide-based screening assays in primary epigenetic discovery programs.