

**Comparison of
metabolomic profiles
of androgen-
dependent and
refractory prostate
cancer cells using
liquid
chromatography-mass
spectrometry (LC-MS)**

ASMS 2010

Nagireddy Putluri¹, Shaiju Vareed¹, Theodore Sana², Steven Fischer², Gagan Thangjam¹, Vasanta Putluri¹, Ali Shojaie⁵, Judith Giri³, George Michailidis⁵, Martha K Terris⁴, and Arun Sreekumar¹

¹*Cancer Center, Medical College of Georgia, Augusta, GA;* ²*Metabolomics Laboratory Application Group, Agilent Technologies, Santa Clara, CA;* ³*Department of Biostatistics,, Medical College of Georgia, Augusta, GA;*

³*Pathology, Biomedical and Radiological Technologies, Medical College of Georgia, Augusta, GA,* ⁴*Department of Urology, Charlie Norwood Veteran Affairs Medical Center, Augusta, GA,* ⁵ *Department of Statistics, University of Michigan, Ann Arbor, MI*



Comparison of metabolomic profiles of androgen-dependent and refractory prostate cancer cells using liquid chromatography-mass spectrometry (LC-MS)

Introduction

Prostate cancer (PC) is the second most prevalent cancer among American men which is primarily treated by androgen ablation therapy. Although a number of patients respond to this regimen, a significant subset fail and the tumor invariably progresses into a hormone refractory metastatic state, which is lethal. Notably, clinicians lack tools that would allow them to select patients who respond to androgen ablation. Earlier we had reported the first unbiased metabolomic signature for localized prostate cancer. Advancing further, we attempt to delineate the subset of metabolome in prostate cancer which is regulated by androgen-action. In this study we examine global metabolomic alteration regulated by androgen in prostate derived cell lines using a validated unbiased LC-MS platform. In addition to revealing potential prognostic markers, these metabolomic profiles highlight the biochemical machinery altered in prostate cancer by androgen action.

Experimental

Instrumentation: An Agilent 1200 SL HPLC system with binary pump and degasser, well plate auto-sampler with thermostat controlled column compartment, and an Agilent 6520 Accurate-Mass QTOF mass spectrometer with dual ESI source operated in the positive and negative ionization modes. Dynamic mass axis calibration was achieved by continuous infusion of a reference mass solution using an isocratic pump connected to a dual sprayer feeding into an electrospray ionization source.

HPLC Method: Flow rate 0.6 ml/min, column temperature 60°C. Solvent A is Water + 0.2% Acetic Acid, solvent B is Methanol + 0.2% Acetic Acid. A 13 minute gradient from 2% B to 98% B and a 6 minute hold at 98% B with a stop time at 19 minutes, followed by 5 minutes post-run column re-equilibration. An Agilent Zorbax SB-Aq (2.1 X 50mm, 1.8 μ m) analytical column and an Agilent Zorbax-SB-C8 Rapid Resolution Cartridge (2.1X 30mm, 3.5 μ m) guard column were used for separation of the metabolome.

Experimental

Global Metabolomic Profiling of Prostate Cancer

Figure 1: Experimental overview





Figure 2: Reproducibility of profiling platform for spiked standards (n=12 compounds) across 7 independent replicates A) positive ionization mode B) negative ionization mode


LC-MS was used to profile metabolites in prostate derived cell lines

- 22Rv1, LNCaP and VCaP: androgen responsive (prostate cancer),
- DU145 and PC3: androgen non-responsive (prostate cancer)
- RWPE (benign).

Using reverse phase positive and negative mode ionizations 6747 and 2288 metabolites respectively, were detected across the prostate-derived cell lines

Results and Discussion

Figure 3: Heat map of compounds detected in at least 25% of cell lines examined (n= 17 including replicates)



Results and Discussion




Figure 6: Distribution of altered metabolites across the three sample classes

204 metabolites were found only in the prostate cancer cell lines, of which 52 were specific to androgen responsive cells and 114 were detected solely in androgen refractory cells. 141 metabolites were seen in benign cells only.

Figure 8: Heat map representation of compounds significantly altered across the three classes A) positive mode (n=311) and B) negative mode (n=171). In total 482 were significantly altered between benign, androgen-responsive and refractory cells.





Figure 7: Principal component analysis using compounds significantly altered across the three classes A) 311 positive mode (n=311) and B) negative mode (n=171)

List of named compounds significantly ($p < 0.005$, $FC > 2$) altered in prostate cancer by androgen action

1. Kynurinine
2. Guanine
3. Cartisone
4. Carnitine
5. Thymidine
6. Hippuric acid
7. DL-2-Indole-3-Lacticacid
8. DL-3 Phenyllactic acid
9. 3-Hydroxy-3-Methylvalerate



Future Work:

Validate androgen regulated metabolome in prostate cancer tissue specimens and delineate their mechanistic role in prostate cancer progression. The latter will be examined using Oncomine Concept Maps (OCM) and their associated genes using the Kyoto Encyclopedia of Genes and Genomes (KEGG, www.kegg.com).

Ref: Sreekumar, A et al. Metabolomic profiles delineate potential role for sarcosine in prostate cancer progression. Nature. 2009 Feb 12;457(7231):910-4.

Funding: 1R01CA133458-01 (AS), 1 R03 CA139489-01 (AS) and RCA145444A (AS and GM), Georgia Cancer Coalition (AS)