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Compound parameters

- ☒ Fragmentor (V)
 - From: 50 To: 200 Step Size: 10
 - ☐ Optimize from defined range
 - ☐ Optimize from method setpoint Step size: 10
- ☒ Collision Energy
 - From: 0 To: 80 Step Size: 4
 - ☐ Find product ions from defined range
 - ☐ Optimize from defined range
 - ☐ Optimize from method setpoint Step size: 4

Source parameters

- ☒ Determine parameter range from method source setpoints
 - ☒ Gas Temperature (°C)
 - ☒ Gas Flow (L/min)
 - ☒ Sheath Gas Temperature (°C)
 - ☒ Sheath Gas Flow (L/min)
 - ☒ Capillary Voltage (V)
 - ☐ Nebulizer (psi)
 - ☒ Nozzle Voltage (V)

What type of optimization do you wish to perform?

Select one of the optimization modes

Guided Method Optimization:

The user is guided through a series of steps to optimize a method. A systematic pause for user input and data review is required between each major phase of the optimization.

- Creates the most suitable MRM transitions from a list of chemical formulas or precursor ions
- Optimizes and fine-tunes MRM specific parameters
- Optimizes and fine-tunes ion source and front-end ion optics parameters
- Optionally updates compound and method database

Fully Automated Method Optimization:

Optimization is carried out in a fully automated manner. All user input for each phase is defined prior to the start of the optimization. Data review is only available once the optimization has completed.

- Creates the most suitable MRM transitions from a list of chemical formulas or precursor ions
- Optimizes and fine-tunes MRM specific parameters
- Optimizes and fine-tunes ion source and front-end ion optics parameters
- Optionally updates compound and method database

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	Status	Method	Data File	Sample Type	iReflex Type
1	Completed	ESDemo_MRM method.m	Blank_1.d	Blank	Carryover
2	Completed iReflex	ESDemo_MRM method.m	Blank_1-CarryoverBlank-001.d	Blank	Carryover
3	Completed	ESDemo_MRM method.m	Sample_1.d	Sample	No iReflex Workflow
4	Completed	ESDemo_MRM method.m	Sample_2.d	Sample	No iReflex Workflow
5	Completed	ESDemo_MRM method.m	Blank_2.d	Blank	Carryover
6	Completed iReflex	ESDemo_MRM method.m	Blank_2-CarryoverBlank-001.d	Blank	Carryover
7	Completed iReflex	ESDemo_MRM method.m	Blank_2-CarryoverBlank-002.d	Blank	Carryover
8	Completed	ESDemo_MRM method.m	Sample_3.d	Sample	No iReflex Workflow

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	Status	Method	Data File	Sample Type	Inj Vol (μl)	iReflex Type	iReflex Acq. Method	Red. Inj. Vol.
1	Completed	Above cal range.m	Blank_before run.d	Blank	As Method	No iReflex Workflow	Above cal range.m	
2	Completed	Above cal range.m	QC_before run.d	QC	As Method			
3	Completed	Above cal range.m	sample outside cal range_5.d	Sample	3	Above Cal. Range Append	Above cal range.m	0.5
4	Completed	Above cal range.m	sample inside cal range_5.d	Sample	1	Above Cal. Range Append	Above cal range.m	0.5
5	Completed iReflex	Above cal range.m	sample outside cal range_5-Fi-Blank.d	Blank	As Method	No iReflex Workflow		
6	Completed iReflex	Above cal range.m	sample outside cal range_5-iReflex.d	Sample	0.5	No iReflex Workflow		削減されたサンプル量による再注入

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	Status	Method	Data File	Sample Type	iReflex Type	iReflex Acq. Method	Red. Inj. Vol.	First Tier Quant Method	Second Tier Quant Method
1	Completed	FS1.m	BlankData_B1.d	Blank	No iReflex Workflow	FS3.m		Quant FS method 1.m	Quant FS method tier2.m
2	Completed	FS1.m	QCData_Q1.d	QC					
3	Completed	FS1.m	Sampledata_S1.d	Sample	Fast Screening Insert	FS3.m		Quant FS method 1.m	Quant FS method tier2.m
4	Completed	FS1.m	Sampledata_S2.d	Sample	Fast Screening Insert	FS3.m		Quant FS method 1.m	Quant FS method tier2.m
5	Completed iReflex	FS3.m	Sampledata_S2-iReflex-Blank.d	Blank	No iReflex Workflow				D:\Projects\Secondary injection\Methods\Quant FS method tier2.m
6	Completed iReflex	FS3.m	Sampledata_S2-iReflex.d	Sample	No iReflex Workflow				D:\Projects\Secondary injection\Methods\Quant FS method tier2.m
7	Completed iReflex	FS1.m	Sampledata_S2-Fi-Blank.d	Blank	No iReflex Workflow				D:\Projects\Secondary injection\Methods\Quant FS method tier2.m
8	Completed	FS1.m	Sampledata_S3.d	Sample	Fast Screening Insert	FS3.m		Quant FS method 1.m	Quant FS method tier2.m
9	Completed	FS1.m	BlankData_B2.d	Blank	No iReflex Workflow	FS3.m		Quant FS method 1.m	Quant FS method tier2.m

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MS Checktune Report - G6475A			
Instrument Information			
Model	G6475A	Checktune Date	2022-08-25T16:08:07-07:00
Serial Number	SQ22229001	SWFW Version	3.0.14248.134
Ion Source	AJS ESI	Ionization Mode	ESI
Last Autotune Date	2022-08-29T11:36:41-07:00	Last Tuned By	SYSTEM (SYSTEM)
Vacuum Pressure	2.25E+0 [R] (Torr); 2.90E-5 [H] (Torr)	Overall Result	Passed
Positive Ion Mode			
MS1 Peak Width Unit, Scan Speed Normal	Result	Passed	
MS2 Peak Width Unit, Scan Speed Normal	Result	Passed	
MS1 Peak Width Narrow, Scan Speed Normal	Result	Passed	
MS2 Peak Width Narrow, Scan Speed Normal	Result	Passed	
MS1 Peak Width Wide, Scan Speed Normal	Result	Passed	
MS2 Peak Width Wide, Scan Speed Normal	Result	Passed	
MS1 Peak Width Widest, Scan Speed Normal	Result	Passed	
MS2 Peak Width Widest, Scan Speed Normal	Result	Passed	
MS2 Scan Speed Fast	Result	Passed	
MS2 Scan Speed Ultra	Result	Passed	
MS1 Lag Factor	Result	Passed	
MS2 Lag Factor	Result	Passed	
Gain	Result	Passed	
Negative Ion Mode			
MS1 Peak Width Unit, Scan Speed Normal	Result	Passed	
MS2 Peak Width Unit, Scan Speed Normal	Result	Passed	
MS1 Peak Width Narrow, Scan Speed Normal	Result	Passed	
MS2 Peak Width Narrow, Scan Speed Normal	Result	Passed	
MS1 Peak Width Wide, Scan Speed Normal	Result	Passed	
MS2 Peak Width Wide, Scan Speed Normal	Result	Passed	
MS1 Peak Width Widest, Scan Speed Normal	Result	Passed	
MS2 Peak Width Widest, Scan Speed Normal	Result	Passed	
MS2 Scan Speed Fast	Result	Passed	
MS2 Scan Speed Ultra	Result	Passed	
MS1 Lag Factor	Result	Passed	
MS2 Lag Factor	Result	Passed	
Gain	Result	Passed	

Maintenance		
Early Maintenance Feedback Counters		
Autotune	Alert threshold (Days)	Expires on
☒ Enable	31	04-Jun-2022
	Set to factory value	Set to custom value
Checktune	Alert threshold (Days)	Expires on
☒ Enable	7	13-May-2022
	Set to factory value	Set to custom value
Injection count	Alert threshold (Counts)	Remaining (Counts)
☒ Enable	10000	972
Reset counter	Set to factory value	Set to custom value
Diverter valve switches	Alert threshold (Counts)	Remaining (Counts)
☒ Enable	10000	9894
Reset counter	Set to factory value	Set to custom value
Rough pump	Alert threshold (Days)	Expires on
☒ Enable	365	13-Apr-2023
Reset timer	Set to factory value	Set to custom value
Gas filter	Alert threshold (Days)	Expires on
☒ Enable	365	13-Apr-2023
Reset timer	Set to factory value	Set to custom value
Detector health	Detector lifetime remaining (%) 94.8%	
Nebulizer status		
Capillary status		
Spray stability status		

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1. Robustness of the Agilent Ultivo Triple Quadrupole LC/MS with the ESI Source, Agilent Technologies, Inc. 2019, 5994-0671EN.
2. Metzger, Benedikt; Willmann, Lucas; Simplify your Method Development using the Agilent 1260 Infinity Prime LC System and InfinityLab LC/MSD iQ, Agilent Technologies, Inc. 2020, 5994-2124EN.
3. Sartain, Mark; Sosienski, Theresa; Yang, Dan-Hui Dorothy; Robustness of the Agilent Ultivo Triple Quadrupole LC/MS for Routine Analysis in Food Safety, Agilent Technologies, Inc. 2017, 5991-8741EN.

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© Agilent Technologies, Inc. 2022