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
1. H. Stephen Stoker, "General, Organic and Biological Chemistry." 5th Edition Publisher: Boston, MA: Houghton Mifflin Co., 2006.
2. Anna Glinko, Michael J. Bozym, Michelle L. Owens, Karyn Usher, and Ronald E. Majors, 「Agilent ZORBAX Eclipse Plus カラムを用いた水溶性ビタミンの逆相 HPLC 分離」アジレント・テクノロジー、資料番号 5989-9313JAJP, August 6, 2008.
3. William Long and Anne E. Mack, "Fast Analysis of Environmental Phenols with Poroshell 120 EC-C18 Columns," Agilent Technologies publication 5990-6156EN, August 31, 2010.
4. Agilent Technical Note, "Transfer of Methods between Poroshell 120 EC-C18 and ZORBAX Eclipse Plus C18 Columns," Agilent Technologies publication 5990-6588EN, February 9, 2011.
5. Anne E. Mack and William J. Long "Fast, Low Pressure analysis of food and beverage additives using a Superficially Porous Agilent Poroshell 120 EC-C18 Column." Agilent Technologies publication 5990-6082EN, July 8, 2010.
6. William J. Long and Anne E. Mack, "Reduce Tubing Volume to Optimize Column Performance," Agilent Technologies publication 5990-4964EN, November 16, 2009.
7. Udo Huber, Analysis of Water-Soluble Vitamins by HPLC, 12/98, Agilent Technologies publication 5968-2971EN.
8. Long, W. and Mack, A. "Fast Analysis of Sulfa Drugs using the Agilent 1100 Series LC with Agilent Poroshell 120 EC-C18 columns," Agilent Technologies publication 5990-5572EN, 2010.
9. Detlef Wilhelm, "Determination of water soluble vitamins with the Agilent 1120 Compact LC after method development with the Agilent 1200 Series Rapid Resolution LC system and back transfer," Agilent Technologies publication 5990-4379EN, June 15, 2010.
10. Dorothy J. Phillips, Mark Capparella, Uwe D. Neue, Zoubair El Fallah, A new small particle packing for faster analysis with high resolution Journal of Pharmaceutical and Biomedical Analysis 15 (1997) 1389-1395.

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© Agilent Technologies, Inc., 2011