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
1. Steger-Hartmann, T. et al. (2002). Investigations into the environmental fate and effects of iopromide (ultravist), a widely used iodinated x-ray contrast medium. *Water Research*, 36(1), 266–274.
2. Ternes, T. A. & Hirsch, R. (2000). Occurrence and behavior of x-ray contrast media in sewage facilities and the aquatic environment. *Environmental Science & Technology*, 34(13), 2741–2748.
3. Carballa, M., et al. (2004). Behavior of pharmaceuticals, cosmetics and hormones in a sewage treatment plant. *Water Research*, 38(12), 2918–2926.
4. Kormos, J. L., Schulz, M. & Ternes, T. A. (2011). Occurrence of iodinated x-ray contrast media and their biotransformation products in the urban water cycle. *Environmental Science & Technology*, 45(20), 8723–8732.
5. Drewes, J. E., Fox, P. & Jekel, M. (2001). Occurrence of iodinated x-ray contrast media in domestic effluents and their fate during indirect potable reuse. *Journal of Environmental Science and Health Part A — Toxic/Hazardous Substances & Environmental Engineering*, 36(9), 1633–1645.
6. Drewes, J. E., et al. (2003). Fate of pharmaceuticals during ground water recharge. *Ground Water Monitoring and Remediation*, 23(3), 64–72.
7. Snyder, S. A., et al. (2007). Role of membranes and activated carbon in the removal of endocrine disruptors and pharmaceuticals. *Desalination*, 202(1–3), 156–181.
8. Dickenson, E. R. V. et al. (2011). Indicator compounds for assessment of wastewater effluent contributions to flow and water quality. *Water Research*, 45(3), 1199–1212.
9. Duirk, S. E., et al. (2011). Formation of toxic iodinated disinfection byproducts from compounds used in medical imaging. *Environmental Science & Technology*, 45(16), 6845–6854.
10. Schulz, M. et al. (2008). Transformation of the x-ray contrast medium iopromide in soil and biological wastewater treatment. *Environmental Science & Technology*, 42(19), 7207–7217.
11. Richardson, S. D. et al. (2008). Occurrence and mammalian cell toxicity of iodinated disinfection byproducts in drinking water. *Environmental Science & Technology*, 42(22), 8330–8338.
12. Richardson, S. D. et al. (2007). Occurrence, genotoxicity, and carcinogenicity of regulated and emerging disinfection byproducts in drinking water: A review and roadmap for research. *Mutation Research-Reviews in Mutation Research*, 636(1–3), 178–242.
13. Smith, E. M. et al. (2010). Comparison of byproduct formation in waters treated with chlorine and iodine: relevance to point-of-use treatment. *Environmental Science & Technology*, 44(22), 8446–8452.
14. Putschew, A., Schittko, S. & Jekel, M. (2001). Quantification of triiodinated benzene derivatives and x-ray contrast media in water samples by liquid chromatography-electrospray tandem mass spectrometry. *Journal of Chromatography A*, 930(1–2), 127–134.
15. Hirsch, R. et al. (2000). A sensitive method for the determination of iodine containing diagnostic agents in aqueous matrices using LC-electrospray-tandem-MS detection. *Fresenius Journal of Analytical Chemistry*, 366(8), 835–841.
16. Ternes, T. A. (2001). Analytical methods for the determination of pharmaceuticals in aqueous environmental samples. *Trac-Trends in Analytical Chemistry*, 20(8), 419–434.
17. Vanderford, B. J. et al. (2011). Assessment of sample preservation techniques for pharmaceuticals, personal care products, and steroids in surface and drinking water. *Analytical and Bioanalytical Chemistry*, 399(6), 2227–2234.

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© Agilent Technologies, Inc. 2012