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
## 2) 染色体解析、Fluorescence *in situ* hybridization (FISH)

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

Distinguishment of Morphologic and immunophenotypic type in 4 groups		DLBCL	DLBCL/BL	Burkitt (BL)	Lymphoblastic (LBL)
Morphology	Cell-size	Large (twice of lymphocytes)	small/medium-sized or mixture of large	small/medium-sized	small/medium-sized
	Nuclei	oval/round/varied	oval/round/varied	round	oval/round
	nucleoli	medium-sized/large/varied central/paracentral	medium-sized/large/varied paracentral basophilic	medium-sized paracentral basophilic	indistinct/multiple
	Chromatin	coarse/fine	finely clumped	finely clumped	dispersed/condensed
	Cytoplasm	amphophilic/basophilic scanty/abundant	basophilic scanty	basophilic scanty	light-blue/bluish grey scanty
	Growth pattern	diffuse/partial involvement	diffuse monotonous	diffuse monotonous cohesive	lipid vacuoles diffuse
	Mitotic features	sometimes	numerous	numerous	sometimes/numerous
Immunophenotype	Starry-sky	rare	frequent	frequent	focal
	Apoptosis	rare	frequent	frequent	focal
	TdT	-	-	-	+
	CD10	-/+	+-	+	+
	Bcl6	-/+	+-	+	-
	CCDN1	-+	-	-	-
	Bcl2	-/+	-+	-	-/+
Ki67 >90%/homogenous	rare	common	always	common	
Ki67 <90%/heterogeneous	common	sometimes	No	rare	

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文献：Matthiesen SH, Hansen CM. Fast and non-toxic *in situ* hybridization without blocking of repetitive sequences. *PLoS One*. 2012;7(7):e40675.

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
- 1) Swerdlow SH, et al. World Health Organization Classification of Tumors: Tumors of Haematopoietic and lymphoid tissues. *International Agency for Research on Cancer (IARC) Press*, Lyon France 2017.
- 2) Jaffe ES, et al. Hematopathology. Second edition *Elsevier* Philadelphia 2017
- 3) 磯部泰司 三浦偉久男：染色体・遺伝子検査：日本臨床 72: 418-429, 2014
- 4) 大島孝一：悪性リンパ腫：検査方法、診断、分類、コンパニオン診断まで *SRL 宝函* 2017 37(2) 11-21
- 5) Swerdlow SH, Campo E, Harris NL et al. (eds) WHO classification of tumor, Tumours of haematopoietic and lymphoid tissues. Lyon 2008.
- 6) Green TM, Young KH, Visco C, Xu-Monette ZY, Orazi A, Go RS, Nielsen O, Gadeberg OV, Mourits-Andersen T, Frederiksen M, Pedersen LM, Møller MB. Immunohistochemical double-hit score is a strong predictor of outcome in patients with diffuse large B-cell lymphoma treated with rituximab plus cyclophosphamide, doxorubicin, vincristine, and prednisone. *J Clin Oncol.* 2012 30:3460-7.
- 7) Johnson NA, Slack GW, Savage KJ, Connors JM, Ben-Neriah S, Rogic S, Scott DW, Tan KL, Steidl C, Sehn LH, Chan WC, Iqbal J, Meyer PN, Lenz G, Wright G, Rimsza LM, Valentino C, Brunhoeber P, Grogan TM, Brazier RM, Cook JR, Tubbs RR, Weisenburger DD, Campo E, Rosenwald A, Ott G, Delabie J, Holcroft C, Jaffe ES, Staudt LM, Gascoyne RD. Concurrent expression of MYC and BCL2 in diffuse large B-cell lymphoma treated with rituximab plus cyclophosphamide, doxorubicin, vincristine, and prednisone. *J Clin Oncol.* 2012 30:3452-9.
- 8) Hu S, Xu-Monette ZY, Tzankov A, Green T, Wu L, Balasubramanyam A, Liu WM, Visco C, Li Y, Miranda RN, Montes-Moreno S, Dybkaer K, Chiu A, Orazi A, Zu Y, Bhagat G, Richards KL, Hsi ED, Choi WW, Zhao X, van Krieken JH, Huang Q, Huh J, Ai W, Ponzoni M, Ferreri AJ, Zhou F, Slack GW, Gascoyne RD, Tu M, Variakojis D, Chen W, Go RS, Piris MA, Møller MB, Medeiros LJ, Young KH. MYC/BCL2 protein coexpression contributes to the inferior survival of activated B-cell subtype of diffuse large B-cell lymphoma and demonstrates high-risk gene expression signatures: a report from The International DLBCL Rituximab-CHOP Consortium Program. *Blood.* 2013 121:4021-31
- 9) Tomita N, Tokunaka M, Nakamura N, Takeuchi K, Koike J, Motomura S, Miyamoto K, Kikuchi A, Hyo R, Yakushijin Y, Masaki Y, Fujii S, Hayashi T, Ishigatsubo Y, Miura I. Clinicopathological features of lymphoma/leukemia patients carrying both BCL2 and MYC translocations. *Haematologica.* 2009 94:935-43.
- 10) Maruyama D, Watanabe T, Maeshima AM, Nomoto J, Taniguchi H, Azuma T, Mori M, Munakata W, Kim SW, Kobayashi Y, Matsuno Y, Tobinai K. Modified cyclophosphamide, vincristine, doxorubicin, and methotrexate (CODOX-M)/ifosfamide, etoposide, and cytarabine (IVAC) therapy with or without rituximab in Japanese adult patients with Burkitt lymphoma (BL) and B cell lymphoma, unclassifiable, with features intermediate between diffuse large B cell lymphoma and BL. *Int J Hematol.* 2010 92:732-43.
- 11) Clark Schneider KM, Banks PM, Collie AM, Lanigan CP, Manilich E, Durkin LM, Hill BT, Hsi ED. Dual expression of MYC and BCL2 proteins predicts worse outcomes in diffuse large B-cell lymphoma. *Leuk Lymphoma.* 2015 23:1-9

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© Agilent Technologies, Inc. 2019