



## Blue-White Color Screening

Blue-white color screening for recombinant plasmids is available when transforming this host strain (containing the *lacI<sup>q</sup>ZAM15* gene on the F' episome) with a plasmid that provides  $\alpha$ -complementation (e.g. Stratagene's pBluescript<sup>®</sup> II vector). When *lacZ* expression is induced by IPTG in the presence of the chromogenic substrate X-gal, colonies containing plasmids with inserts will be white, while colonies containing plasmids without inserts will be blue. If an insert is suspected to be toxic, plate the cells on media without X-gal and IPTG. Color screening will be eliminated, but lower levels of the potentially toxic protein will be expressed in the absence of IPTG.

### Critical Success Factors and Troubleshooting

**Aliquoting Cells:** Keep the cells on ice at all times during aliquoting. It is essential that the microcentrifuge tubes that the cells will be aliquoted into are placed on ice before the cells are thawed and that the cells are aliquoted directly into the pre-chilled tubes.

**Cuvette Gap Width:** Use a cuvette with a 0.1-cm gap to maximize the transformation efficiency and to minimize the possibility of arcing. A cuvette with a 0.2-cm gap is not recommended because the transformation efficiency is lower and the possibility of arcing is higher.

**Quantity and Volume of DNA:** The greatest efficiency is obtained from the transformation of 1  $\mu$ l of 0.01 ng/ $\mu$ l DNA per 40  $\mu$ l of cells. The volume of DNA may be increased to up to 4  $\mu$ l but the transformation efficiency may be reduced and the possibility of arcing may be increased if the DNA solution contains salts. A greater number of colonies may be obtained by increasing the amount of DNA added to the cells, although the overall efficiency may be lower.

**Ionic Strength of DNA Solution:** The sample DNA to be electroporated must be in a low-ionic-strength buffer, such as TE buffer or water. DNA samples containing too much salt will cause arcing at high voltage, possibly damaging both the sample and the machine.

**Reduction of Satellite Colonies:** Due to the high concentration of bacteria in electroporation, the transformation plates should be incubated at 37°C for less than 24 hours. If using an incubation period >24 hours, Stratagene recommends selecting for transformants on LB plates containing tetracycline (30  $\mu$ g/ml) in addition to the antibiotic used to select for the transforming plasmid. Including tetracycline helps suppress satellite colony growth and does not affect the transformation efficiency. In addition, all colonies that grow in the presence of tetracycline contain an F' episome, improving blue-white color selection.

**Plating the Transformation Mixture:** If plating <100  $\mu$ l of cells, pipet the cells into a 200- $\mu$ l pool of SOC medium and then spread the mixture with a sterile spreader. If plating  $\geq$ 100  $\mu$ l, the cells can be spread on the plates directly. Tilt and tap the spreader to remove the last drop of cells.

### Preparation of Media and Reagents

#### SOB Medium (per Liter)

20.0 g of tryptone  
5.0 g of yeast extract  
0.5 g of NaCl  
Add dH<sub>2</sub>O to a final volume of 1 liter and then autoclave  
Add 10 ml of filter-sterilized 1 M MgCl<sub>2</sub> and 10 ml of filter-sterilized 1 M MgSO<sub>4</sub> prior to use

#### SOC Medium (per 100 ml)

**Prepare immediately before use**  
2 ml of filter-sterilized 20% (w/v) glucose or 1 ml of filter-sterilized 2 M glucose  
SOB medium (autoclaved) to a final volume of 100 ml

#### LB Agar (per Liter)

10 g of NaCl  
10 g of tryptone  
5 g of yeast extract  
20 g of agar  
Add deionized H<sub>2</sub>O to a final volume of 1 liter  
Adjust pH to 7.0 with 5 N NaOH and then autoclave  
Pour into petri dishes (~25 ml/100-mm plate)

#### LB-Ampicillin Agar (per Liter)

1 liter of LB agar, autoclaved and cooled to 55°C  
Add 10 ml of 10 mg/ml filter-sterilized ampicillin  
Pour into petri dishes (~25 ml/100-mm plate)

#### TE Buffer

10 mM Tris-HCl (pH 7.5)  
1 mM EDTA

### Plates for Blue-White Color Screening

Prepare the LB agar and when adding the antibiotic, also add 5-bromo-4-chloro-3-indolyl- $\beta$ -D-galactopyranoside (X-gal) to a final concentration of 80  $\mu$ g/ml [prepared in dimethylformamide (DMF)] and isopropyl-1-thio- $\beta$ -D-galactopyranoside (IPTG) to a final concentration of 20 mM (prepared in sterile dH<sub>2</sub>O). Alternatively, 100  $\mu$ l of 10 mM IPTG and 100  $\mu$ l of 2% X-gal may be spread on solidified LB agar plates 30 minutes prior to plating the transformations. (For consistent color development across the plate, pipet the X-gal and the IPTG into a 100- $\mu$ l pool of SOC medium and then spread the mixture across the plate. Do not mix the IPTG and the X-gal before pipetting them into the pool of SOC medium because these chemicals may precipitate.)

### Limited Product Warranty

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

### Endnotes

\*U.S. Patent Nos. 6,338,965 and 6,040,184.

\*\*U.S. Patent Nos. 6,017,748, 5,552,314 and equivalent foreign patents.

pBluescript<sup>®</sup> and SURE<sup>®</sup> are registered trademarks of Stratagene in the United States.

**For in vitro use only. This certificate is a declaration of analysis at the time of manufacture.**