

QuikHyb Hybridization Solution

INSTRUCTION MANUAL

Catalog #201220 (250 ml) and #201221 (1 liter)

Revision A

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

QUIKHYB HYBRIDIZATION SOLUTION

CONTENTS

Materials Provided	1
Storage Conditions	1
Introduction	2
Radioactive Hybridization with QuikHyb Hybridization Solution	2
Prehybridization	2
Probe Preparation and Hybridization	3
Washing the Membrane.....	4
Detection	4
Stripping the Membrane for Reuse.....	4
Appendix I: Hybridization and Melting Temperatures	5
Calculating the Hybridization Temperature	5
Calculating the Melting Temperature ¹	5
Preparation of Media and Reagents	6
Reference	6
Endnotes	6
MSDS Information	6
Quick-Reference Protocol	8
Radioactive Hybridization.....	8

QuikHyb Hybridization Solution

MATERIALS PROVIDED

Material provided	Quantity	
	Catalog #201220	Catalog #201221
QuikHyb hybridization solution ^a	One 250-ml bottle	One 1-liter bottle

^a Store at 4°C. Do not freeze.

Note *The molar concentration of NaCl in QuikHyb hybridization solution is 0.5 M.*

STORAGE CONDITIONS

QuikHyb hybridization solution: 4°C

Revision A

© Agilent Technologies, Inc. 2009.

INTRODUCTION

Nucleic acid hybridization of a nucleic acid probe to a target sequence immobilized on a solid support such as a nylon membrane is a powerful technique routinely used in the analysis of gene structure and expression and in the diagnosis of disease.

Factors that affect the rate and efficiency of nucleic acid binding to a target sequence include hybridization time, temperature, concentration of the target sequence and concentration of the probe. Other variables which improve the hybridization reaction are the use of rate enhancers and blocking agents to suppress nonspecific hybridization.

QuikHyb hybridization solution* was developed for use in Southern, Northern, slot-blot and plaque-lift formats using randomly labeled radioactive and nonradioactive nucleic acid probes. This easy-to-use solution reduces the time required for hybridization from the customary 12–24 hours down to 1–2 hours, resulting in sensitive detection of target sequences.

RADIOACTIVE HYBRIDIZATION WITH QUIKHYB HYBRIDIZATION SOLUTION

Prehybridization

Notes *Gently invert the QuikHyb hybridization solution immediately before use in order to ensure a homogenous mixture.*

Prior to prehybridization (after crosslinking or baking), briefly dip the membrane in deionized water to remove excess salt from the membrane.

The membrane should be covered with a thin film of the QuikHyb solution at all times.

1. Use a minimum of 33 μ l of the QuikHyb hybridization solution/cm² of the blot. If this volume is not sufficient to cover the membrane at all times during the prehybridization and hybridization, add more solution until it does.

Minimum Total Volume

2 ml when using 50-ml conical tubes

10 ml when using heat-sealable bags or roller bottles

* Patent Pending.

2. For *double-stranded probes*, prehybridize the membrane in QuikHyb solution at 68°C for 10–20 minutes.

For *oligonucleotide probes and riboprobes*, calculate the melting temperature (T_m) (see *Appendix I: Hybridization and Melting Temperatures* for the mathematical formula). Prehybridize the membrane in QuikHyb solution for 10–20 minutes at 5–10°C below the T_m .

Probe Preparation and Hybridization

Add the labeled probe to 100 µl of 10 mg/ml sonicated salmon sperm DNA.

Suggested Probe Concentration

2.5×10^6 total counts/2 ml of hybridization solution

Specific Activity of the Probe

10^8 cpm/µg or greater

Double-Stranded Probes

Boil the double-stranded probe and salmon sperm DNA mixture for 2 minutes and then add the mixture to the prehybridization solution. Due to the viscosity of the QuikHyb hybridization solution, the probe may not mix well in the large glass tubes (roller bottles) used in hybridization ovens. Perform the following steps to ensure the probe is mixed well with the prehybridization solution.

1. Remove 1 ml of the prehybridization solution from the glass tube and mix the solution with the probe.
2. Place this mixture back in the glass tube and hybridize in QuikHyb solution at 68°C for 1 hour.

Oligonucleotide Probes and Riboprobes

Add the probe and salmon sperm DNA mixture directly to the prehybridization solution without boiling. Due to the viscosity of the QuikHyb hybridization solution, the probe may not mix well in the large glass tubes (roller bottles) used in hybridization ovens. Perform the following steps to ensure the probe is mixed well with the prehybridization solution.

1. Remove 1 ml of the prehybridization solution from the glass tube and mix the solution with the probe.
2. Place this mixture back in the glass tube and hybridize in QuikHyb solution at 5–10°C below the T_m for 1 hour (see *Appendix I: Hybridization and Melting Temperatures*).

Washing the Membrane

Perform the following washes (gentle agitation is required) for the double-stranded probes, oligonucleotide probes and riboprobes:

1. Wash twice for 15 minutes at room temperature with a 2× SSC buffer and 0.1% (w/v) SDS wash solution[§] (~2.5 ml/cm²).
2. Wash once for 30 minutes at 60°C with a 0.1× SSC buffer and 0.1% (w/v) SDS wash solution[§] for a high-stringency wash.

Detection

Wrap the membrane in plastic wrap and place the wrapped membrane on Kodak® X-OMAT® AR film with an intensifying screen at –80°C. Expose the film overnight.

Stripping the Membrane for Reuse

1. Heat the 0.1× SSC buffer and 0.1% (w/v) SDS wash solution to boiling.
2. In a glass dish, pour the 0.1× SSC buffer and 0.1% (w/v) SDS wash solution over the membrane and wash the membrane twice for 15 minutes.
3. Proceed with the prehybridization step for the next hybridization or store the blot in plastic wrap, desiccated.

[§] See *Preparation of Media and Reagents*.

APPENDIX I: HYBRIDIZATION AND MELTING TEMPERATURES

Calculating the Hybridization Temperature

The hybridization temperature for oligonucleotide probes equals 5–10°C below the T_m . For best results with oligonucleotide probes, perform hybridizations for 1 hour.

Calculating the Melting Temperature¹

Notes *The first method below overestimates the T_m of hybrids involving longer nucleotides. The second formula works only for Na^+ concentrations of ≤ 1 M.*

The molar concentration of NaCl in QuikHyb hybridization solution is 0.5 M.

Oligonucleotides Shorter than 18 Bases

$$T_m = 2^\circ\text{C}(\text{A} + \text{T}) + 4^\circ\text{C}(\text{G} + \text{C})$$

Oligonucleotides 14 Bases and Longer (up to 60–70 Nucleotides)

$$T_m = 81.5 + 16.6(\log_{10}[\text{Na}^+]) + 0.41(\%G + C) - (600/N)$$

where N is the chain length.

PREPARATION OF MEDIA AND REAGENTS

20× SSC Stock Solution 175.3 g of NaCl (3 M final concentration) 88.2 g of sodium citrate-trisodium salt (300 mM final concentration) Adjust the pH to 7.0 with HCl Add distilled water (dH ₂ O) to 1 liter	20% (w/v) SDS Stock Solution Dissolve 20 g of SDS in 90 ml of distilled water Mix well and heat to 68°C if necessary Then add dH ₂ O to a final volume of 100 ml
2× SSC Buffer and 0.1% (w/v) SDS Wash Solution 100 ml of 20× SSC buffer 5 ml of 20% (w/v) SDS dH ₂ O to 1 liter	0.1× SSC Buffer and 0.1% (w/v) SDS Wash Solution 5 ml of 20× SSC buffer 5 ml of 20% (w/v) SDS dH ₂ O to 1 liter

REFERENCE

1. Sambrook, J., Fritsch, E. F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

ENDNOTES

Kodak® and X-OMAT® are registered trademarks of Eastman Kodak Company.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.

STRATAGENE

An **Agilent Technologies** Division

QuikHyb Hybridization Solution

Catalog #201220 (250 ml) and #201221 (1 liter)

QUICK-REFERENCE PROTOCOL

Radioactive Hybridization

- ♦ Prehybridize for 15 minutes
- ♦ Denature the double-stranded probe
- ♦ Add probe to the prehybridizing blot
- ♦ Hybridize for 1 hour
- ♦ Wash twice for 15 minutes at room temperature with a $2\times$ SSC buffer and 0.1% (w/v) SDS wash solution
- ♦ Wash once for 30 minutes at higher stringency temperature with a $0.1\times$ SSC buffer and 0.1% (w/v) SDS wash solution
- ♦ Expose to Kodak® X-OMAT® AR film at -80°C overnight