



# **BL21-Gold Competent Cells, BL21-Gold(DE3) Competent Cells, and BL21- Gold(DE3)pLysS Competent Cells**

## **Instruction Manual**

**Catalog #230130 (BL21-Gold Competent Cells),  
#230132 [BL21-Gold(DE3) Competent Cells], and  
#230134 [BL21-Gold(DE3)pLysS Competent Cells]**  
Revision B

**Research Use Only. Not for Use in Diagnostic Procedures.**

230130-12

## LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

## ORDERING INFORMATION AND TECHNICAL SERVICES

### Email

techservices@agilent.com

### World Wide Web

www.genomics.agilent.com

### Telephone

| Location                 | Telephone   |
|--------------------------|---|
| United States and Canada | 800 227 9770  |
| Austria                  | 01 25125 6800   |
| Benelux                  | 02 404 92 22  |
| Denmark                  | 45 70 13 00 30  |
| Finland                  | 010 802 220   |
| France                   | 0810 446 446  |
| Germany                  | 0800 603 1000   |
| Italy                    | 800 012575  |
| Netherlands              | 020 547 2600  |
| Spain                    | 901 11 68 90  |
| Sweden                   | 08 506 4 8960   |
| Switzerland              | 0848 8035 60  |
| UK/Ireland               | 0845 712 5292   |
| All Other Countries      | Please visit <a href="http://www.genomics.agilent.com">www.genomics.agilent.com</a> and click <b>Contact Us</b> |

# **BL21-Gold Competent Cells, BL21-Gold(DE3) Competent Cells, and BL21-Gold(DE3)pLysS Competent Cells**

## **CONTENTS**

|  |           |
|--|-----------|
| <b>Materials Provided</b> .....  | <b>1</b>  |
| <b>Storage Conditions</b> .....  | <b>1</b>  |
| <b>Academic and Nonprofit Laboratory Assurance Letter</b> .....            | <b>2</b>  |
| <b>Introduction</b> .....  | <b>3</b>  |
| Host Strains and Genotypes.....  | 3         |
| Features .....   | 4         |
| <b>Transformation Guidelines</b> .....                                     | <b>5</b>  |
| Storage Conditions .....   | 5         |
| Use of 14-ml BD Falcon Polypropylene Round-Bottom Tubes.....               | 5         |
| Aliquoting Cells .....   | 5         |
| Quantity of DNA Added .....  | 5         |
| Length of the Heat Pulse .....   | 5         |
| <b>Transformation Protocol</b> .....                                       | <b>6</b>  |
| <b>Induction of Target Protein Using IPTG</b> .....                        | <b>7</b>  |
| <b>Induction of Target Protein by Infection with Lambda CE6</b> .....      | <b>8</b>  |
| Growth and Maintenance of High-Titer Bacteriophage Lambda CE6 Stocks ..... | 8         |
| Phage Amplification.....   | 9         |
| Induction of Target Protein by Infection with Lambda CE6.....              | 9         |
| <b>Troubleshooting</b> .....   | <b>10</b> |
| <b>Preparation of Media and Reagents</b> .....                             | <b>11</b> |
| <b>References</b> .....  | <b>12</b> |
| <b>Endnotes</b> .....  | <b>12</b> |
| <b>MSDS Information</b> .....  | <b>12</b> |

# BL21-Gold Competent Cells, BL21-Gold(DE3) Competent Cells, and BL21-Gold(DE3)pLysS Competent Cells

## MATERIALS PROVIDED

| Materials provided   | Tube color | Catalog number     |                    |                    | Efficiency                               |
|--|------------|--------------------|--------------------|--------------------|--|
|  |            | #230130            | #230132            | #230134            | (cfu/ $\mu$ g of pUC18 DNA) <sup>a</sup> |
| BL21-Gold competent cells  | Purple     | 10 $\times$ 0.1 ml | —                  | —                  | $\geq 1 \times 10^8$                     |
| BL21-Gold(DE3) competent cells                                     | Red-orange | —                  | 10 $\times$ 0.1 ml | —                  | $\geq 1 \times 10^8$                     |
| BL21-Gold(DE3)pLysS competent cells                                | White      | —                  | —                  | 10 $\times$ 0.1 ml | $\geq 1 \times 10^8$                     |
| pUC18 control plasmid (0.1 ng/ $\mu$ l in TE buffer <sup>b</sup> ) | —          | 10 $\mu$ l         | 10 $\mu$ l         | 10 $\mu$ l         | —  |

<sup>a</sup> These competent cell efficiencies are guaranteed when cells are used according to the specifications outlined in this instruction manual.

<sup>b</sup> See *Preparation of Media and Reagents*.

## STORAGE CONDITIONS

**Competent Cells:  $-80^{\circ}\text{C}$ . Do not store the competent cells in liquid nitrogen.**

## **ACADEMIC AND NONPROFIT LABORATORY ASSURANCE LETTER**

The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy and is protected by U.S. patents assigned to Brookhaven Science Associates (BSA). BSA will grant a nonexclusive license for use of this technology, including the enclosed materials, based on the following assurances:

1. These materials are to be used for noncommercial research purposes only. A separate license is required for any commercial use, including the use of these materials for research purposes or production purposes by any commercial entity. Information about commercial licenses may be obtained from the Office of Intellectual Property and Industrial Partnerships, Brookhaven National Laboratory, Bldg. 475D, Upton, New York, 11973 [telephone (631) 344-7134].

2. No materials that contain the cloned copy of T7 gene 1, the gene for T7 RNA polymerase, may be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of this license and agrees to be bound by its terms. This limitation applies to strains BL21-Gold, BL21-Gold(DE3), BL21-Gold(DE3)pLysS included in this kit and any derivatives you may make of them.

You may refuse this license by returning the enclosed materials unused. By keeping or using the enclosed materials, you agree to be bound by the terms of this license.

### **Commercial Entities Outside of the US**

The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy and is protected by U.S. Patents assigned to Brookhaven Science Associates (BSA). To protect its patent properties BSA requires commercial entities doing business in the United States, its Territories or Possessions to obtain a license to practice the technology. This applies for in-house research use of the T7 system as well as commercial manufacturing using the system. Commercial entities outside the U.S. that are doing business in the U.S., must also obtain a license in advance of purchasing T7 products. Commercial entities outside the U.S. that are using the T7 system solely for in-house research need not obtain a license if they do no business in the United States. However all customers, whether in the U.S. or outside the U.S. must agree to the terms and conditions in the Assurance Letter which accompanies the T7 products. Specifically, no materials that contain the cloned copy of T7 gene 1, the gene for T7 RNA polymerase, may be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of the assurance letter and agrees to be bound by its terms. This limitation applies to strains BL21-Gold, BL21-Gold(DE3), BL21-Gold(DE3)pLysS included in this kit and any derivatives you may make of them.

To obtain information about licensing, please contact the Office of Intellectual Property and Industrial Partnerships, Brookhaven National Laboratory, Building 475D, Upton, NY 11973 [telephone: 631-344-7134; Fax: 631-344-3729].

## INTRODUCTION

---

BL21-Gold competent cells, BL21-Gold(DE3) competent cells, and BL21-Gold(DE3)pLysS competent cells are improved versions of BL21 competent cells.<sup>3</sup> These high-performance competent cells provide increased transformation efficiency and produce high-quality miniprep DNA.

These BL21-Gold-derived expression strains are ideal for performing protein expression studies that utilize the T7 RNA polymerase promoter to direct high-level expression. Derived from *Escherichia coli* B, these expression strains naturally lack the Lon protease, which can degrade recombinant proteins. In addition, these strains are engineered to be deficient for a second protease, the OmpT protein.

These BL21-Gold-derived expression strains incorporate major improvements over the original BL21 strain. The BL21-Gold strains feature the Hte phenotype present in the highest efficiency Agilent competent cell strain, XL10-Gold.<sup>4</sup> The presence of the Hte phenotype increases the transformation efficiency of the BL21-Gold cells to  $>1 \times 10^8$  cfu/ $\mu$ g of pUC18 DNA. In addition, the gene that encodes endonuclease I (*endA*), which rapidly degrades plasmid DNA isolated by most miniprep procedures, is inactivated. These two improvements allow direct cloning of many protein expression constructs.

All three of the BL21-Gold-derived expression strains are resistant to tetracycline. In addition, the BL21-Gold(DE3)pLysS strain is resistant to chloramphenicol.

### Host Strains and Genotypes

| Host strain                             | Reference | Genotype   |
|---|-----------|--|
| BL21-Gold strain <sup>o</sup>           | 1         | <i>E. coli</i> B F <sup>-</sup> <i>ompT</i> <i>hsdS</i> (r <sub>B</sub> <sup>-</sup> m <sub>B</sub> <sup>-</sup> ) <i>dcm</i> <sup>+</sup> Tet <sup>r</sup> <i>gal</i> <i>endA</i> Hte   |
| BL21-Gold(DE3) strain <sup>o</sup>      | 1         | <i>E. coli</i> B F <sup>-</sup> <i>ompT</i> <i>hsdS</i> (r <sub>B</sub> <sup>-</sup> m <sub>B</sub> <sup>-</sup> ) <i>dcm</i> <sup>+</sup> Tet <sup>r</sup> <i>gal</i> $\lambda$ (DE3) <i>endA</i> Hte                           |
| BL21-Gold(DE3)pLysS strain <sup>o</sup> | 1         | <i>E. coli</i> B F <sup>-</sup> <i>ompT</i> <i>hsdS</i> (r <sub>B</sub> <sup>-</sup> m <sub>B</sub> <sup>-</sup> ) <i>dcm</i> <sup>+</sup> Tet <sup>r</sup> <i>gal</i> $\lambda$ (DE3) <i>endA</i> Hte [pLysS Cam <sup>r</sup> ] |

<sup>o</sup> This strain, a derivative of *E. coli* B, is a general protein expression strain that lacks both the Lon protease and the OmpT protease, which can degrade proteins during purification.<sup>2</sup> The Dcm methylase, naturally lacking in *E. coli* B, is inserted.

## Features

We offer three different competent cell lines that provide varying levels of expression control with T7 promoter-driven vectors such as the pCAL vectors<sup>‡</sup> and the pET vectors.<sup>‡</sup> BL21-Gold(DE3) is an all-purpose strain for high-level protein expression and easy induction. The BL21-Gold(DE3)pLysS strain provides tighter control of protein expression for expression of toxic proteins. When used with the CE6 bacteriophage, the BL21-Gold cells provide the tightest control of protein expression. Table I illustrates the features of the BL21-Gold expression strains for protein expression.

<sup>‡</sup> Available separately from Agilent. Visit <http://www.genomics.agilent.com> for details.

**Table I**

| Expression strain                   | Induction Method   | Advantages                               | Disadvantages  |
|-------------------------------------|--|--|--|
| BL21-Gold competent cells           | Infection with lambda bacteriophage CE6  | Tightest control of uninduced expression | Induction is not as efficient as DE3 derivatives<br>Induction (infection) process is more cumbersome |
| BL21-Gold(DE3) competent cells      | Isopropyl-1-thio-β-D-galactopyranoside (IPTG) induction of T7 polymerase from <i>lacUV5</i> promoter | High-level expression                    | Leaky expression of T7 polymerase can lead to uninduced expression of potentially toxic proteins     |
| BL21-Gold(DE3)pLysS competent cells | Isopropyl-1-thio-β-D-galactopyranoside (IPTG) induction of T7 polymerase from <i>lacUV5</i> promoter | Ease of induction                        | Slight inhibition of induced expression compared with BL21-Gold(DE3)                                 |

# TRANSFORMATION GUIDELINES

---

**Important** For optimal transformation efficiency, please read the guidelines outlined in the following sections before proceeding with the Transformation Protocol.

## Storage Conditions

The competent cells are very sensitive to slight variations in temperature. Storing the competent cells at the bottom of a  $-80^{\circ}\text{C}$  freezer directly from the dry ice shipping container is required in order to prevent a loss of transformation efficiency. Transferring tubes from one freezer to another may also result in a loss of efficiency. The transformation efficiency of the competent cells is guaranteed when the competent cells are used according to the specifications outlined in this instruction manual.

## Use of 14-ml BD Falcon Polypropylene Round-Bottom Tubes

It is important to use 14-ml BD Falcon polypropylene round-bottom tubes (BD Biosciences Catalog #352059) for the *Transformation Protocol* because the critical incubation period during the heat pulse is optimized specifically for the thickness and shape of these tubes.

## Aliquoting Cells

Store the competent cells on ice at all times while aliquoting. It is essential to place the 14-ml polypropylene tubes on ice before the competent cells are thawed and to aliquot the competent cells directly into the prechilled 14-ml polypropylene tubes. It is also important to use 100  $\mu\text{l}$  of competent cells/transformation. Using an inadequate volume of competent cells results in lower transformation efficiencies.

## Quantity of DNA Added

Greatest efficiencies (i.e., transformants/microgram of DNA) are observed when adding 1  $\mu\text{l}$  of DNA at a concentration of 0.1 ng/ $\mu\text{l}$  per 100  $\mu\text{l}$  of competent cells. Although the overall transformation efficiency may be lower, a greater number of colonies will be obtained when transforming up to 50 ng.

## Length of the Heat Pulse

Optimal transformation efficiencies are observed when transformation reactions are heat-pulsed for 20–25 seconds. Transformation efficiencies decrease sharply when the duration of the heat pulse is  $<20$  seconds or  $>25$  seconds.



## TRANSFORMATION PROTOCOL

---

1. Thaw the competent cells on ice.

**Note** *Store the competent cells **on ice at all times** while aliquoting. It is essential that the 14-ml BD Falcon polypropylene round-bottom tubes are placed on ice before the competent cells are thawed and that 100  $\mu$ l of competent cells are aliquoted directly into each **prechilled** polypropylene tube.*

2. Gently mix the competent cells. Aliquot 100  $\mu$ l of the competent cells into the appropriate number of prechilled 14-ml BD Falcon polypropylene round-bottom tubes. Prepare an additional 100- $\mu$ l aliquot of cells for use as a transformation control.
3. Add 1–50 ng of DNA to each transformation reaction and swirl gently. For the control transformation reaction, add 1  $\mu$ l of the pUC18 control plasmid to a separate 100- $\mu$ l aliquot of the competent cells and swirl gently.
4. Incubate the reactions on ice for 30 minutes.
5. Preheat SOC medium<sup>§</sup> in a 42°C water bath for use in step 8.
6. Heat-pulse each transformation reaction in a 42°C water bath for 20 seconds. **The duration of the heat pulse is critical for optimal transformation efficiencies (see *Length of the Heat Pulse*).**
7. Incubate the reactions on ice for 2 minutes.
8. Add 0.9 ml of preheated (42°C) SOC medium to each transformation reaction and incubate the reactions at 37°C for 1 hour with shaking at 225–250 rpm.
9. Concentrate the cells from the experimental transformation by centrifugation (200  $\times$  g for 3–5 minutes) and plate the entire transformation reaction (using a sterile spreader) onto a single LB agar plate<sup>§</sup> that contains the appropriate antibiotic.<sup>||</sup>

For the pUC18 control transformation, first place a 195- $\mu$ l pool of SOC medium on an LB–ampicillin agar plate.<sup>§</sup> Add 5  $\mu$ l of the control transformation reaction to the pool of SOC medium and then spread the mixture using a sterile spreader.<sup>||</sup>

10. Incubate the plates overnight at 37°C.

<sup>§</sup> See *Preparation of Media and Reagents*.

<sup>||</sup> When spreading the transformation reactions, tilt and tap the spreader to remove the last drop of cells.

## Transformation Summary for the pUC18 Control Plasmid

| Host strain                         | Plating quantity | Expected number of colonies | Efficiency (cfu/ $\mu$ g of pUC18 DNA) |
|-------------------------------------|------------------|-----------------------------|--|
| BL21-Gold competent cells           | 5 $\mu$ l        | >50                         | $\geq 1 \times 10^8$                   |
| BL21-Gold(DE3) competent cells      | 5 $\mu$ l        | >50                         | $\geq 1 \times 10^8$                   |
| BL21-Gold(DE3)pLysS competent cells | 5 $\mu$ l        | >50                         | $\geq 1 \times 10^8$                   |

## INDUCTION OF TARGET PROTEIN USING IPTG

The following induction protocol is a general guide for expression of genes under the control of IPTG-inducible promoters on an analytical scale (1 ml of induced culture). Most commonly, this protocol is used to analyze protein expression of individual transformants when using BL21(DE3) host strains in combination with plasmids containing T7 promoter constructs (e.g. pET vectors). Expression cassettes under the control of the *trp/lac* hybrid promoter, *tac*, can be also induced using this protocol. In the case of *tac* promoter constructs, non-DE3 lysogen strains can be employed as hosts.

**Note** *The transformation procedure described above will produce varying numbers of colonies depending on the transformation efficiency obtained for the expression plasmid. It is prudent to test more than one colony as colony-to-colony variations in protein expression are possible.*

1. Inoculate 1-ml aliquots of LB broth (see *Preparation of Media and Reagents*) containing the antibiotic required to maintain the expression plasmid with single colonies from the transformation. Shake at 220–250 rpm at 37°C overnight.

**Note** *For the BL21-Gold(DE3)pLysS host strain, the overnight culture must contain chloramphenicol at a final concentration of 50  $\mu$ g/ml in addition to the antibiotic required to maintain the expression plasmid. Chloramphenicol serves to maintain the pACYC-based plasmid carrying the T7 lysozyme gene derivative.*

2. The next morning, pipet 50  $\mu$ l of each culture into fresh 1-ml aliquots of LB broth containing no selection antibiotics. Incubate these cultures with shaking at 220–250 rpm at 37°C for 2 hours.
3. Pipet 100  $\mu$ l of each of the cultures into clean microcentrifuge tubes and place the tubes on ice until needed for gel analysis. These will serve as the non-induced control samples.
4. To the rest of the culture in each tube add IPTG to a final concentration of 1 mM. Incubate with shaking at 220–250 rpm at 37°C for 2 hours.

**Note** *These values for IPTG concentration and induction time are starting values only and may require optimization depending on the gene expressed.*

5. After the induction period, place the cultures on ice.
6. Pipet 20  $\mu$ l of each of the induced cultures into clean microcentrifuge tubes. Add 20  $\mu$ l of 2 $\times$  SDS gel sample buffer<sup>§</sup> to each microcentrifuge tube.
7. Mix the non-induced samples held on ice to resuspend the cells. Pipet 20  $\mu$ l from each tube into a clean microcentrifuge tube. Add 20  $\mu$ l of 2 $\times$  SDS gel sample buffer to each of the 20- $\mu$ l aliquots of cells.
8. Heat all tubes to 95°C for 5 minutes. Load the associated non-induced and induced samples in adjacent lanes for analysis by SDS-PAGE. Stain the protein gel with Coomassie<sup>®</sup> Brilliant Blue stain.

## **INDUCTION OF TARGET PROTEIN BY INFECTION WITH LAMBDA CE6**

---

Expression of genes under the control of the T7 promoter (e.g. genes in pET vectors) can be achieved in non-DE3 lysogen host strains (e.g. BL21) if the strain harboring the expression plasmid is subsequently infected with lambda CE6. Lambda CE6 expresses T7 polymerase, which in turn drives the transcription of the gene downstream of the T7 promoter. The following protocols describe the growth and maintenance of lambda CE6 and the use of lambda CE6 for infecting host strains. We offer the Lambda CE6 Induction Kit (Catalog #235200) for use in protein expression protocols which incorporate CE6 infection.

### **Growth and Maintenance of High-Titer Bacteriophage Lambda CE6 Stocks**

1. Inoculate 5 ml of modified\* NZY broth<sup>§</sup> with a single colony of LE392 host cells. Shake overnight at 37°C at 220–250 rpm.
2. Centrifuge the overnight culture for 15 minutes at 1700–2000  $\times$  g at 4°C. Resuspend the cells in 10 mM MgSO<sub>4</sub> to a final OD<sub>600</sub> of 0.5.
3. Combine 250  $\mu$ l of cells (at OD<sub>600</sub> = 0.5) with 1  $\times$  10<sup>6</sup> pfu of lambda CE6 in 14-ml BD Falcon polypropylene round-bottom tubes in triplicate. Incubate at 37°C for 15 minutes.
4. Add 3 ml of melted NZY top agar<sup>§</sup> to each cell suspension and plate on warm agarose plates. <sup>§</sup> Incubate the plates overnight at 37°C.

\* NZY broth to be used for lambda infection protocols should be supplemented with maltose at a final concentration of 0.2%. Add 1 ml of 20% maltose solution (filter-sterilized) per 100 ml of NZY broth to achieve the correct final concentration of maltose.

<sup>§</sup> See *Preparation of Media and Reagents*.

5. Flood each plate with 5 ml of SM solution (see *Preparation of Media and Reagents*) and rock the plates for 2 hours at room temperature.
6. Remove the SM solution (which contains the lambda CE6) from each plate and pool the volumes in a 50-ml conical tube.
7. Centrifuge the SM solution at  $1700\text{--}2000 \times g$  for 15 minutes at  $4^{\circ}\text{C}$ .
8. Remove the supernatant and determine the titer of the solution.
9. Store the lambda CE6 stock at  $4^{\circ}\text{C}$ .

## Phage Amplification

If the titer drops over time, or if more phage are needed, grow up LE392 cells in 10 ml of medium and add bacteriophage lambda CE6 at a multiplicity of infection of 1:1000 (CE6-to-cell ratio). Continue growing the culture at  $37^{\circ}\text{C}$  for 5–6 hours and spin down the cellular debris. Titer of the supernatant should be  $\geq 5.0 \times 10^9$  pfu/ml. For general information regarding phage amplification, see reference 5.

## Induction of Target Protein by Infection with Lambda CE6

**Note** *This protocol is designed for induction in 50-ml culture volumes. If induction of a larger volume of culture is desired, it will be necessary to increase the volume of the overnight culture in step 1. The increased volume of overnight culture is necessary to achieve the required cell density ( $A_{600} \leq 1$ ) in the larger volume of broth the following day..*

1. Inoculate 5 ml of modified\* NZY broth containing the antibiotic required to maintain the expression plasmid with a single colony of BL21 cells (not a DE3 lysogen) harboring the expression plasmid. Shake overnight at  $37^{\circ}\text{C}$  at 200–250 rpm.
2. In the morning, centrifuge 1.0 ml of the overnight culture, resuspend the cells in 1.0 ml of fresh modified\* NZY broth, and pipet the resuspended cells into a flask containing 50 ml of fresh modified\* NZY broth (no selection antibiotics).
3. Record the  $A_{600}$  of the diluted culture. It should be  $\leq 0.1$ . If the  $A_{600}$  is  $> 0.1$ , use more fresh modified\* NZY broth to dilute the culture to  $A_{600} \leq 0.1$ . If the  $A_{600}$  is  $< 0.1$ , the time required to reach an  $A_{600}$  of 0.3 (in step 4) will be extended.
4. Grow the culture to an  $A_{600}$  of 0.3 and then add glucose to a final concentration of 4 mg/ml (e.g. 1.0 ml of a 20% glucose solution to the 50-ml culture).

\* NZY broth to be used for lambda infection protocols should be supplemented with maltose at a final concentration of 0.2%. Add 1 ml of 20% maltose solution (filter-sterilized) per 100 ml of NZY broth to achieve the correct final concentration of maltose.

5. Grow the culture to an  $A_{600}$  of 0.6–1.0 and then add  $MgSO_4$  to a final concentration of 10 mM (e.g. 500  $\mu$ l of a 1.0 M solution of  $MgSO_4$  to the 50-ml culture).
6. Remove a portion of the culture to serve as the uninduced control and infect the rest with bacteriophage lambda CE6 at a multiplicity of infection (MOI) of 5–10 particles per cell. (To optimize induction, cultures may be split into 3 or 4 aliquots and infected with varying dilutions of bacteriophage lambda CE6. The subsequent induction can be monitored by SDS-PAGE or by a functional assay, if available.)
7. Grow the culture for 2–3 hours after infection with lambda CE6.
8. Remove 5–20  $\mu$ l of the culture for determination by SDS-PAGE, and harvest the remaining culture by centrifugation. Store the pellets at  $-70^\circ C$ .

**Note** *If induction will be monitored using Coomassie stain, silver stain, or another nonspecific protein stain, we recommend running a control of CE6-infected BL21 cells harboring the plasmid without a cloned insert.*

## TROUBLESHOOTING

| Observation         | Suggestion  |
|---------------------|---|
| Plasmid instability | Inserts containing repeated sequence(s) or secondary structure may undergo rearrangement if established in the BL21-Gold series of competent cells as they are <i>recA</i> <sup>+</sup> . To minimize this possibility, establish the insert in a recombination-deficient host strain, such as SURE competent cells <sup>‡</sup> ( <i>recB</i> <sup>-</sup> and <i>recJ</i> <sup>-</sup> ) or XL1-Blue competent cells <sup>‡</sup> ( <i>recA</i> <sup>-</sup> ), prior to expression in BL21-Gold cells. |
| Clone toxicity      | More tightly controlled induction may be achieved by infecting BL21-Gold cells with bacteriophage CE6 than by IPTG-induction of BL21-Gold (DE3) cells. When using IPTG induction, uninduced expression levels are reduced when using the BL21-Gold (DE3) $\rho$ LysS host strain, compared to the BL21-Gold (DE3) host strain. See <i>Features</i> in the <i>Introduction</i> .   |

<sup>‡</sup> Available separately from Agilent. Visit <http://www.genomics.agilent.com> for details.

## PREPARATION OF MEDIA AND REAGENTS

|   |   |
|---|---|
| <p><b>LB Broth (per Liter)</b><br/>           10 g of NaCl<br/>           10 g of tryptone<br/>           5 g of yeast extract<br/>           Add dH<sub>2</sub>O to a final volume of 1 liter<br/>           Adjust pH to 7.0 with 5 N NaOH<br/>           Autoclave<br/>           Pour into petri dishes (~25 ml/100-mm plate)</p>   | <p><b>LB Agar (per Liter)</b><br/>           10 g of NaCl<br/>           10 g of tryptone<br/>           5 g of yeast extract<br/>           20 g of agar<br/>           Add deionized H<sub>2</sub>O to a final volume of 1 liter<br/>           Adjust pH to 7.0 with 5 N NaOH<br/>           Autoclave<br/>           Cool to 55°C, and then add antibiotic, if required<br/>           Pour into petri dishes (~25 ml/100-mm plate)</p> |
| <p><b>LB–Ampicillin Agar (per Liter)</b><br/>           1 liter of LB agar, autoclaved<br/>           Cool to 55°C<br/>           Add 10 ml of 10-mg/ml filter-sterilized ampicillin<br/>           Pour into petri dishes (~25 ml/100-mm plate)</p>  | <p><b>SOC Medium (per 100 ml)</b><br/> <b>Note</b> <i>This medium should be prepared immediately before use</i><br/>           1 ml of a 2 M filter-sterilized glucose solution<br/> <b>or</b> 2 ml of 20% (w/v) glucose<br/>           SOB medium to a final volume of 100 ml<br/>           Filter sterilize</p>  |
| <p><b>SOB Medium (per Liter)</b><br/>           20.0 g of tryptone<br/>           5.0 g of yeast extract<br/>           0.5 g of NaCl<br/>           Add deionized H<sub>2</sub>O to a final volume of 1 liter<br/>           Autoclave<br/>           Add 10 ml of filter-sterilized 1 M MgCl<sub>2</sub> and 10 ml of filter-sterilized 1 M MgSO<sub>4</sub> prior to use</p> | <p><b>2× SDS gel sample buffer</b><br/>           100 mM Tris-HCl (pH 6.5)<br/>           4% SDS (electrophoresis grade)<br/>           0.2% bromophenol blue<br/>           20% glycerol<br/><br/> <b>Note</b> <i>Add dithiothreitol to a final concentration in the 2× buffer of 200 mM prior to use. This sample buffer is useful for denaturing, discontinuous acrylamide gel systems only.</i></p>                                     |
| <p><b>TE Buffer</b><br/>           10 mM Tris-HCl (pH 7.5)<br/>           1 mM EDTA</p>   |   |

(table continues on following page)

(table continues from previous page)

|  |   |
|--|---|
| <b>SM Solution</b><br>5 g of NaCl<br>2 g of MgSO <sub>4</sub> · 7H <sub>2</sub> O<br>50 ml of 1 M Tris-HCl (pH 7.5)<br>5 ml 2% gelatin<br>Add deionized H <sub>2</sub> O to a final volume of 1 liter<br>Adjust the pH to 7.5<br>Autoclave   | <b>NZY Broth, (per Liter)</b><br>5 g of NaCl<br>2 g of MgSO <sub>4</sub> · 7H <sub>2</sub> O<br>5 g of yeast extract<br>10 g of NZ amine (casein hydrolysate)<br>Add deionized H <sub>2</sub> O to a final volume of 1 liter<br>Adjust the pH to 7.5 with NaOH<br>Autoclave |
| <b>Agarose Plates (per Liter)</b><br>Melt 20 g of agarose in 500 ml of deionized H <sub>2</sub> O<br>Add the following:<br>5 g of NaCl<br>5 g of yeast extract<br>10 g of tryptone<br>Add deionized H <sub>2</sub> O to a final volume of 1 liter<br>Autoclave<br>Pour into petri dishes (~25 ml/100-mm plate) | <b>NZY Top Agar (per Liter)</b><br>Prepare 1 liter of NZY broth<br>Add 0.7% (w/v) agarose<br>Autoclave<br>Prior to use, melt agar in microwave, then hold at 48°C in water bath   |

## REFERENCES

1. Weiner, M. P., Anderson, C., Jerpseth, B., Wells, S., Johnson-Browne, B. *et al.* (1994) *Strategies* 7(2):41–43.
2. Grodberg, J. and Dunn, J. J. (1988) *J Bacteriol* 170(3):1245-53.
3. Jerpseth, M., Jerpseth, B., Briester, L. and Greener, A. (1998) *Strategies* 11(1):3-4.
4. Jerpseth, B., Callahan, M. and Greener, A. (1997) *Strategies* 10(2):37–38.
5. Sambrook, J., Fritsch, E. F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

## ENDNOTES

Coomassie® is a registered trademark of Imperial Chemical Industries.

## MSDS INFORMATION

Material Safety Data Sheets (MSDSs) are provided online at <http://www.genomics.agilent.com>. MSDS documents are not included with product shipments.