



## TRANSFORMATION GUIDELINES AND TROUBLESHOOTING

**Storage Conditions:** Supercompetent cells are very sensitive to even small variations in temperature and must be stored at the bottom of a  $-80^{\circ}\text{C}$  freezer. Transferring tubes from one freezer to another may result in a loss of efficiency.

**Use of  $\beta$ -Mercaptoethanol ( $\beta$ -ME) and NZY<sup>+</sup> broth:** For the highest transformation efficiencies ( $1 \times 10^9$  cfu/ $\mu\text{g}$  pUC18 DNA), use the optimized transformation protocol that includes treatment of the cells with  $\beta$ -ME (see step 4) and use of NZY<sup>+</sup> as the outgrowth medium (see steps 1 and 10). When the  $\beta$ -ME step is omitted from the protocol and SOC is used as the outgrowth medium, expect transformation efficiencies of approximately 90% of those obtained using the optimized protocol. The  $\beta$ -ME mix provided is diluted and ready to use. Stratagene cannot guarantee results with  $\beta$ -ME from other sources.

**Quantity and Volume of DNA:** The greatest efficiency is obtained from the transformation of 1  $\mu\text{l}$  of 0.01 ng/ $\mu\text{l}$  of supercoiled DNA per aliquot of cells. A greater number of colonies may be obtained by transforming up to 50 ng DNA, although the resulting efficiency may be lower. The volume of the DNA solution added to the reaction may be increased to up to 5  $\mu\text{l}$ , but the transformation efficiency may be reduced.

**Heat Pulse Duration and Temperature:** Optimal transformation efficiency is observed when cells are heat-pulsed at  $54^{\circ}\text{C}$  for 60 seconds in the tubes provided. Efficiency decreases sharply when either the temperature or the duration of the heat pulse is altered.

**Blue-White Color Screening:** Blue-white color screening for recombinant plasmids is available when transforming the host strain with a plasmid that provides  $\alpha$ -complementation (e.g. Stratagene's pBluescript<sup>®</sup> II). When lacZ expression is induced by IPTG in the presence of the chromogenic substrate X-gal, colonies containing plasmids with inserts will remain white, while colonies containing plasmids without inserts will be blue. When performing blue-white color screening, incubate the LB agar plates containing IPTG and X-gal at  $37^{\circ}\text{C}$  for at least 17 hours to allow color development. The blue color can be enhanced by subsequent incubation of the plates for two hours at  $4^{\circ}\text{C}$ . If an insert is suspected to be toxic, plate the cells on media without X-gal and IPTG. Color screening will be eliminated, but lower levels of the potentially toxic protein will be expressed in the absence of IPTG.

## PREPARATION OF MEDIA AND REAGENTS

<b>LB Agar (per Liter)</b> 10 g of NaCl 10 g of tryptone 5 g of yeast extract 20 g of agar Add deionized H <sub>2</sub> O to a final volume of 1 liter Adjust pH to 7.0 with 5 N NaOH Autoclave Pour into petri dishes (~25 ml/100-mm plate)	<b>SOB Broth (per Liter)</b> 20.0 g of tryptone 5.0 g of yeast extract 0.5 g of NaCl Add deionized H <sub>2</sub> O to a final volume of 1 liter Autoclave Add 10 ml of filter-sterilized 1 M MgCl <sub>2</sub> and 10 ml of filter-sterilized 1 M MgSO <sub>4</sub> prior to use
<b>NZY<sup>+</sup> Broth (per Liter)</b> 10 g of NZ amine (casein hydrolysate) 5 g of yeast extract 5 g of NaCl Add deionized H <sub>2</sub> O to a final volume of 1 liter Adjust to pH 7.5 using NaOH Autoclave Add the following filter-sterilized supplements prior to use: 12.5 ml of 1 M MgCl <sub>2</sub> 12.5 ml of 1 M MgSO <sub>4</sub> 20 ml of 20% (w/v) glucose (or 10 ml of 2 M glucose)	<b>SOC Broth (per 100 ml)</b> <i>Note This medium should be prepared immediately before use.</i>  2 ml of filter-sterilized 20% (w/v) glucose or 1 ml of filter-sterilized 2 M glucose SOB medium (autoclaved) to a final volume of 100 ml
	<b>LB-Ampicillin Agar (per Liter)</b> 1 liter of LB agar, autoclaved Cool to $55^{\circ}\text{C}$ Add 10 ml of 10-mg/ml filter-sterilized ampicillin Pour into petri dishes (~25 ml/100-mm plate)

## Preparation of Agar Plates for Blue-White Color Screening

To prepare plates for blue-white screening, prepare LB agar as indicated above. When adding the antibiotic, also add 5-bromo-4-chloro-3-indolyl- $\beta$ -D-galactopyranoside (X-gal) to a final concentration of 80  $\mu\text{g}/\text{ml}$  [prepared in dimethylformamide (DMF)] and isopropyl-1-thio- $\beta$ -D-galactopyranoside (IPTG) to a final concentration of 20 mM (prepared in sterile dH<sub>2</sub>O). Alternatively, 100  $\mu\text{l}$  of 10 mM IPTG and 100  $\mu\text{l}$  of 2% X-gal may be spread on solidified LB agar plates 30 minutes prior to plating the transformations. (For consistent color development across the plate, pipet the X-gal and the IPTG into a 100- $\mu\text{l}$  pool of SOC medium and then spread the mixture across the plate. Do not mix the IPTG and the X-gal before pipetting them into the pool of SOC medium because these chemicals may precipitate.)

## LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

## ENDNOTES

pBluescript<sup>®</sup>, SoloPack<sup>®</sup> and XL10-Gold<sup>®</sup> are registered trademarks of Stratagene in the United States.

**Stratagene Technical Services**  
USA/Canada (Toll-free) 800 894 1304  
Europe (Toll-free) 00800 7400 7400  
  
**Email** tech\_services@stratagene.com  
**World Wide Web** www.stratagene.com