

Lambda FIX II/Xho I Partial Fill-In Vector Kit

INSTRUCTION MANUAL

Catalog #248211 (*Xho* I-Partial Fill-In), #248612 (*Xho* I/Gigapack III Gold Packaging Extract), and #248712 (*Xho* I/Gigapack III XL Packaging Extract)

Revision A

For In Vitro Use Only
248211-12

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

Lambda FIX II/Xho I Partial Fill-In Vector Kit

CONTENTS

Materials Provided	1
Storage Conditions	1
Introduction	2
Preparation of Host Strains	3
Host Strain Genotypes.....	3
Growing and Maintaining the Host Strains.....	3
Preparing a –80°C Bacterial Glycerol Stock.....	4
Preparing the Insert	4
Ligating the Insert	5
Packaging	6
General Information.....	6
Gigapack III XL Packaging Extract.....	6
Packaging Instructions.....	6
Titering the Packaging Reaction	9
Amplifying the Library	10
Performing Plaque Lifts	11
Hybridizing and Screening	12
Rapid Restriction Mapping	13
Troubleshooting	14
Preparation of Media and Reagents	15
References	16
Endnotes	16
MSDS Information	16

Lambda FIX II/Xho I Partial Fill-In Vector Kit

MATERIALS PROVIDED

Materials provided	Quantity		
	Catalog #248211	Catalog #248612	Catalog #248712
Lambda FIX II vector predigested with Xho I ^a (1 µg/µl in TE buffer)	10 µg	10 µg	10 µg
pMF/BamH I test insert (12 kb) ^b (0.5 µg/µl in TE buffer)	5 µg	5 µg	5 µg
Host strains ^c			
XL1-Blue MRA strain	0.5-ml bacterial glycerol stock	0.5-ml bacterial glycerol stock	0.5-ml bacterial glycerol stock
XL1-Blue MRA (P2) strain	0.5-ml bacterial glycerol stock	0.5-ml bacterial glycerol stock	0.5-ml bacterial glycerol stock
Gigapack III Gold-11 packaging extract	—	11 × 25 µl	—
Gigapack III XL-11 packaging extract	—	—	11 × 25 µl
λcl857 Sam7 wild-type lambda control DNA	—	1.05 µg	1.05 µg
VCS257 host strain ^{c,d}	—	1 ml bacterial glycerol stock	1 ml bacterial glycerol stock

^a The first two nucleotides of the Xho I site are filled in, leaving 3'-CT-5' overhangs

^b The first two nucleotides of the BamH I site are filled in, leaving 3'-AG-5' overhangs.

^c For host strain handling instructions, please see *Preparation of Host Strains*.

^d The VCS257 host strain, included for plating the λcl857 Sam7 positive control, is a derivative of DP50 *supF* and should be used only when plating the packaged test DNA. The control DNA used with Gigapack III Gold packaging extract requires a *supF* mutation in the bacterial host to plate efficiently.

STORAGE CONDITIONS

Lambda FIX II Vector and pMF/BamH I test insert: Aliquot and store at –20°C. Do not pass through more than two freeze-thaw cycles. For short-term storage, store at 4°C for 1 month.

Bacterial Host Strains: –80°C; avoid repeated freeze-thaw cycles.

Gigapack III Packaging Extracts: Gigapack III packaging extract is very sensitive to slight variations in temperature. Store the packaging extracts at the bottom of a –80°C freezer directly from the dry ice shipping container in order to prevent a loss of packaging efficiency. Transferring tubes from one freezer to another may also result in a loss of efficiency. **Do not allow the packaging extracts to thaw!** Do not store the packaging extracts in liquid nitrogen as the tubes may explode.

λcl857 Sam7 wild-type lambda control DNA: Shipped frozen; store immediately at –80°C.

INTRODUCTION

The Lambda FIX II vector is a replacement vector used for cloning large fragments of genomic DNA (see Figures 1 and 2). The Lambda FIX II system takes advantage of *spi* (sensitive to P2 inhibition) selection. Lambda phages containing active *red* and *gam* genes are unable to grow on host strains that contain P2 phage lysogens. Lambda phages without these genes are able to grow on strains lysogenic for P2 such as XL1-Blue MRA (P2), a P2 lysogen of XL1-Blue MRA. The *red* and *gam* genes in the Lambda FIX II DNA are located on the stuffer fragment; therefore, the wild-type Lambda FIX II phage cannot grow on XL1-Blue MRA (P2). When the stuffer fragment is replaced by an insert, the recombinant Lambda FIX II vector becomes Red⁻/Gam⁻, and the phage is able to grow on the P2 lysogenic strain. Therefore, by plating the library on the XL1-Blue MRA (P2) strain, only recombinant phages are allowed to grow. The strain XL1-Blue MRA is also provided as a control strain and later for growth of the recombinant when the selection is no longer necessary. The unique arrangement of the polylinker for the Lambda FIX II vector permits the isolation of the insert and flanking T3 and T7 bacteriophage promoters as an intact cassette by digestion with *Not* I. T3 and T7 promoters flanking the insertion sites can be used to generate end-specific RNA probes for use in chromosomal walking and restriction mapping.




Figure 1 Map of the Lambda FIX II replacement vector.



Figure 2 Multiple cloning site sequence of the Lambda FIX II replacement vector.

PREPARATION OF HOST STRAINS

Host Strain Genotypes

Host strain	Genotype
XL1-Blue MRA strain	$\Delta(mcrA)183 \Delta(mcrCB-hsdSMR-mrr)173 endA1 supE44 thi-1 gyrA96 relA1 lac$
XL1-Blue MRA (P2) strain	XL1-Blue MRA (P2 lysogen)

Growing and Maintaining the Host Strains

The bacterial host strains are shipped as bacterial glycerol stocks. For the appropriate media, please refer to the following table:

Host strain	Agar plates for bacterial streak	Medium for bacterial glycerol stock	Medium for bacterial cultures for titering phage (final concentration)
XL1-Blue MRA strain ^a	LB	LB	LB with 0.2% (w/v) maltose–10 mM MgSO ₄
XL1-Blue MRA (P2) strain ^a	LB	LB	LB with 0.2% (w/v) maltose–10 mM MgSO ₄
VCS257 strain ^b	LB	LB	LB with 0.2% (w/v) maltose–10 mM MgSO ₄

^a The XL1-Blue MRA and XL1-Blue MRA (P2) host strains are modified to enhance the stability of clones containing methylated DNA; in addition, these strains enhance the stability of nonstandard DNA structures.

^b For use with Gigapack III packaging extract and wild-type control only. Supplied with Gigapack III packaging extract.

On arrival, prepare the following from the bacterial glycerol stock using the appropriate media as indicated in the previous table:

Note *The host strains may thaw during shipment. The vials should be stored immediately at –20° or –80°C, but most strains remain viable longer if stored at –80°C. It is also best to avoid repeated thawing of the host strains in order to maintain extended viability.*

1. Revive the stored cells by scraping off splinters of solid ice with a sterile wire loop.
2. Streak the splinters onto a LB agar plate (see *Preparation of Media and Reagents*).
3. Incubate the plate overnight at 37°C.
4. Seal the plate with Parafilm® laboratory film and store the plate at 4°C for up to 1 week.
5. Restreak the cells onto a fresh plate every week.

Preparing a -80°C Bacterial Glycerol Stock

1. In a sterile 50-ml conical tube, inoculate 10 ml of appropriate liquid medium with one colony from the plate. Grow the cells to late log phase.
2. Add 4.5 ml of a sterile glycerol–liquid medium solution (5 ml of glycerol + 5 ml of appropriate medium) to the bacterial culture from step 1. Mix well.
3. Aliquot into sterile centrifuge tubes (1 ml/tube).

This preparation may be stored at -20°C for 1–2 years or at -80°C for more than 2 years.

PREPARING THE INSERT

Note *In order to ligate your digested insert DNA to the Lambda FIX II vector, it is important to fill-in the first two nucleotides of the BamH I-compatible site. This strategy makes it virtually impossible for the insert to religate to itself, thus lowering the background and the multiple inserts. Our Klenow Fill-In Kit (catalog #200410) provides the components required to perform partial fill-in reactions to prepare genomic DNA for insertion into the Lambda FIX II vector. An example of a standard nucleotide fill-in reaction is as follows:*

Assume that 50 μg of the insert DNA is digested with either *Mbo* I, *Sau* 3A or *Bam*H I.

Prepare the following reaction mixture in a microcentrifuge tube:

50 μg of digested insert DNA
30 μl of 10 \times fill-in buffer
5 μl of 10 mM dGTP
5 μl of 10 mM dATP
15 U of Klenow
Water to a final volume of 300 μl

Incubate at room temperature for 15 minutes.

Add 33 μl of 10 \times STE buffer
Add 100 μl of 1 \times STE buffer
Extract with phenol–chloroform
Repeat the phenol–chloroform extraction
Ethanol precipitate
Spin down the DNA
Lyophilize the pellet
Resuspend the pellet in 25 μl of Tris-EDTA

LIGATING THE INSERT

Note *In all ligations, the final glycerol content should be less than 5% (v/v). Do not exceed 5% (v/v) glycerol. Due to the high molecular weight of the lambda vector, the contents may be very viscous. It is important to microcentrifuge the contents of the lambda vector tube briefly at $11,000 \times g$, then gently mix the solution by stirring with a yellow pipet tip prior to pipetting.*

Prepare the following reaction mixture in a microcentrifuge tube:

1.0 μl of the Lambda FIX II predigested with *Xho* I (1 μg)
0.8 μl of pMF/*Bam*H I insert (0.4 μg)
0.5 μl of 10 \times ligase buffer (see *Preparation of Media Reagents*)
0.5 μl of 10 mM rATP (pH 7.5)
2 U of T4 DNA ligase
Water up to a final volume of 5 μl

Incubate the ligation at 4°C overnight.

When ligating inserts containing the DNA of interest, use a volume of insert DNA up to 2.5 μl . Use an equal molar ratio of the *Bam*H I- *Mbo* I-, or *Sau* 3A-digested insert with the Lambda FIX II vector. The Lambda FIX II vector can accommodate inserts ranging from 9–23 kb. If ligating a 20,000-bp insert to the vector, use 0.4 μg of insert for every 1 μg of vector. If the insert used is free from contaminants and contains a high percentage of ligatable ends, expect about 5×10^5 – 1×10^7 recombinant plaques when using high-efficiency packaging extracts, such as Gigapack III Plus or Gigapack III Gold packaging extracts.

PACKAGING

General Information

Packaging extracts are used to package recombinant lambda phage with high efficiency. The single-tube format of Gigapack III packaging extract simplifies the packaging procedure and increases the efficiency and representation of libraries constructed from highly methylated DNA. Each packaging extract is restriction minus (HsdR⁻ McrA⁻ McrBC⁻ McrF⁻ Mrr⁻) to optimize packaging efficiency and library representation. When used in conjunction with restriction-deficient plating cultures, Gigapack III packaging extract improves the quality of DNA libraries constructed from methylated DNA.^{1,2,3,4}

Optimal packaging efficiencies are obtained with lambda DNAs that are concatemeric. Ligations should be carried out at DNA concentrations of 0.2 µg/µl or greater, which favors concatemers and not circular DNA molecules that only contain one *cos* site. DNA to be packaged should be relatively free from contaminants. *Polyethylene glycol (PEG), which is contained in some ligase buffers, can inhibit packaging.* The volume of DNA added to each extract should be between 1 and 4 µl. To obtain the highest packaging efficiency [i.e., the number of plaque-forming units per microgram (pfu/µg) of DNA], package 1 µl of the ligation reaction and never more than 4 µl. Increased volume (i.e., >4 µl) yields more plaque-forming units per packaging reaction, but fewer plaque-forming units per microgram of DNA.

DNA that is digested with restriction enzymes and religated packages less efficiently (by a factor of 10–100) than uncut lambda DNA. For example, uncut wild-type lambda DNA packages with efficiencies exceeding 1×10^9 pfu/µg of vector when using a Gigapack III packaging extract. However, predigested vector, when ligated to a test insert, yields $\sim 5 \times 10^6$ – 1×10^7 recombinant plaques/µg of vector.

Gigapack III XL Packaging Extract

Gigapack III XL packaging extract is an *in vitro* packaging extract that preferentially size selects for extra large inserts, while maintaining the highest packaging efficiencies commercially available. This extract is specifically designed for use in generating genomic libraries. For example, a 20-kb insert will be packaged with a 95% higher efficiency than a 14-kb insert when using replacement vectors such as the Lambda FIX II vector. This simplifies library construction and reduces losses of precious genomic samples by eliminating the need to size select genomic inserts on sizing columns or sucrose gradients.

Packaging Instructions

For optimal packaging efficiency, package 1 µl of the ligation and never more than 4 µl. For further selection of large inserts, use the Gigapack III XL packaging extract, a size-selective packaging extract.

Preparing the Host Bacteria

Note *Prepare an overnight culture of the VCS257 strain (see the table in Preparation of Host Strains) prior to performing the protocol for the positive wild-type lambda DNA control (see Positive Wild-Type Lambda DNA Control for the Gigapack III Packaging Extract).*

1. Streak the bacterial glycerol stock onto the appropriate agar plates (see the table in *Preparation of Host Strains*). Incubate the plates overnight at 37°C.
2. Inoculate an appropriate medium, supplemented with 10 mM MgSO₄ and 0.2% (w/v) maltose, with a single colony.
3. Grow at 37°C, shaking for 4–6 hours (do not grow past an OD₆₀₀ of 1.0). Alternatively, grow overnight at 30°C, shaking at 200 rpm.

Note *The lower temperature keeps the bacteria from overgrowing, thus reducing the number of nonviable cells. Phage can adhere to nonviable cells resulting in a decreased titer.*

4. Pellet the bacteria at 500 × g for 10 minutes.
5. Gently resuspend the cells in half the original volume with sterile 10 mM MgSO₄.
6. Dilute the cells to an OD₆₀₀ of 0.5 with sterile 10 mM MgSO₄.

Note *The bacteria should be used immediately following dilution.*

Packaging Protocol

Note *Polyethylene glycol, which is present in some ligase buffers, can inhibit packaging.*

1. Remove the appropriate number of packaging extract vials from a –80°C freezer and place the extract vials on dry ice.
2. Quickly thaw the packaging extract by holding the tube between your fingers until the contents of the tube just begins to thaw.
3. Add the experimental DNA **immediately** (1–4 µl containing 0.1–1.0 µg of ligated DNA) to the packaging extract.
4. Stir the tube with a pipet tip to mix well. **Gentle** pipetting is allowable provided that air bubbles are not introduced.
5. Spin the tube quickly (for 3–5 seconds), if desired, to ensure that all contents are at the bottom of the tube.

6. Incubate the tube at room temperature (22°C) for 2 hours. **Do not exceed 2 hours.**

Note *The highest efficiency occurs between 90 minutes and 2 hours. Efficiency may drop dramatically during extended packaging times.*

7. Add 500 µl of SM buffer (see *Preparation of Media and Reagents*) to the tube.
8. Add 20 µl of chloroform and mix the contents of the tube gently.
9. Spin the tube briefly to sediment the debris.
10. The supernatant containing the phage is ready for titering. The supernatant may be stored at 4°C for up to 1 month.

Testing the Efficiency of the Gigapack III Packaging Extract with the Wild-Type Lambda Control DNA (Optional)

Use the following procedure to test the efficiency of the Gigapack III packaging extract with the λ cI857 *Sam7* wild-type lambda control DNA:

1. Thaw the frozen wild-type lambda control DNA on ice and gently mix after thawing.
2. Using 1 µl of the wild-type lambda control DNA (~0.2 µg), proceed with steps 1–10 in the *Packaging Protocol*.

Note *Because of the high titer achieved with the wild-type lambda control DNA, we recommend stopping the control packaging reaction with 1 ml of SM buffer. This should make the plaques easier to count.*

3. Prepare two consecutive 10^{-2} dilutions of the packaging reaction from step 10 in the *Packaging Protocol* in SM buffer. (The final dilution is 10^{-4} .)
4. Add 10 µl of the 10^{-4} dilution to 200 µl of the VCS257 host strain. (This strain is recommended for plating the wild-type lambda control DNA only.) Incubate at 37°C for 15 minutes. Add 3 ml of NZY top agar (see *Preparation of Media and Reagents*), melted and cooled to ~48°C, and quickly pour the dilution onto dry, pre-warmed NZY agar plates.
5. Incubate the plates for at least 12 hours at 37°C.
6. Count the plaques. Approximately 400 plaques should be obtained on the 10^{-4} dilution plate when the reaction is stopped with 1 ml of SM buffer.

TITERING THE PACKAGING REACTION

1. Streak the bacterial glycerol stock onto the appropriate agar plates (see the table in *Preparation of Host Strains*). Incubate the plates overnight at 37°C.
2. Inoculate an appropriate medium, supplemented with 10 mM MgSO₄ and 0.2% (w/v) maltose, with a single colony.
3. Grow at 37°C, shaking for 4–6 hours (do not grow past an OD₆₀₀ of 1.0). Alternatively, grow overnight at 30°C, shaking at 200 rpm.

Note *The lower temperature keeps the bacteria from overgrowing, thus reducing the number of nonviable cells. Phage can adhere to nonviable cells resulting in a decreased titer.*

4. Pellet the bacteria at 500 × g for 10 minutes.
5. Gently resuspend the cells in half the original volume with sterile 10 mM MgSO₄.
6. Dilute the cells to an OD₆₀₀ of 0.5 with sterile 10 mM MgSO₄.

Note *The bacteria should be used immediately following dilution.*

7. Prepare dilutions of the final packaged reaction in SM buffer. Add 1 µl of the final packaged reaction to 200 µl of host cells diluted in 10 mM MgSO₄ to an OD₆₀₀ of 0.5. If desired, also add 1 µl of a 1:10 dilution of the packaged reaction in SM buffer to 200 µl of host cells.
8. Incubate the phage and the bacteria at 37°C for 15 minutes to allow the phage to attach to the cells.
9. Add 3 ml of NZY top agar (48°C) and plate immediately on pre-warmed NZY agar plates.

- Count the plaques and determine the titer in plaque-forming units per milliliter (pfu/ml).

Both recombinant and nonrecombinant phage will grow on XL1-Blue MRA, but only recombinant phage will grow on XL1-Blue MRA (P2). Plaques should be visible after 8–12 hours of incubation at 37°C.

Note *Historically, the host strain LE392 has been used with this vector; however, E. coli restriction systems in this strain have a significant negative effect on the efficiency of DNA cloning and the ability to generate libraries representative of the gene population. The strains provided with this kit, XL1-Blue MRA and XL1-Blue MRA (P2), are McrA⁻, McrB⁻ and Mrr⁻; these modifications have been demonstrated to cause up to a 10-fold increase in the yield of recombinant phage containing methylated DNA. In addition, these strains have been further modified to enhance the stability of nonstandard DNA structures.*

AMPLIFYING THE LIBRARY

It is usually desirable to amplify libraries prepared in lambda vectors to make a large, stable quantity of a high-titer stock of the library. However, more than one round of amplification is not recommended, since slower growing clones may be significantly underrepresented.

Day 1

- Prepare the host strains as outlined in *Preparation of Host Strains*.

Day 2

- Dilute the cells to an OD₆₀₀ of 0.5 in 10 mM MgSO₄. Use 600 µl of cells at an OD₆₀₀ of 0.5/150-mm plate.
- Combine aliquots of the packaged mixture or library suspension containing ~5 × 10⁴ pfu of bacteriophage with 600 µl of host cells at an OD₆₀₀ of 0.5 in 14-ml BD Falcon polypropylene tubes. To amplify 1 × 10⁶ plaques, use a total of 20 aliquots (expect 5 × 10⁴ plaques/150-mm plate for each aliquot).

Note *Do not add more than 300 µl of phage/600 µl of cells.*

- Incubate the tubes containing the phage and host cells for 15 minutes at 37°C.
- Mix 6.5 ml of NZY top agar, melted and cooled to ~48°C, with each aliquot of infected bacteria and spread evenly onto a freshly poured 150-mm NZY bottom agar plate.

6. Incubate the plates at 37°C for 6–8 hours. Do not allow the plaques to get larger than 1–2 mm. On completion, the plaques should be touching.
7. Overlay the plates with ~8–10 ml of SM buffer. Store the plates at 4°C overnight (with *gentle* rocking if possible). This allows the phage to diffuse into the SM buffer.

Day 3

8. Recover the bacteriophage suspension from each plate and pool it into a sterile polypropylene container. Rinse the plates with an additional 2 ml of SM buffer and pool. Add chloroform to a 5% (v/v) final concentration. Mix well and incubate for 15 minutes at room temperature.
9. Remove the cell debris by centrifugation for 10 minutes at 500 × g.
10. Recover the supernatant and transfer it to a sterile polypropylene container. If the supernatant appears cloudy or has a high amount of cell debris, repeat steps 8 and 9. If the supernatant is clear, add chloroform to a 0.3% (v/v) final concentration and store at 4°C. We recommend storing aliquots of the amplified library in 7% (v/v) DMSO at –80°C.
11. Check the titer of the amplified library using host cells and serial dilutions of the library. (Assume ~10⁹–10¹¹ pfu/ml.)

PERFORMING PLAQUE LIFTS

1. Titer the library to determine the concentration (prepare fresh host cells to use in titrating and in screening).
2. Plate on large 150-mm agar plates (≥2-day-old) to 50,000 pfu/plate with 600 µl of host cells at an OD₆₀₀ of 0.5/plate and 6.5 ml of NZY top agar/plate. (Use 20 plates to screen 1 × 10⁶ plaques.)
3. Incubate the plates at 37°C for ~8 hours.
4. Chill the plates for 2 hours at 4°C to prevent the top agar from sticking to the nitrocellulose membrane.

Note Use forceps and wear gloves for the following steps.

5. Transfer the plaques onto a nitrocellulose membrane for 2 minutes. Use a needle to prick through the agar for orientation. (If desired, waterproof ink in a syringe needle may be used.)

If making duplicate nitrocellulose membranes, allow the second membrane to transfer for ~4 minutes.

Note *Pyrex® dishes are convenient for the following steps. All solutions should be at room temperature.*

- a. Denature the nitrocellulose membrane after lifting by submerging the membrane in a 1.5 M NaCl and 0.5 M NaOH denaturation solution for 2 minutes.

Note *If using charged nylon, wash with gloved fingertips to remove the excess top agar.*

- b. Neutralize the nitrocellulose membrane for 5 minutes by submerging the membrane in a 1.5 M NaCl and 0.5 M Tris-HCl (pH 8.0) neutralization solution.

- c. Rinse the nitrocellulose membrane for no more than 30 seconds by submerging the membrane in a 0.2 M Tris-HCl (pH 7.5) and 2× SSC buffer solution (see *Preparation of Media and Reagents*).

6. Blot briefly on Whatman® 3MM paper.
7. Crosslink the DNA to the membranes using the autocrosslink setting on the Stratalinker UV crosslinker (120,000 μJ of UV energy) for ~30 seconds. Alternatively, oven bake at 80°C for ~1.5–2 hours.
8. Store the stock agar plates of the transfers at 4°C to use after screening.

HYBRIDIZING AND SCREENING

Following the preparation of the membranes for hybridization, perform prehybridization, probe preparation, hybridization, and washes for either oligonucleotide probes or double-stranded probes and then expose the membranes to film as outlined in standard methodology texts.^{5,6} Following these procedures, perform secondary and tertiary screenings also as outlined in the standard methodology texts.^{5,6} After an isolate is obtained, refer to Sambrook *et al.*⁶ for suggested phage miniprep and maxiprep procedures.

RAPID RESTRICTION MAPPING

The insertion sites of the Lambda FIX II vector are flanked by T3 and T7 promoters, which permit the generation of end-specific hybridization probes. End-specific probes can be made once a recombinant clone containing an insert is isolated. In addition, the Lambda FIX II vector has unique *Not* I sites flanking the RNA promoters, which permits the excision from the lambda vector of insert DNA with the T3 and T7 promoter sequences as an intact fragment.

TROUBLESHOOTING

Observations	Suggestions
Packaging efficiency is too low	Storage conditions are inappropriate. Gigapack III packaging extract is very sensitive to slight variations in temperature. Store the packaging extracts at the bottom of a -80°C freezer and avoid transferring tubes from one freezer to another
	Do not allow the packaging extracts to thaw
	Avoid use of ligase buffers containing PEG, which can inhibit packaging
	The concentration of DNA may be too low. Ligate at DNA concentrations of $0.2\ \mu\text{g}/\mu\text{l}$ or greater and package between 1 and $4\ \mu\text{l}$ of the ligation reaction
Neither a bacterial lawn nor plaques is observed on the plate	The proteins contained within the packaging extract are diluted to suboptimal concentration when combined with an excessive volume of the ligation reaction mixture. Never package $>4\ \mu\text{l}$ of the ligation reaction
	Chloroform, which is added after packaging to prevent bacterial contamination, if present while titering, inhibits growth of the lawn of bacterial host cells. Be sure to spin down the chloroform completely prior to removing an aliquot for titering

PREPARATION OF MEDIA AND REAGENTS

<p>10× Fill-In Buffer 60 mM Tris-HCl (pH 7.5) 60 mM NaCl 60 mM MgCl₂ 0.5% gelatin 10 mM dithiothreitol (DTT)</p>	<p>LB Agar (per Liter) 10 g of NaCl 10 g of tryptone 5 g of yeast extract 20 g of agar Add deionized H₂O to a final volume of 1 liter Adjust pH to 7.0 with 5 N NaOH Autoclave Pour into petri dishes (~25 ml/100-mm plate)</p>
<p>LB Broth (per Liter) 10 g of NaCl 10 g of tryptone 5 g of yeast extract Add deionized H₂O to a final volume of 1 liter Adjust to pH 7.0 with 5 N NaOH Autoclave</p>	
<p>NZY Agar (per Liter) 5 g of NaCl 2 g of MgSO₄ · 7H₂O 5 g of yeast extract 10 g of NZ amine (casein hydrolysate) 15 g of agar Add deionized H₂O to a final volume of 1 liter Adjust the pH to 7.5 with NaOH Autoclave Pour into petri dishes (~80 ml/150-mm plate)</p>	<p>NZY Broth (per Liter) 5 g of NaCl 2 g of MgSO₄ · 7H₂O 5 g of yeast extract 10 g of NZ amine (casein hydrolysate) Add deionized H₂O to a final volume of 1 liter Adjust the pH to 7.5 with NaOH Autoclave</p>
<p>NZY Top Agar (per Liter) Prepare 1 liter of NZY broth Add 0.7% (w/v) agarose Autoclave</p>	<p>TE Buffer 10 mM Tris-HCl (pH 7.5) 1 mM EDTA</p>
<p>20× SSC Buffer (per Liter) 175.3 g of NaCl 88.2 g of sodium citrate 800.0 ml of deionized H₂O Adjust to pH 7.0 with a few drops of 10 N NaOH Add deionized H₂O to a final volume of 1 liter</p>	<p>SM Buffer (per Liter) 5.8 g of NaCl 2.0 g of MgSO₄ · 7H₂O 50.0 ml of 1 M Tris-HCl (pH 7.5) 5.0 ml of 2% (w/v) gelatin Add deionized H₂O to a final volume of 1 liter Autoclave</p>
<p>5× Transcription Buffer 200 mM Tris-HCl (pH 8.0) 40 mM MgCl₂ 10 mM spermidine 250 mM NaCl</p>	<p>10× Ligation Buffer 500 mM Tris-HCl (pH 7.5) 70 mM MgCl₂ 10 mM dithiothreitol (DTT)</p>

REFERENCES

1. Kohler, S. W., Provost, G. S., Kretz, P. L., Dyaico, M. J., Sorge, J. A. *et al.* (1990) *Nucleic Acids Res* 18(10):3007-13.
2. Kretz, P. L., Reid, C. H., Greener, A. and Short, J. M. (1989) *Nucleic Acids Res* 17(13):5409.
3. Kretz, P. L., Kohler, S. W. and Short, J. M. (1991) *J Bacteriol* 173(15):4707-16.
4. Kretz, P. L. and Short, J. M. (1989) *Strategies in Molecular Biology* 2(2):25-26.
5. Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G. *et al.* (1987). *Current Protocols in Molecular Biology*. John Wiley & Sons, New York.
6. Sambrook, J., Fritsch, E. F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

ENDNOTES

Pyrex® is a registered trademark of Corning Glass Works.
Whatman® is a registered trademark of Whatman Ltd.
BD Falcon is a trademark of BD Biosciences.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.