

Cloned *Pfu* DNA Polymerase AD

Catalog #600353, 600355, and 600357

600353-12, Revision A.01

MATERIALS PROVIDED

Materials provided	Quantity		
	Catalog #600353	Catalog #600355	Catalog #600357
Cloned <i>Pfu</i> DNA Polymerase AD	100 U	500 U	1000 U
10× Cloned <i>Pfu</i> Reaction Buffer AD	1 ml	2 × 1 ml	4 × 1 ml

Storage: Store at –20°C upon receipt.

INTRODUCTION

Pfu DNA polymerase AD* features a reformulated buffer system for increased economy with the same high performance as our original *Pfu* DNA polymerase, a proofreading DNA polymerase isolated from *Pyrococcus furiosus*.¹⁻³ This enzyme is an ideal choice for a variety of applications requiring high-fidelity DNA synthesis by the polymerase chain reaction (PCR), including cloning, gene expression, and site-directed mutagenesis. Successful PCR using *Pfu* DNA polymerase is readily performed requiring only slight modifications from PCR protocols optimized with *Taq* DNA polymerase.

OPTIMIZATION PARAMETERS (50 µL REACTION VOLUME)

Parameter	Genomic or Vector DNA Targets ≤2 kb	cDNA Targets ≤2 kb
Extension time	1–2 minutes per kb	2 minutes per kb
Cloned <i>Pfu</i> DNA polymerase AD	1.25–2.5 U	1.25–2.5 U
Input template	25–100 ng genomic DNA; 0.5–50 ng vector DNA	0.5–1 µl cDNA (from cDNA synthesis reaction)
Primers (each)	50–125 ng (0.1–0.5 µM)	50–125 ng (0.1–0.5 µM)
dNTP concentration	100–250 µM each dNTP	100–250 µM each dNTP
Mg ²⁺ concentration	2 mM (obtained in 1.0× dilution of provided 10× reaction buffer)	3 mM (supplement reaction with 1mM MgSO ₄)
Denaturing temperature	95°C	95°C
Extension temperature	72°C	68°C

PCR PROTOCOL

The reaction conditions given here are for amplification of a typical single-copy chromosomal target of ≤2 kb. See the *Optimization Parameters* section for guidelines on amplifying vector DNA or cDNA targets. The reaction conditions are for one reaction and must be adjusted for multiple samples. The final volume of each reaction is 50 µl.

1. Add the components *in order* into sterile thin-walled PCR tubes while mixing gently.

Reaction Mixture for a Typical Single-Copy Chromosomal Locus PCR Amplification (≤2 kb)

Component	Amount per reaction
Distilled water (dH ₂ O)	40.1 µl
10× Cloned <i>Pfu</i> Reaction Buffer AD ^a	5.0 µl
dNTP mix (25 mM each dNTP)	0.4 µl
DNA template (100 ng/µl)	1.0 µl
Primer #1 (100 ng/µl)	1.25 µl ^b
Primer #2 (100 ng/µl)	1.25 µl ^b
Cloned <i>Pfu</i> DNA Polymerase AD (2.5 U/µl)	1.0 µl ^c
Total reaction volume	50 µl

^a The 10× buffer provides a final 1× Mg²⁺ concentration of 2 mM. When amplifying cDNA, add Mg²⁺ to a final 1× concentration of 3 mM. (For example, Mg²⁺ concentration may be adjusted to 3 mM in the final 50-µl reaction volume by adding 2 µl of a PCR-grade 25 mM MgSO₄ solution and reducing the amount of dH₂O to 38.1 µl.)

^b The recommended primer concentration of 0.1–0.5 µM corresponds to 50–125 ng for a typical 18- to 25-mer oligonucleotide primer in a 50-µl reaction volume.

^c The amount of *Pfu* DNA polymerase AD varies depending on the length of the PCR target. The standard amount for PCR targets up to 2 kb in length is 1 µl (2.5 U).

- Perform PCR using optimized cycling conditions. Suggested cycling parameters for cloned *Pfu* DNA polymerase AD-based PCR are provided below. (Optimized cycling parameters are not necessarily transferable between thermal cyclers. Consult the instrument manufacturer's recommendations if further optimization of cycling parameters is required.) Analyze the PCR amplification products on a 0.7–1.0% (w/v) agarose gel.

PCR Cycling Program for Cloned *Pfu* DNA Polymerase AD

Single-Block Temperature Cyclers

Segment	Number of cycles	Temperature	Duration
1	1	94–98°C ^a	1 minute
2	25–30	94–98°C ^a	30 seconds
		Primer $T_m - 5^\circ\text{C}$	30 seconds
		72°C ^b	1–2 minute per kb
3	1	72°C ^b	10 minutes

^a Denaturing temperatures above 95°C are recommended only for GC-rich templates.

^b An extension temperature of 68°C is recommended for cDNA target amplification.

Stratagene RoboCycler Temperature Cycler

Segment	Number of cycles	Temperature	Duration
1	1	94–98°C ^a	1 minute
2	25–30	94–98°C ^a	1 minute
		Primer $T_m - 5^\circ\text{C}$	1 minute
		72°C ^b	1–2 minute per kb
3	1	72°C ^b	10 minutes

^a Denaturing temperatures above 95°C are recommended only for GC-rich templates.

^b An extension temperature of 68°C is recommended for cDNA target amplification.

TROUBLESHOOTING AND APPLICATION NOTES

- Low yield:** If PCR product yields are lower than expected, optimize the extension time, amount of DNA template (excess DNA can inhibit PCR), and amount of *Pfu* DNA polymerase AD. Optimize the PCR program denaturation time (typically 30–60 seconds at 95°C is sufficient; prolonged denaturation steps may damage the template DNA) and the annealing temperature. Extraneous salts contributed by primer or template DNA solutions may inhibit the PCR reaction. Use high-quality, gel purified primers and highly-purified template DNA.
- Multiple bands:** The annealing temperature may require optimization. Typically annealing temperatures will range between 55°C and 72°C. Try adding Perfect Match PCR Enhancer (Stratagene Catalog #600129) to improve specificity. Redesign primers.
- PCR adjuncts and cosolvents:** Including a cosolvent, such as 1–10% (v/v) DMSO or 5–20% glycerol, or an adjunct, such as 1.25–10% formamide or 0.01–1 U of Perfect Match PCR Enhancer (Stratagene Catalog #600129), in the PCR reaction may increase performance for some targets/cycling programs.
- Amplification of targets >2 kb:** For amplification of PCR targets >2 kb in length, we recommend using *PfuTurbo* DNA Polymerase AD (Catalog #600255, #600257, and #600259). This enzyme features equivalent fidelity with increased target-length capability compared to *Pfu* DNA polymerase AD.
- PCR cloning:** If generating PCR fragments for cloning applications, use the StrataClone Blunt PCR Cloning Kit (Catalog #240207) or another blunt PCR cloning strategy. *Pfu* DNA polymerase does **not** exhibit terminal deoxynucleotidyltransferase (TdT) activity, which is characterized by the addition of nontemplate-directed nucleotide(s) at the 3' end of PCR-generated fragments.

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

LIMITED LABEL LICENSE FOR CLONED *Pfu* DNA POLYMERASE PRODUCTS

This product is covered by the claims of one or more of the following U.S. Patents: 5,545,552; 5,556,772; 5,866,395; 5,948,663; 6,183,997; 6,444,428; 6,489,150; 6,734,293; 7,045,328, and patents pending. Purchase of this product conveys to the purchaser only the non-transferable right under these patents to use the product **for research use only by the purchaser**. No rights are granted to the purchaser hereunder to sell, modify for resale or otherwise transfer this product, either alone or as a component of another product, to any third party. Agilent reserves all other rights, and this product may not be used in any manner other than as provided herein. For information on obtaining a license to use this product for purposes other than research, please contact Agilent Technologies, Stratagene Product Division Business Development, 11011 North Torrey Pines Road, La Jolla, California 92037. Phone (858) 535-5400.

REFERENCES

- Cha, R. S. and Tilly, W. G. (1995). In PCR Primer, C. W. Dieffenbach and G. S. Dveksler (Eds.). Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York.
- Flaman, J. M., Frebourg, T., Moreau, V., Charbonnier, F., Martin, C. et al. (1994) Nucleic Acids Res 22(15):3259-60.
- Lundberg, K. S., Shoemaker, D. D., Adams, M. W., Short, J. M., Sorge, J. A. et al. (1991) Gene 108(1):1-6.

ENDNOTES

* U.S. Patent Nos. 7,045,328; 6,489,150; 5,948,663; 5,866,395; 5,545,552 and patents pending.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.