HaloPlex

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
The screenshot shows the SureDesign website interface. On the left, there is a banner with a blue DNA double helix graphic and the text "Collaborate to impact research around the world". On the right, there is a "Sign In" form highlighted with an orange border. The form includes a text input field for the email address, a password input field labeled "Password", a "Forgot Password?" link, a "Sign In" button, a "Register" button, and a language dropdown menu set to "English".

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The image displays two screenshots of the SureDesign web application interface, illustrating the steps to access the HaloPlex design tool.

Top Screenshot: The interface shows the 'Start Designing' section with various design tools. The 'HaloPlex' icon is highlighted with a yellow box. A yellow arrow points from the instruction text to the 'Show Advanced Options' checkbox, which is also highlighted with a yellow box.

Bottom Screenshot: The 'Create Designs' menu is open, showing a list of design tools. The 'HaloPlex' option is highlighted with a yellow box. A yellow arrow points from the instruction text to the 'Show Advanced Options' checkbox, which is also highlighted with a yellow box.

The interface includes a navigation bar with 'Home', 'Find Designs', and 'Create Designs' tabs. The 'Messages (0)' section on the right indicates no messages are present. The bottom section shows 'Designs: In Progress' and 'Designs: Recent and Favorites'.

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The screenshot shows the SureDesign web application interface. At the top, there's a navigation bar with 'SureDesign' logo and buttons for 'Home', 'Find Designs', and 'Create Designs'. Below this, a sub-navigation bar shows 'Create Designs' and 'HaloPlex' tabs. On the right, there's a 'Show Advanced Options' checkbox (checked), a shopping cart icon with 'Cart (0)', and links for 'Logout', 'Settings', and 'Help - Create Designs'. The main content area is titled '1. Choose Wizard' and '2. Choose Design or Probegroup'. Under '1. Choose Wizard', there are two radio buttons: 'Standard' and 'Advanced'. The 'Advanced' option is selected and highlighted with a yellow box. Below it, a description states: 'The advanced wizard allows you to create a design composed of multiple probegroups. You can select new probes for defined target regions; or you can include probegroups that already exist in your account.' Under '2. Choose Design or Probegroup', there are two radio buttons: 'Create Design' and 'Create Probegroup'. The 'Create Design' option is selected and highlighted with a yellow box. Below it, a description states: 'A design is a set of one or more probegroups that is manufactured as a single library for use in a SureSelect DNA or HaloPlex protocol.' At the bottom right, there is a blue 'Continue' button highlighted with a yellow box.

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The screenshot shows the SureDesign web application interface. The top navigation bar includes 'SureDesign' and 'Help - Add/Review Content'. The main content area is titled 'HaloPlex (Advanced)' and 'Add/Review Content'. On the left, a sidebar shows a progress bar with 'Define Design' (checked) and 'Add/Review Content' (active). Below this, options include 'Tile Genes or Regions', 'Select Existing Probes', and 'Finalize'. The main panel has a 'Select Content Addition Method' section with two radio buttons: 'Design new probes by tiling genes or regions' (selected and highlighted with a yellow box) and 'Select probes from an existing Design or Probegroup'. Below the second option is a 'Select all probes' dropdown. At the bottom left, a 'HaloPlex Design' summary box displays details: Name: HaloPlexdemo, Species: H. sapiens (hg19), Category: HaloPlex, Target Regions: # Regions: NA, Size: NA, Sequenced Regions: # Amplicons: NA, Size: NA, Price Tier: NA, Coverage: NA. It includes buttons for 'UCSC View', 'BAM', and 'Download'. At the bottom right, there are 'Cancel', 'Back', and 'Next' buttons, with the 'Next' button highlighted by a yellow box.

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SureDesign Help - Tile Genes or Regions

HaloPlex (Advanced) **Tile Genes or Regions**

Define Design ☒
Add/Review Content ☒
▶ **Tile Genes or Regions**
Select Existing Probes
Finalize

HaloPlex Design
Name: HaloPlexdemo
Species: H. sapiens (hg19)
Category: HaloPlex

Target Regions
Regions: NA
Size: NA

Sequenced Regions
Amplicons: NA
Size: NA
Price Tier: NA
Coverage: NA

UCSC View BAM
Download

Target Summary
• 14 Target IDs resolved to 10 targets comprising 174 regions.
• 4 Target IDs were not found.

Target Details [View targets in UCSC](#)

Target ID	# Regions	Base Pairs	Position
ABL	Not found		
BCL2ALPHA	Not found		
BCL2BETA	Not found		
CRK-II	Not found		
AKT2	16	2104	chr19:40739769-40771184
APC	19	9221	chr5:112043405-112179833
BCL3	9	1605	chr19:45252038-45262882
BCR	24	4421	chr22:23523138-23657719
BRCA1	26	6363	chr17:41197685-41277212
BRCA2	28	10983	chr13:32890588-32972917
CBL	17	3135	chr11:119077118-119177081
CCND1	5	988	chr11:69456072-69466060
CDK4	8	1156	chr12:58142298-58145510
CSE1R	22	3443	chr5:149433622-149466000

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SureDesign Help - Tile Genes or Regions

HaloPlex (Advanced) Tile Genes or Regions

Define Design ✓
Add/Review Content ✓
▶ Tile Genes or Regions
Select Existing Probes
Finalize

* Probegroup Name: HaloPlexdemo_1

Selection Parameters

For the Stringency parameter, Agilent recommends the Balanced option. Selecting the more stringent option (Maximize Specificity) may decrease design coverage but minimize the capture of off-target genomic regions. Selecting the less stringent option (Maximize Coverage) may increase design coverage but can significantly increase the capture of off-target genomic regions, particularly if your target regions include UTRs, introns, or other non-coding regions. Consequently, you will need to increase sequencing in order to achieve your desired read-depth. Please read the help pages for more details.

Stringency: Balanced

Optimize for FFPE Samples: ☐

Reset

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Balanced
Maximize Specificity
Balanced
Maximize Coverage

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Tile Genes or Regions
Select Existing Probes
Finalize

Select Content Addition Method

☐ Design new probes by tiling genes or regions
☒ Select probes from an existing Design or Probegroup

Select all probes
Select all probes
Select probes for specific genes or regions
Select probes by Amplicon ID

HaloPlex Design
Name: HaloPlexdemo
Species: H. sapiens (hg19)
Category: HaloPlex

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SureDesign Help - Select Probe Source

HaloPlex (Advanced) Select Existing Probes Select Probe Source

Define Design ☒ Add/Review Content Tile Genes or Regions Select Existing Probes Finalize

HaloPlex Design
 Name: HaloPexdemo
 Species: H. sapiens (hg19)
 Category: HaloPlex
Target Regions
 # Regions: 201
 Size: 40.539 kbp
Sequenced Regions
 # Amplicons: 2227
 Size: 118.754 kbp
 Price Tier: Tier 1
 Coverage: 99.24 %
 UCSC View BAM Download

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HaloPlex (Advanced) Select Existing Probes Select probes for specific genes or regions

Define Design ✓
Add/Review Content ✓
Tile Genes or Regions
▶ Select Existing Probes
Finalize

HaloPlex Design
Name: HaloPexdemo
Species: H. sapiens (hg19)
Category: HaloPlex

Target Regions
Regions: 201
Size: 40,539 kbp

Sequenced Regions
Amplicons: 2227
Size: 118,754 kbp
Price Tier: Tier 1
Coverage: 99.24 %
UCSC View BAM
Download

Targets:
Enter gene IDs, gene symbols, or accessions.
SureDesign will map these targets to genomic locations
using the parameters below, then select probes from the
defined probe source that match the locations.
SMAD4
NM_005359
ENST00000342988
CCDS11950
4089

Enter UCSC browser or BED coordinates
chr18:48573407-48573675
chr18 48573406 48573675

Enter coordinates with identifiers
chr18:48575655-48575704 MyGene1
chr18 48577713 48577795 MyGene2 # exon 1
chr18 48578994 48579032 MyGene2 # exon 2

Enter cytoband
7p22.1

Upload
Example
Clear

*** Databases**
☒ RefSeq
☒ Ensembl
☒ CCDS
☒ Gencode
☒ VEGA
☒ SNP
☒ CytoBand

Regions of Interest
☒ Coding Exons
☐ Coding Exons + UTRs ☒ 5' UTR ☒ 3' UTR
☐ Entire Transcribed Region

Include Flanking Bases:
3': 10 bp 5': 10 bp
☐ Allow Synonyms

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HaloPlex (Advanced) Select Existing Probes Select probes for specific genes or regions

Define Design ✓
Add/Review Content ✓
Tile Genes or Regions
Select Existing Probes
Finalize

Probegroup Name HaloPexdemo_1

Target Summary

- 3 Target IDs resolved to 3 targets comprising 83 regions.
- 0 Target IDs were not found.

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SureDesign Help - Select probes for specific genes or regions

HaloPlex (Advanced) Select Existing Probes **Select probes for specific genes or regions**

Define Design ☒
 Add/Review Content ☒
 Tile Genes or Regions
Select Existing Probes
 Finalize

HaloPlex Design
 Name: HaloPexdemo
 Species: H. sapiens (hg19)
 Category: HaloPlex

Target Regions
 # Regions: 201
 Size: 40,539 kbp

Sequenced Regions
 # Amplicons: 2227
 Size: 118,754 kbp
 Price Tier: Tier 1
 Coverage: 99.24 %
 UCSC View BAM
 Download

Targets:
 # Enter gene IDs, gene symbols, or accessions.
 # SureDesign will map these targets to genomic locations
 # using the parameters below, then select probes from the
 # defined probe source that match the locations.
 SMAD4
 NM_005359
 ENST00000342988
 CCDS11950
 4089
 # Enter UCSC browser or BED coordinates
 chr18:48573407-48573675
 chr18 48573406 48573675
 # Enter coordinates with identifiers
 chr18:48575655-48575704 MyGene1
 chr18 48577713 48577795 MyGene2 # exon 1
 chr18 48578994 48579032 MyGene2 # exon 2
 # Enter cytoband
 7p22.1

Upload
 Example
 Clear

*** Databases**
☒ RefSeq
☒ Ensembl
☒ CCDS
☒ Gencode
☒ VEGA
☒ SNP
☒ CytoBand

Regions of Interest
☒ Coding Exons
☐ Coding Exons + UTRs ☒ 5' UTR ☒ 3' UTR
☐ Entire Transcribed Region
 Include Flanking Bases:
 3': 10 bp 5': 10 bp
☐ Allow Synonyms

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SureDesign Help - Select probes for specific genes or regions

HaloPlex (Advanced) Select Existing Probes **Select probes for specific genes or regions**

Define Design ☒
 Add/Review Content ☒
 Tile Genes or Regions
Select Existing Probes
 Finalize

Probegroup Name

Target Summary

- 40 Target IDs resolved to 40 targets comprising 777 regions.
- 0 Target IDs were not found.

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SureDesign Help - Add/Review Content

HaloPlex (Advanced) Add/Review Content

Define Design

▶ Add/Review Content

Tile Genes or Regions

Select Existing Probes

Finalize

Select Content Addition Method

☐ Design new probes by tiling genes or regions

☒ Select probes from an existing Design or Probegroup

Select probes by Amplicon ID

- Or -

☐ Finalize this design

Probegroup Summary

Name	# Targets	Target Size	# Amplicons	Actions
HaloPlexdemo_1	Processing			

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SureDesign Help - Add/Review Content

HaloPlex (Advanced) Add/Review Content

Define Design ☒ Add/Review Content

Tile Genes or Regions

Select Existing Probes

Finalize

Select Content Addition Method

☒ Design new probes by tiling genes or regions

☐ Select probes from an existing Design or Probegroup

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
Note:

Tier 1 designs must contain up to 15K amplicons.

Tier 1 Plus designs must contain 15K to 20K amplicons.

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<Design ID>_ Report.txt
<ProbeGroup名>_ Report.txt

	A	B	C	D	E	F
1	# File Summary					
2	File Type: HaloPlex Design Report					
3	Created By: HaloPlex Advanced Design Wizard					

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10	Design Name: HaloPlexdemo					
11	Design ID:					
12	Design Category: HaloPlex					
13	Species: H. sapiens (UCSC hg19, GRCh37, February 2009)					
14	Platform: Illumina					
15	Read Length: 150 bp					
16	# Probegroup Summary					
17	Number of Probegroups: 2					
18	Probegroup 1: HaloPlexdemo_1					
19	Probegroup 2: HaloPlexdemo_2					

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<Design ID>_ targets.txt
<ProbeGroup名>_ targets.txt

	A	B	C	D	E	F
1	# File Type: Input Targets File					
2	# Timestamp:					
3	# User:					
4	# Created By: HaloPlex Advanced Design Wizard					

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
	A	B	C	D	E	F	G	H
1	# File Summary							
2	File Type: HaloPlex Design Report							
3	Created By: HaloPlex Advanced Design Wizard							
4	User:							
5	Workgroup:							
6	Folder:							
7	Timestamp:							
8								
9	# Design Summary							
10	Design Name: HaloPlexdemo							
11	Design ID: 04818-1672989845							
12	Design Category: HaloPlex							
13	Species: H. sapiens (UCSC hg19, GRCh37, February 2009)							
14	Platform: Illumina							
15	Read Length: 150 bp							
16								
17	# Probegroup Summary							
18	Number of Probegroups: 2							
19	Probegroup 1 : HaloPlexdemo_2							
20	Probegroup 2 : HaloPlexdemo_1							
21								
22	# Target Summary							
23	Target Region Size: 9.653 kbp							
24								
25	# Amplicon Summary							
26	Total Amplicons: 996							
27	Total Target Bases Analyzable: 9.375 kbp							
28	Total Sequenceable Design Size: 34.674 kbp							
29	Target Coverage: 97.12 %							
30	Recommended Minimum Sequencing per Sample: 6.934 Mbp							
31	Pricing: Illumina Tier 1 (Probe Region Size = 1 - 500 kbp; up to 15K probes)							

	A	B	C	D	E	F	G	H	I
34	# Target and Probe Details								
35	# TargetID: The identifier entered in the Targets list.								
36	# Interval: The genomic interval of the target.								
37	# Regions: The number of regions within this target.								
38	# Size: The total size (in base pairs) of the regions.								
39	# Database(s): The databases in which this target was found.								
40	# Coverage: Bases overlapped by probes (extended +/- 100 bp) to represent likely capture.								
41	# High Coverage: Number of regions where analyzable amplicon overlap >= 90%.								
42	# Low Coverage: Number of regions where analyzable amplicon overlap < 90								
43									
44	TargetID	Interval	Regions	Size	Databases	Coverage	HighCover	LowCoverage	
45	IL10	chr1:2069	5	637	CCDS, Ens	100	5	0	
46	IL11	chr19:558	5	700	CCDS, Ens	99.86	5	0	
47	IL13	chr5:1319	4	560	CCDS, Ens	100	4	0	
48	IL1A	chr2:1135	6	936	CCDS, Ens	100	6	0	
49	IL1B	chr2:1135	6	930	CCDS, Ens	100	6	0	
50	IL2	chr4:1233	4	542	CCDS, Ens	100	4	0	
51	IL3	chr5:1313	5	559	CCDS, Ens	100	5	0	
52	IL4	chr5:1320	5	663	CCDS, Ens	100	5	0	
53	IL5	chr5:1318	5	547	CCDS, Ens	100	5	0	
54	IL6	chr7:2276	4	1007	CCDS, Ens	100	4	0	
55	IL7	chr8:7964	7	743	CCDS, Ens	100	7	0	
56	IL8	chr4:7460	4	384	CCDS, Ens	100	4	0	
57	IL9	chr5:1352	5	535	CCDS, Ens	100	5	0	
58	chr1:1340	chr1:1340	1	105	NA	100	1	0	
59	chr1:3245	chr1:3245	1	95	NA	100	1	0	
60	chr1:6956	chr1:6956	1	115	NA	100	1	0	
61	chr1:6647	chr1:6647	1	42	NA	100	1	0	
62	chr1:7214	chr1:7214	1	339	NA	18.29	0	1	
63	chr1:3246	chr1:3246	1	95	NA	100	1	0	
64	chr1:3247	chr1:3247	1	119	NA	100	1	0	

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- Target Regions
- Amplicons
- Sequenceable Regions
- Covered Regions
- Missed Regions
- Analyzable Regions

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SureDesign

HaloPlex (Advanced) **Add/Review Content**

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☐ Select probes from an existing Design or Probe Set

Select all probes

- Or -

☐ Finalize this design

Probe group Summary

Name	# Targets
HaloPlexdemo_2	65
HaloPlexdemo_1	7

HaloPlex Design

Name: HaloPlexdemo_2

Species: H. sapiens (hg19)

Category: HaloPlex

Target Regions

Regions: 72

Size: 9.653 kbp

Sequenced Regions

Amplicons: 996

Size: 34.674 kbp

Price Tier: Tier 1

Coverage: 97.12 %

UCSC View **BAM**

Download

UCSC Genome Browser on Human (GRCh37/hg19)

move <<< << < > >> >>> zoom in 1.5x 3x 10x base zoom out 1.5x 3x 10x 100x

multi-region chr1:206,941,970-206,945,790 3,821 bp. gene, chromosome range, search terms, help pages, see example go examples

chr1 (q32.1) 1p31.1 1q12 q41 43 44

Scale chr1: 206,942,500| 206,943,000| 206,943,500| 206,944,000| 206,944,500| 206,945,000| 206,945,500| hg19

Target Regions

Covered

Missed

Agilent HaloPlex - HaloPlexdemo - Target regions of interest given as input to probe selection

Agilent HaloPlex - HaloPlexdemo - Genomic regions expected to be amplified

Agilent HaloPlex - HaloPlexdemo - Target regions of interest that will not be amplified

Agilent HaloPlex - HaloPlexdemo - Amplicons

Reference Assembly Fix Patch Sequence Alignments

Reference Assembly Alternate Haplotype Sequence Alignments

UCSC Genes (RefSeq, GenBank, CCDS, Rfam, tRNAs & Comparative Genomics)

IL10 IL11

RefSeq Curated

RefSeq genes from NCBI

Sequences

SNPs

Publications: Sequences in Scientific Articles

OMIM Gene Phenotypes - Dark Green Can Be Disease-causing

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Select Content Addition Method

☒ Design new probes by tiling genes or regions
☐ Select probes from an existing Design or Probegroup

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SureDesign HaloPlex Design Report

General Information

Design Information

Design Name: HaloPlexdemo
Design ID: [redacted]
Design Category: HaloPlex
Species: H. sapiens (UCSC hg19, GRCh37, February 2009)
Platform: Illumina
Read Length: 150 bp

Probegroup Summary
Number of Probegroups: 2
Probegroup 1 : HaloPlexdemo_1
Probegroup 2 : HaloPlexdemo_2

Target Summary
Target Region Size : 9.653 kbp

Amplicon Summary
Total Amplicons: 996
Total Target Bases Analyzable: 9.375 kbp
Total Sequenceable Design Size: 34.674 kbp
Target Coverage: 97.12 %
Recommended Minimum Sequencing per Sample: 6.934 Mbp
Pricing: Illumina Tier 1 (Probe Region Size = 1 - 500 kbp; up to 15K probes)
P/N: G9901C/G9901B

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<Design ID>_ Report.txt
<ProbeGroup名>_ Report.txt

Report.txt	
1	# File Summary
2	File Type: HaloPlex Design Report
3	Created By: HaloPlex Advanced Design Wizard
4	Usage:
10	Design Name: HaloPlexdemo
11	Design ID:
12	Design Category: HaloPlex
13	Species: H. sapiens (UCSC hg19, GRCh37, February 2009)
14	Platform: Illumina
15	Read Length: 150 bp
16	
17	# Probegroup Summary
18	Number of Probegroups: 2
19	Probegroup 1 : HaloPlexdemo_1
20	Probegroup 2 : HaloPlexdemo_2
21	

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<Design ID>_ targets.txt
<ProbeGroup名>_ targets.txt

Targets.txt	
1	# File Type: Input Targets File
2	# Timestamp:
3	# User:
4	# Created By: HaloPlex Advanced Design Wizard
10	chr1:324500-324594
11	chr1:324682-324776
12	chr1:324779-324897
13	chr1:664795-664836
14	chr1:721405-721743
15	IL1A
16	IL1B
17	IL2
18	IL3
19	IL4
20	IL5
21	IL6

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


The screenshot shows the 'Order' tab in the SureDesign HaloPlex interface. The form includes fields for Design Name, Design ID, Quantity (set to 1), Sequencing Platform (Illumina), Automation Kit (unchecked), Sample size (set to 48), and Design Part Number (G9901C). A callout box with three bullet points explains the settings: Quantity (number of kits), Automation Kit (select if used), and Sample Size (number of samples per kit, with an example of 48 samples). At the bottom are 'Add to Cart', 'Request a Quote', and 'Cancel' buttons. A note at the bottom left states that HaloPlex designs are not compatible with HaloPlex^{HS}.

SureDesign Help - Order

HaloPlex **Order**

Design Name: [blurred]
Design ID: [blurred]
★ Quantity: 1
Sequencing Platform: Illumina
Automation Kit: ☐
Sample size: 48
Design Part Number: G9901C

Note: HaloPlex designs are not compatible with HaloPlex^{HS}. Please re-design this library as a HaloPlex^{HS} library if you wish to upgrade your platform.



Workgroup Name: [blurred]
Company: [blurred]
Company: [blurred]
You are currently logged in as [blurred]
SureDesign is the Agilent purchase portal for registered users.
To check-out your order, click the 'Add to Cart' button. Your order will be transferred to the Agilent Genomics page. You must be logged in to see your discounted prices.

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SureDesign: Request for quote

We have received your request-for-quote for the following HaloPlex kit:

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差出人 ; Agilent SureDesign [suredesign_noreply@agilent.com]

タイトル ; SureDesign: Request for quote for HaloPlex kit received

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