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
The screenshot shows the SureDesign website interface. On the left, there is a banner with a blue DNA double helix graphic and the text "Collaborate to impact research around the world". On the right, there is a "Sign In" form. The form includes a text input field for the email address, a password input field labeled "Password", a "Forgot Password?" link, a "Sign In" button, a "Register" button, and a language dropdown menu set to "English". The entire sign-in form is highlighted with a yellow border.

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The screenshot displays the SureDesign web application interface. At the top, the 'SureDesign' logo is on the left, and navigation links for 'Home', 'Find Designs', and 'Create Designs' are in the center. On the right, there are links for 'Cart (0)', 'Logout', 'Settings', and 'Help - Home'. Below the navigation bar, the 'Start Designing' section features four icons: 'SureSelect DNA' (highlighted with a yellow box), 'SureSelect RNA', 'SureSelect All-In-One', and 'OneSeq'. A dropdown menu is open under 'Create Designs', showing 'SureSelect DNA' (highlighted with a yellow box), 'SureSelect RNA', 'OneSeq', 'HaloPlex', 'CGH', 'CustomFISH', 'SureGuide', and 'SurePrint OLS'. To the right of the 'Start Designing' section, there is a 'Messages (0)' box. Below the 'Start Designing' section, there is a 'Designs: In Progress' section. The main content area shows the 'Create Designs' wizard for 'SureSelect DNA'. The first step is '1. Choose Wizard', with 'Advanced' selected (highlighted with a yellow box). The second step is '2. Choose Design or Probegroup', with 'Create Design' selected (highlighted with a yellow box). The wizard provides detailed instructions for each step.

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SureDesign

Help - Define Design

SureSelect DNA (Advanced) Define Design

Define Design

Add/Review Content

Select Optimized Probes

Tile Genes or Regions

Upload Probes

Select Existing Probes

Finalize

* Design Name:

* Species: Select ⓘ

Build:

* Create In: Select

Hybridization: SureSelect XT HS2 / XT HS / XT LI / QXT

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SureDesign Help - Select Optimized Probes

SureSelect DNA (Advanced) Add/Review Content **Select Optimized Probes**

Define Design ☒ Add/Review Content ☒ **Select Optimized Probes** Tile Genes or Regions Upload Probes Select Existing Probes Finalize

SureSelect DNA Design
Name: Oncogene
Species: H. sapiens (hg19)
Category: SureSelect DNA
Hybridization: 90 Minutes
Probes
Probes: NA
Size: NA
Price Tier: [i](#) NA

*** Targets:**
Enter gene IDs, gene symbols, or accessions.
SureDesign will map these identifiers to genomic locations
using the parameters below, then select probes from the
reference genome that match the locations.
SMAD4
NM_005359
ENST00000342988
CCDS11950
4089
Enter UCSC browser or BED coordinates
chr18:48573407-48573675
chr18 48573406 48573675
Enter coordinates with identifiers
chr18:48575655-48575704 MyGene1
chr18 48577713 48577795 MyGene2 # exon 1
chr18 48578994 48579032 MyGene2 # exon 2
Enter cytoband
7p22.1

*** Databases**
☒ RefSeq
☒ Ensembl
☒ CCDS
☒ Gencode
☒ VEGA
☒ SNP
☒ CytoBand

Regions of Interest
☒ Coding Exons
☐ Coding Exons + UTRs ☒ 5' UTR ☒ 3' UTR
☐ Entire Transcribed Region
Include Flanking Bases:
3' : 10 bp 5' : 10 bp
☐ Allow Synonyms

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Download

Target Summary

13 Target IDs resolved to 13 targets comprising 206 regions.

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SureDesign Help - Select Optimized Probes

SureSelect DNA (Advanced) > Add/Review Content > **Select Optimized Probes**

Define Design ☒
Add/Review Content ☒
Select Optimized Probes
Tile Genes or Regions
Upload Probes
Select Existing Probes
Finalize

SureSelect DNA Design
Name: Oncogene
Species: H. sapiens (hg19)
Category: SureSelect DNA
Hybridization: 90 Minutes

Probes
Probes: NA
Size: NA
Price Tier: NA

Target Summary

- 13 Target IDs resolved to 13 targets comprising 206 regions.
- 0 Target IDs were not found.

Target Details View targets in UCSC

| Target ID | # Regions | Base Pairs | Position |
|-----------|-----------|------------|--------------------------|
| ARHGAP9 | 21 | 3013 | chr12:57866290-57882414 |
| C9ORF24 | 7 | 1124 | chr9:34379064-34397641 |
| CASK | 29 | 3554 | chrX:41379663-41782251 |
| COL27A1 | 62 | 6913 | chr9:116918221-11702985 |
| CWH43 | 16 | 2424 | chr4:48988437-49063917 |
| ELL | 13 | 2168 | chr19:18555552-18632875 |
| FFAR3 | 1 | 1061 | chr19:35849783-35850843 |
| GART | 22 | 3559 | chr21:34876421-34911631 |
| IL1RL1 | 11 | 2078 | chr2:102954715-102968391 |
| IL23R | 11 | 2193 | chr1:67633794-67724821 |
| KRTAP5-7 | 1 | 518 | chr11:71238337-71238854 |
| NKX2-6 | 2 | 946 | chr8:23559954-23564121 |
| OCIAD1 | 10 | 972 | chr4:48833529-48862789 |

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Optimized performance XT HS / XT HS2 / XT LI / QXT (Recommended)

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
The screenshot shows the SureDesign web interface for 'Add/Review Content'. The left sidebar lists steps: Define Design (checked), Add/Review Content (active), Select Optimized Probes, Tile Genes or Regions, Upload Probes, Select Existing Probes, and Finalize. The main area shows 'Select Content Addition Method' with four radio button options: 'Select Optimized Probes', 'Design new probes by tiling genes or regions' (selected and highlighted with an orange box), 'Upload probes', and 'Select probes from an existing Design or Probegroup'. Below these is a 'Select all probes' dropdown. At the bottom right, the 'Next' button is highlighted with an orange box. The bottom left shows 'UCSC View' and 'Download' buttons. The top right has a 'Help - Add/Review Content' link.

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Select Content Addition Method

☐ Select Optimized Probes

☒ Design new probes by tiling genes or regions

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SureDesign Help - Tile Genes or Regions

SureSelect DNA (Advanced) Add/Review Content **Tile Genes or Regions**

Define Design ☒
 Add/Review Content ☒
 Select Optimized Probes
Tile Genes or Regions
 Upload Probes
 Select Existing Probes
 Finalize

SureSelect DNA Design
 Name: Oncogene
 Species: H. sapiens (hg19)
 Category: SureSelect DNA
 Hybridization: 90 Minutes

Probes
 # Probes: NA
 Size: NA
 Price Tier: [i](#) NA

UCSC View Download

*** Targets:**
 # Enter gene IDs, gene symbols, or accessions.
 # SureDesign will map these identifiers to genomic locations
 # using the parameters below, then select probes from the
 # reference genome that match the locations.
 SMAD4
 NM_005359
 ENST00000342988
 CCDS11950
 4089
 # Enter UCSC browser or BED coordinates
 chr18:48573407-48573675
 chr18 48573406 48573675
 # Enter coordinates with identifiers
 chr18:48575655-48575704 MyGene1
 chr18 48577713 48577795 MyGene2 # exon 1
 chr18 48578994 48579032 MyGene2 # exon 2
 # Enter cytoband
 7p22.1

Upload
 Example

Clear

*** Databases**
☒ RefSeq
☒ Ensembl
☒ CCDS
☒ Gencode
☒ VEGA
☒ SNP
☒ CytoBand

Regions of Interest
☒ Coding Exons
☐ Coding Exons + UTRs ☒ 5' UTR ☒ 3' UTR
☐ Entire Transcribed Region

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Optimized performance XT / XT2 (Reccomended)

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Optimized performance XT HS / XT HS2 / XT LI / QXT (Reccomended)

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The screenshot shows the 'SureDesign' web interface, specifically the 'Add/Review Content' step. The left sidebar shows a progress bar with 'Define Design' completed and 'Add/Review Content' active. Under 'Add/Review Content', the options are 'Select Optimized Probes', 'Tile Genes or Regions', 'Upload Probes', 'Select Existing Probes', and 'Finalize'. The 'Upload Probes' option is selected. The main area shows 'Select Content Addition Method' with three radio buttons: 'Select Optimized Probes', 'Design new probes by tiling genes or regions', and 'Upload probes' (which is selected and highlighted with an orange box). Below these is a dropdown menu labeled 'Select all probes'. At the bottom, there are buttons for 'UCSC View', 'Download', 'Cancel', 'Back', and 'Next' (which is highlighted with an orange box).

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Define Design ☒ Add/Review Content

Select Optimized Probes
Tile Genes or Regions
Upload Probes
Select Existing Probes
Finalize

Select Content Addition Method

☐ Select Optimized Probes
☐ Design new probes by tiling genes or regions
☒ Upload probes
☐ Select probes from an existing Design or Probegroup

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SureDesign Help - Select Probe Source

SureSelect DNA (Advanced) Select Existing Probes Select Probe Source

Define Design ☒ Add/Review Content ☒ Select Optimized Probes Tile Genes or Regions Upload Probes **Select Existing Probes** Finalize

☒ Design ☐ Probegroup

* Workspace: Agilent Catalog My Workgroup **Agilent Catalog**

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SureDesign Help - Select probes for specific genes or regions

SureSelect DNA (Advanced) Select Existing Probes **Select probes for specific genes or regions**

Define Design ✓

Add/Review Content ✓

Select Optimized Probes

Tile Genes or Regions

Upload Probes

Select Existing Probes

Finalize

SureSelect DNA Design

Name: Oncogene

Species: H. sapiens (hg19)

Category: SureSelect DNA

Hybridization: 90 Minutes

Probes

Probes: NA

Size: NA

Price Tier: [i](#) NA

UCSC View Download

*** Targets:**

Enter gene IDs, gene symbols, or accessions.
SureDesign will map these targets to genomic locations
using the parameters below, then select probes from the
defined probe source that match the locations.

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
SureDesign Help - Add/Review Content

SureSelect DNA (Advanced) Add/Review Content

Define Design

► Add/Review Content

- Select Optimized Probes
- Tile Genes or Regions
- Upload Probes
- Select Existing Probes
- Finalize

SureSelect DNA Design

Name: Oncogene
Species: H. sapiens (hg19)
Category: SureSelect DNA
Hybridization: 90 Minutes

Probes

Probes: NA
Size: NA
Price Tier: NA

UCSC View Download

Select Content Addition Method

☒ Select Optimized Probes

☐ Design new probes by tiling genes or regions

☐ Upload probes

☐ Select probes from an existing Design or Probegroup

Select all probes

- Or -

☐ Finalize this design

Probegroup Summary


| Name | # Targets | Target Size | # Probes | Probe Size | Replicate Count | Actions |
|------------|------------|-------------|----------|------------|-----------------|---------|
| Oncogene_1 | Processing | | | | | |

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| # Target Summary | | | | | | | |
|---|------------|---------|-------|------------|----------|-----------|----------|
| 90 Target IDs resolved to 90 targets comprising 1804 regions. | | | | | | | |
| 0 Target IDs were not found. | | | | | | | |
| Region Size: 299.272 kbp | | | | | | | |
| # Probe Summary | | | | | | | |
| Total Probes: 3120 | | | | | | | |
| Total Probes Size: 157.069 kbp | | | | | | | |
| Recommended Minimum Sequencing per Sample: 31.413 Mbp | | | | | | | |
| # Target Parameters | | | | | | | |
| Databases: RefSeq, Ensembl, CCDS, Gencode, VEGA, SNP, CytoBand | | | | | | | |
| Region: Coding Exons | | | | | | | |
| Region Extension: 10 bases from 3' end and 10 bases from 5' end. | | | | | | | |
| Allow Synonyms: No | | | | | | | |
| # Target and Probe Details | | | | | | | |
| # TargetID: The identifier entered in the Targets list. | | | | | | | |
| # Interval: The genomic interval of the target. | | | | | | | |
| # Regions: The number of regions within this target. | | | | | | | |
| # Size: The total size (in base pairs) of the regions. | | | | | | | |
| # Database(s): The databases in which this target was found. | | | | | | | |
| # Coverage: Bases overlapped by probes (extended +/- 100 bp) to represent likely capture. | | | | | | | |
| # HighCoverage: Number of regions with overlap >= 90%. | | | | | | | |
| # LowCoverage: Number of regions with overlap < 90% | | | | | | | |
| TargetID | Interval | Regions | Size | Databases | Coverage | HighCover | LowCover |
| A1BG | chr19:5885 | 8 | 1860 | Gencode, f | 0 | 0 | 8 |
| A1CF | chr10:5256 | 12 | 2159 | Gencode, f | 0 | 0 | 12 |
| A2ML1 | chr12:8975 | 36 | 5109 | Gencode, f | 0 | 0 | 36 |
| A3GALT2 | chr1:33777 | 5 | 1123 | RefSeq | 0 | 0 | 5 |
| A4GALT | chr22:4308 | 1 | 1357 | Gencode, f | 0 | 0 | 1 |
| A4GNT | chr3:13784 | 2 | 1063 | Gencode, f | 0 | 0 | 2 |
| AACS | chr12:1255 | 18 | 2379 | Gencode, f | 0 | 0 | 18 |
| AADAC | chr3:15153 | 5 | 1312 | Gencode, f | 0 | 0 | 5 |
| AADACL2 | chr3:15145 | 5 | 1306 | Gencode, f | 0 | 0 | 5 |
| AADACL3 | chr1:12776 | 4 | 1137 | Gencode, f | 0 | 0 | 4 |
| AADACL4 | chr1:12704 | 4 | 1304 | Gencode, f | 0 | 0 | 4 |
| AADAT | chr4:17098 | 14 | 1610 | Gencode, f | 0 | 0 | 14 |
| AAED1 | chr9:99404 | 6 | 801 | RefSeq | 0 | 0 | 6 |
| AAK1 | chr2:69688 | 22 | 3490 | Gencode, f | 0 | 0 | 22 |
| AAMDC | chr11:7755 | 3 | 429 | RefSeq | 0 | 0 | 3 |
| AAMP | chr2:21917 | 11 | 1598 | Gencode, f | 0 | 0 | 11 |
| AANAT | chr17:7444 | 4 | 839 | Gencode, f | 0 | 0 | 4 |
| AAR2 | chr20:3487 | 3 | 1215 | RefSeq | 0 | 0 | 3 |
| AARD | chr8:11795 | 2 | 508 | RefSeq | 0 | 0 | 2 |
| AARSD1 | chr17:4116 | 20 | 2448 | Gencode, f | 0 | 0 | 20 |
| AASDH | chr4:57204 | 14 | 3615 | Gencode, f | 0 | 0 | 14 |
| AASDHPPT | chr11:1055 | 6 | 1050 | Gencode, f | 0 | 0 | 6 |
| AATF | chr17:3536 | 12 | 1923 | Gencode, f | 0 | 0 | 12 |
| AATK | chr17:7905 | 15 | 4472 | Gencode, f | 0 | 0 | 15 |
| ABCA10 | chr17:6714 | 40 | 5640 | Gencode, f | 0 | 0 | 40 |
| ABCA11P | chr4:42955 | 6 | 3167 | Gencode | 0 | 0 | 6 |
| ABCA13 | chr7:48211 | 63 | 16531 | Gencode, f | 0 | 0 | 63 |

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SureDesign Help - Add/Review Content

SureSelect DNA (Advanced) **Add/Review Content**

Define Design

Add/Review Content

- Select Optimized Probes
- Tile Genes or Regions
- Upload Probes
- Select Existing Probes
- Finalize

SureSelect DNA Design

Name: Oncogene
Species: H. sapiens (hg19)
Category: SureSelect DNA
Hybridization: 90 Minutes

Probes

Probes: NA
Size: NA
Price Tier: NA

UCSC View Download

Select Content Addition Method

- ☒ Select Optimized Probes
- ☐ Design new probes by tiling genes or regions
- ☐ Upload probes
- ☐ Select probes from an existing Design or Probegroup

Select all probes

- Or -

- ☐ Finalize this design

Probegroup Summary

| Name | # Targets | Target Size | # Probes | Probe Size | Replicate Count | Actions |
|------------|------------|-------------|----------|------------|-----------------|---------|
| Oncogene_1 | Processing | | | | | |

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
Price Tiers:

- Tier 1: 1 - 499 kbp (up to 60,500 probes)
- Tier 2: 0.5 - 2.999 Mbp (up to 121,000 probes)
- Tier 3: 3.0 - 5.999 Mbp (up to 181,500 probes)
- Tier 4: 6.0 - 11.999 Mbp (up to 242,000 probes)
- Tier 5: 12.0 - 24 Mbp (up to 363,000 probes)

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<ProbeGroup名>_Report.txt

```
# File Summary
File Type: SureSelect Design Report
Created By: SureSelect Advanced Design Wizard
```

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<ProbeGroup名>_Targets.txt


```
# File Type: Input Targets File
# Timestamp: 12-Mar-2013
# User: test@agilent.com
# Created By: SureSelect DNA Pr
```

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The screenshot displays the SureDesign software interface during the 'Add/Review Content' step. The sidebar on the left shows the 'Define Design' section with a green checkmark and the 'Add/Review Content' section with a list of options: 'Select Optimized Probes', 'Tile Genes or Regions', 'Upload Probes', 'Select Existing Probes', and 'Finalize'. Below this, the 'SureSelect DNA Design' section shows details for an 'Oncogene' design for 'H. sapiens (hg19)', including 1710 probes and a size of 1,004 kbp. The main panel displays the 'UCSC Genome Browser on Human (GRCh37/hg19)' view of a genomic region (chr12:57,866,038-57,873,633, 7,596 bp). The browser view includes a scale bar, a track of 'Agilent SureSelect DNA - Oncogene - Genomic regions overlapped by probes', and various genomic annotations including 'GLI1', 'ARHGAP9', and 'RefSeq genes from NCBI'. An orange arrow points from the 'UCSC View' button in the sidebar to the UCSC browser window. The bottom bar contains buttons for 'Close Design Wizard', 'Back', and 'Next'.

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Select Content Addition Method

☒ Select Optimized Probes
☐ Design new probes by tiling genes or regions
☐ Upload probes
☐ Select probes from an existing Design or Probegroup

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

SureDesign Home Find Designs Create Designs

Start Designing

- SureSelect DNA
- SureSelect RNA
- SureSelect All-In-One
- OneSeq
- HaloPlex
- CGH

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Define Design ☒ Add/Review Content ☐ Finalize

Select Content Addition Method

☐ Select Optimized Probes
☐ Design new probes by tiling genes or regions
☐ Upload probes
☐ Select probes from an existing Design or Probegroup

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SureDesign Help - Finalize Design

SureSelect DNA (Advanced) Add/Review Content Finalize

Define Design ☒ Add/Review Content ☒ Finalize

You are now ready to finalize this design. After finalizing:

- The design can no longer be modified.
- A design ID will be assigned.
- Design information will be sent to Agilent's manufacturing facility. Manufacturing will not begin until you have submitted a purchase order.

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

The image shows the SureDesign web interface. On the left, there are input fields for Design Name, Design ID, Quantity (set to 1), Automation Kit (radio buttons for None, Bravo, Magnis), Hybridization (90 Minutes), Sequencing Platform (HiSeq), Reagent Kit (XT HS2, index 1-16), Sample size (16), Beads (Included in the kit), Design Part Number (5191-6900), and Reagent Kit Part Number (G9982A). On the right, there are fields for Workgroup Name, Company Name, and Company Number. Below these fields, a message states: "You are currently logged into SureDesign. SureDesign shares a shopping cart with the Agilent Genomics website. To purchase online, you also need to be registered with that site. To check-out, click the cart icon, and you will be transferred to the Agilent Genomics page. You must be logged in to see your discounted prices." At the bottom, there are three buttons: "Add to Cart", "Request a Quote" (highlighted with a yellow border), and "Cancel".

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SureDesign: Request for quote

We have received your request-for-quote for the following SureSelect DNA kit:

Design Details
Design ID:
Design Name: Oncogene
Species: H. sapiens

Kit Details
Sequencing Platform: HiSeq
Reactions: 16
Quantity: 1
Design Part Number: 5190-4806
Reagent Kit Part Number: G9611A

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
SureDesign Help - Define Design

SureSelect DNA (Advanced) Define OneSeq Design

Define OneSeq Design
Select Spike-in Design
Finalize

* Design Name:

* Species: H. sapiens [Select](#)

Build: UCSC hg19, GRCh37, February 2009

* Create In: Agilent [Select](#)

Description:

Keywords:

Cancel Next

SureSelect DNA Design
Name: NA
Species: NA
Category: OneSeq
Size: NA

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SureDesign Help - Select Catalog and Spike-in Design

SureSelect DNA (Advanced) Select Spike-in Design

Define OneSeq Design ☒ **Select Spike-in Design** Finalize

* Catalog Design: OneSeq 1Mb CNV Backbone
* Workspace: Agilent OneSeq CNV Backbone
* Folder: Agilent Catalog

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SureDesign Help - Finalize Design

SureSelect DNA (Advanced) **Finalize**

Define OneSeq Design ✓
Select Spike-in Design ✓
▶ **Finalize**

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