

## Nitrogen Containing Herbicides (EPA Method 507)

**Column:** DB-35

30 m x 0.53 mm I.D., 0.50 µm

**J&W P/N:** 125-1937

**Carrier:** Helium at 38 cm/sec (5 mL/min),  
measured at 150°C

**Oven:** 60°C for 1 min  
60-290°C at 15°/min  
290°C for 5 min

**Injector:** Megabore Direct, 290°C, 1 µL of 3 ng/µL standard  
**Detector:** NPD, 290°C

1. Eptam®
2. Sutan®
3. Vernam®
4. Tillam®
5. Ordram®
6. Treflan®
7. Balan®
8. Ro-Neet®
9. Propachlor
10. Tolban®
11. Propazine
12. Atrazine
13. Simazine
14. Terbacil
15. Sencor®
16. Dual®
17. Paarlán®
18. Prowl®
19. Bromacil
20. Oxadiazon
21. GOAL®
22. Hexazinone

