

Acrylate Impurities on DB-200

Column: DB-200

30 m x 0.53 mm, 1.0 μ m
J&W P/N: 125-2032
Carrier: Helium, 34.5 cm/sec at 35°C
Oven: 35°C for 5 minutes
10°/min to 200°C
Injector: Split, 1:10, 230°C
Detector: FID, 250°C

1. Methanol
 2. Ethanol
 3. Methyl acetate
 4. Ethyl acetate
 5. Methyl acrylate
 6. Methyl propionate
 7. Isobutanol
 8. Butanol
 9. Ethyl acrylate
 10. Ethyl propionate
 11. Methyl methacrylate
 12. Isopropyl acrylate
 13. Isobutyl acetate
 14. Butyl acetate
 15. Isobutyl acrylate
 16. Dibutyl ether
 17. Isobutyl propionate
 18. Butyl acrylate
 19. Butyl propionate
- * Methyl acrylate impurities
• Ethyl acrylate impurities
♦ Butyl acrylate impurities

