

## Hallucinogens

### Column: DB-5ms

30 m x 0.25 mm I.D., 0.25 µm

### J&W P/N: 122-5532

Carrier: Helium at 30 cm/sec, measured at 50°C

Injector: Splitless, 250°C

30 sec purge activation time

Detector: 1 µL of 10-50 ng/µL standard in methanol

MSD, 300°C transfer line

full scan at m/z 40-350

1. 4-Methylaminorex
  2. Mescaline
  3. Dimethyltryptamine
  4. Ketimine
  5. TCP (Tenocyclidine)
  6. PCP
  7. Diethyltryptamine
  8. Bufotenine (peak not shown)
  9. Fentanyl
  10. Ibogaine
  11. LSD
  12. Psilocin (peak not shown)
- \* Impurity

