




Calibrant to use on Single Quad LC/MS

Below is a matrix that shows the type of calibrant to be used with the applicable SQ LC/MS model and source combinations:

		Ion Source				
		ESI	APCI	APPI	MM	AJS
LC/MS Model	G6110A	ES	APCI/APPI	APCI/APPI	APCI/APPI	N/A
	G6120A/B/C	ES	APCI/APPI	APCI/APPI	APCI/APPI	N/A
	G6125B/C	ES	APCI/APPI	APCI/APPI	APCI/APPI	N/A
	G6130A/B/C	ES	APCI/APPI	APCI/APPI	APCI/APPI	ES-L
	G6135B/C	ES	APCI/APPI	APCI/APPI	APCI/APPI	ES-L
	G6140A	ES-L	APCI-L	APCI/APPI	ES-L	N/A
	G6150B	ES-L	APCI-L	APCI/APPI	ES-L	ES-L

Tune Calibrant	Part Number	
ES	G2421-60001	
APCI/APPI	G2432-60001	
ES-L	G1969-85000	
APCI-L	G1969-85010	