

キャピラリGCカラム-フェラル対応表(ロングフェラル)

nominal 内径 [mm]	nominal 外径 [mm]	素材	工場*	Prefixesブランド名	グラファイトフェラル		ベスペルフェラル(SP-1)		85% ベスペル /15% グラファイト (SP-21)		60% ベスペル /40% グラファイト(SP-22)		
					Nom ID	部品番号	Nom ID	部品番号	Nom ID	部品番号 (空焼済)	部品番号	Nom ID	部品番号
0.025	0.15	Fused Silica	MBG	CP, VF, Select	0.4	8010-0301	0.3	8010-0306	0.3	5062-3507	(空焼なし)	0.3	8004-0211
0.05	0.2	Fused Silica	MBG	CP, VF, Select	0.4	8010-0301	0.3	8010-0306	0.3	5062-3507	NA	0.3	8004-0212
0.075	0.26	Fused Silica	MBG	CP, VF, Select	0.4	8010-0301	0.3	8010-0306	0.3	5062-3507	NA	0.3	8004-0211
0.1	0.37	Fused Silica	FOL	DB, HP	0.5	8010-0302	0.4	8010-0307	0.4	5062-3508	8010-0310	0.425	8004-0212
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.425	8004-0212
0.15	0.36	Fused Silica	FOL	DB, HP	0.5	8010-0302	0.4	8010-0307	0.4	5062-3508	8010-0310	0.425	8004-0212
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.425	8004-0212
0.18	0.36	Fused Silica	FOL	DB, HP	0.5	8010-0302	0.4	8010-0307	0.4	5062-3508	8010-0310	0.425	8004-0212
0.2	0.36	Fused Silica	FOL	DB, HP	0.5	8010-0302	0.4	8010-0307	0.4	5062-3508	8010-0310	0.425	8004-0212
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.425	8004-0212
0.21	0.39	Fused Silica	MBG	CP, VF, Select	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.425	8004-0212
0.25	0.36	Fused Silica	FOL	DB, HP	0.5	8010-0302	0.4	8010-0307	0.4	5062-3508	8010-0310	0.425	8004-0212
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.425	8004-0212
0.32	0.45	Fused Silica	FOL	DB, HP	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.5	8004-0214
	0.45	Fused Silica	MBG	CP, VF, Select	0.5	8010-0302	0.5	8010-0308	0.5	5062-3506	8010-0311	0.5	8004-0214
0.32	0.53	CP-SilicaPLOT FS	MBG	CP, VF, Select	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218
0.45	0.69	Fused Silica	FOL	DB, HP	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218
0.53	0.69	Fused Silica	FOL	DB, HP	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218
	0.69	Fused Silica	MBG	CP, VF, Select	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218
0.25	0.5	Ultimetal	MBG	CP, VF, Select	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218
0.32	0.52	Ultimetal	MBG	CP, VF, Select	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218
0.53	0.75	Ultimetal	MBG	CP, VF, Select	0.8	8010-0304	0.8	8010-0309	0.8	5062-3538	8010-0313	0.8	8004-0218

*MBG: Middleburg(旧パリアンGCカラム工場)、FOL: Folsom(アジレントGCカラム工場)

キャピラリGCカラム-フェラル対応表(ショートフェラル)

nominal 内径 [mm]	nominal 外径 [mm]	素材	工場*	Prefixesブランド名	グラファイトフェラル		ベスペルフェラル(SP-1)		85% ベスペル /15% グラファイト (SP-21)	
					Nom ID	部品番号	Nom ID	部品番号	Nom ID	部品番号
0.025	0.15	Fused Silica	MBG	CP, VF, Select	0.4	500-2114	0.26	5062-3515	0.26	5062-3518
0.05	0.2	Fused Silica	MBG	CP, VF, Select	0.4	500-2114	0.26	5062-3515	0.26	5062-3518
0.075	0.26	Fused Silica	MBG	CP, VF, Select	0.4	500-2114	0.37	5062-3515	0.37	5062-3516
0.1	0.37	Fused Silica	FOL	DB, HP	0.5	5080-8853	0.4	5181-3322	0.4	5181-3323
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
0.15	0.36	Fused Silica	FOL	DB, HP	0.5	5080-8853	0.4	5181-3322	0.4	5181-3323
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
0.18	0.36	Fused Silica	FOL	DB, HP	0.5	5080-8853	0.4	5181-3322	0.4	5181-3323
0.2	0.36	Fused Silica	FOL	DB, HP	0.5	5080-8853	0.4	5181-3322	0.4	5181-3323
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
0.21	0.39	Fused Silica	MBG	CP, VF, Select	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
0.25	0.36	Fused Silica	FOL	DB, HP	0.5	5080-8853	0.4	5181-3322	0.4	5181-3323
	0.39	Fused Silica	MBG	CP, VF, Select	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
0.32	0.45	Fused Silica	FOL	DB, HP	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
	0.45	Fused Silica	MBG	CP, VF, Select	0.5	5080-8853	0.475	5062-3513	0.475	5062-3514
0.32	0.53	CP-SilicaPLOT FS	MBG	CP, VF, Select	0.8	500-2118	0.74	5062-3511	0.74	5062-3512
0.45	0.69	Fused Silica	FOL	DB, HP	0.8	500-2118	0.74	5062-3511	0.74	5062-3512
0.53	0.69	Fused Silica	FOL	DB, HP	0.8	500-2118	0.74	5062-3511	0.74	5062-3512
	0.69	Fused Silica	MBG	CP, VF, Select	0.8	500-2118	0.74	5062-3511	0.74	5062-3512
0.25	0.5	Ultimetal	MBG	CP, VF, Select	0.8	500-2118	0.74	5062-3511	0.74	5062-3512
0.32	0.52	Ultimetal	MBG	CP, VF, Select	0.8	500-2118	0.74	5062-3511	0.74	5062-3512
0.53	0.75	Ultimetal	MBG	CP, VF, Select	0.8	500-2118				

*MBG: Middleburg(旧パリアンGCカラム工場)、FOL: Folsom(アジレントGCカラム工場)

キャピラリGCカラム-PTV注入口(グラフパック2Mフェラル)対応表

nominal 内径 [mm]	nominal 外径 [mm]	素材	工場*	Prefixesブランド名	GRAPHPACK 2M フェラル		GRAPHPACK 2M 注入口アダプタ		
					Nom ID	部品番号	対応外径	部品番号	部品名
0.025	0.15	Fused Silica	MBG	CP, VF, Select			0.31	5182-9754	Graphpack-2M connector 0.2mm
0.05	0.2	Fused Silica	MBG	CP, VF, Select	0.3	5021-7136	0.31	5182-9754	Graphpack-2M connector 0.2mm
0.075	0.26	Fused Silica	MBG	CP, VF, Select	0.3	5021-7136	0.31	5182-9754	Graphpack-2M connector 0.2mm
0.1	0.37	Fused Silica	FOL	DB, HP	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
	0.39	Fused Silica	MBG	CP, VF, Select	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
0.15	0.36	Fused Silica	FOL	DB, HP	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
	0.39	Fused Silica	MBG	CP, VF, Select	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
0.18	0.36	Fused Silica	FOL	DB, HP	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
0.2	0.36	Fused Silica	FOL	DB, HP	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
	0.39	Fused Silica	MBG	CP, VF, Select	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
0.21	0.39	Fused Silica	MBG	CP, VF, Select	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
0.25	0.36	Fused Silica	FOL	DB, HP	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
	0.39	Fused Silica	MBG	CP, VF, Select	0.4	5021-7137	0.45	5182-9761	Graphpack-2M connector 0.32/0.25mm
0.32	0.45	Fused Silica	FOL	DB, HP	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
	0.45	Fused Silica	MBG	CP, VF, Select	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
0.45	0.69	Fused Silica	FOL	DB, HP	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
0.53	0.69	Fused Silica	FOL	DB, HP	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
	0.69	Fused Silica	MBG	CP, VF, Select	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
0.25	0.5	Ultimetal	MBG	CP, VF, Select	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
0.32	0.52	Ultimetal	MBG	CP, VF, Select	0.7	5021-7134	0.7	5182-9762	Graphpack-2M connector 0.53mm
0.53	0.75	Ultimetal	MBG	CP, VF, Select					

*MBG: Middleburg(旧パリアンGCカラム工場)、FOL: Folsom(アジレントGCカラム工場)