

Selected Oxygenates on GS-OxyPLOT

Column: GS-OxyPLOT
10m x 0.53 mm I.D.
Agilent P/N: 115-4912
Carrier: Helium at 41 cm/sec
Oven: 150°C Isothermal
Injector: Split, 1:40, 250°C
Detector: FID, 290°C

1. n-Dodecane
2. Methyl t-butyl ether
3. n-Tridecane
4. Iso-Buteraldehyde
5. n-Tetradecane
6. Methanol
7. Acetone
8. n-Pentadecane

