

以下の製品は2012年8月31日をもって終売とさせていただきます。

型番	品名
200180	BacterioMatchReporterStrainCompCells
200190	BacterioMatchII ScreengReprtrCompCell
200192	BacterioMatchIVldtnRprtrCompCell
200195	BacterioMatchIIElectrocompReprtrCell
201101	PhiX174/HaeIIIMarker,100ug
201104	Lambda/HindIIIMarker,2x500ug
201107	Lambda/HindIII-PhiX174/HaeIIIMrkr,100ug
201109	Lambda/HindIIIMarker,500u
211186	PCR-ScriptAmpCloningKitw/ElectroporationCompetentCells
211188	PCR-ScriptAmpCloningKit,10Rx
211189	PCR-ScriptAmpCloningKit,50Rx
211190	PCR-ScriptAmpCloningKit,25Rx
211192	PCR-ScriptCamCloningKit
211199	pCMV-Script PCR Cloning Kit
240064	pTRGXRTargetPlasmid
240065	BacterioMatchTwo-hybridSystemVector
240066	pTRFTargetPlasmid,Supercoiled
240067	pBTPlasmid,Supercoiled
300141	ProteinaseK,500mg
302010	Humanbeta-ActinRT-PCRPrimer
302047	HumanG6PDRT-PCRPrimer
302055	Humana-actinRT-PCRPrimer
302110	Mousebeta-ActinRT-PCRPrimer
600152	SP6RNAPolymerase,10,000U
720130	CntrlRT-PCRPrmrs,4GenesVaryAbundLvl,Hmn
929900	AgilentCloningPCRComboPack
929901	AgilentSite-MutagenesisComboPack
929902	AgilentGeneralRTPCRComboPack
929903	AgilentGenotypingPCRComboPack
929904	AgilentGeneralPCRComboPack