

以下の製品は2012年12月31日をもって終売とさせていただきます。

型番	品名
540006	Kit,TotalRNA,Brain,Human,3x25ug
540008	TotalRNA,Brain,Cerebellum,Hmn,25ug
540010	Kit,TotalRNA,Colon,Human,3x25ug
540012	Kit,TotalRNA,Heart,Human,3x25ug
540014	Kit,TotalRNA,Kidney,Human,3x25ug
540018	Kit,TotalRNA,Liver,Hmn,3x25ug
540020	Kit,TotalRNA,Lung,Human,3x25ug
540022	Kit,TotalRNA,LymphNode,Hmn,3x25ug
540024	Kit,TotalRNA,Pancreas,Hmn,3x25ug
540026	Kit,TotalRNA,Placenta,Hmn,3x25ug
540030	Kit,TotalRNA,SkeletalMuscle,Hmn,3x25ug
540032	Kit,TotalRNA,Skin,Human,3x25ug
540036	Kit,TotalRNA,Spleen,Human,3x25ug
540038	Kit,TotalRNA,Stomach,Human,3x25ug
540040	Kit,TotalRNA,Thyroid,Human,3x25ug
540042	TotalRNA,Stomach,Diseased,Hmn,3x25ug
540044	Kit,TotalRNA,Uterus,Human,3x25ug
540046	Kit,TotalRNA,Breast,Hmn,Female,3x25ug
540048	Kit,TotalRNA,Larynx,Human,3x25ug
540050	Kit,TotalRNA,Testes,Human,3x25ug
540052	Kit,TotalRNA,Vulva,Human,3x25ug
540054	TotalRNA,BrainStem,Human,3x25ug
540056	Kit,TotalRNA,Omentum,Human,3x25ug
540058	Kit,TotalRNA,FallopianTube,Hmn,3x25ug
540070	Kit,TotalRNA,Rectum,Human,3x25ug
540072	Kit,TotalRNA,Ovary,Human,3x25ug
540074	Kit,TotalRNA,Cervix,Human,3x25ug
540088	Kit,TotalRNA,HeLa-S3Cells,Hmn,3x25ug
540094	Kit,TotalRNA,MDA-MB-453Cells,3x25ug
540104	Kit,TotalRNA,K-562Cells,3x25ug
540106	Kit,TotalRNA,U937Cells,3x25ug
540108	Kit,TotalRNA,JurkatCells,3x25ug
540118	TtlRNA,Brain,FrontalCortex,Hmn,3x25ug
540122	Kit,TotalRNA,CavalVein,Hmn,3x25ug
540124	Kit,TotalRNA,Colon,Ascndng,Hmn,3x25ug
540126	Kit,TotalRNA,Colon,Descndng,Hmn,3x25ug
540128	Kit,TotalRNAParotid,Human,3x25ug
540132	Kit,TotalRNA,Duodenum,Human,3x25ug
540134	Kit,TotalRNA,AdrenalGland,Hmn,3x25ug
540136	TotalRNA,Brain,Striatum,Hmn,3x25ug
540138	TtlRNA,Brain,OccipitalCortx,Hmn,3x25ug

540142	Kit,TotalRNA,Thymus,Human,3x25ug
540144	TtIRNA,BrainParietalCortex,Hmn,3x25ug
540146	Kit,TotalRNA,Trachea,Human,3x25ug
540148	Kit,TotalRNA,Penis,Human,3x25ug
540150	Kit,TotalRNA,Tongue,Human,3x25ug
540158	Kit,TotalRNA,FetalBrain,Hmn,3x25ug
540174	Kit,TotalRNA,FetalLiver,Hmn,3x25ug
540202	MVPUmbilicalCord,HmnTotalRNA,3x25u
600152	SP6RNAPolymerase,10,000U
600816	BrilliantIIQPCRHighROXMasterMix
600832	BrilliantIISYBRQPCRw/highROX
600833	BrilliantIISYBRQPCRlWROXMstrMx10Pk
600839	BrilliantIIQRT-PCRSYBRw/LwROXMstrMx10pk
600840	BrilliantIIQRT-PCRSYBRw/HghROXMstMx10pk
600841	BrilliantIIQRT-PCRLwROXMstrMx1-Stp10pk
600842	BrilliantIIQRT-PCRHghROXMstrMx1-Stp10pk
211186	Kit,PCR-ScriptAmpw/ElecCompcell,25rx
211188	Kit,PCR-ScriptAmpCloning,10rxn
211189	Kit,PCR-ScriptAmpCloning,50rxn
211190	Kit,PCR-ScriptAmpCloning,25rxn
211192	Kit,PCR-ScriptCamCloning,25rxn
501064	SrfIRestrictionEnzyme,300U
600700	let-7amiRNASpcificForwardPrimer
600701	let-7bmiRNASpcificForwardPrimer
600702	let-7cmiRNASpcificForwardPrimer
600703	let-7dmiRNASpcificForwardPrimer
600704	let-7emiRNASpcificForwardPrimer
600705	let-7fmiRNASpcificForwardPrimer
600706	let-7gmiRNASpcificForwardPrimer
600707	let-7imiRNASpcificForwardPrimer
600708	miR-15amiRNASpcificForwardPrimer
600709	miR-15bmiRNASpcificForwardPrimer
600710	miR-16miRNASpcificForwardPrimer
600711	miR-17-3pmiRNASpcificForwardPrimer
600712	miR-17-5pmiRNASpcificForwardPrimer
600713	miR-19amiRNASpcificForwardPrimer
600714	miR-20amiRNASpcificForwardPrimer
600715	miR-21miRNASpcificForwardPrimer
600716	miR-23amiRNASpcificForwardPrimer
600717	miR-23bmiRNASpcificForwardPrimer
600718	miR-25miRNASpcificForwardPrimer
600719	miR-29amiRNASpcificForwardPrimer
600720	miR-29b-2miRNASpcificForwardPrimer

600721	miR-21c miRNAspecificForwardPrimer
600722	miR-30c miRNAspecificForwardPrimer
600723	miR-23 miRNAspecificForwardPrimer
600724	miR-92-1 miRNAspecificForwardPrimer
600725	miR-98 miRNAspecificForwardPrimer
600726	miR-106a miRNAspecificForwardPrimer
600727	miR-107 miRNAspecificForwardPrimer
600728	miR-122a miRNAspecificForwardPrimer
600729	miR-126 miRNAspecificForwardPrimer
600730	miR-128b miRNAspecificForwardPrimer
600731	miR-140 miRNAspecificForwardPrimer
600732	miR-143 miRNAspecificForwardPrimer
600733	miR-145 miRNAspecificForwardPrimer
600734	miR-146a miRNAspecificForwardPrimer
600735	miR-146b miRNAspecificForwardPrimer
600736	miR-155 miRNAspecificForwardPrimer
600737	miR-188 miRNAspecificForwardPrimer
600738	miR-199a-1 miRNAspecificForwardPrimer
600739	miR-200b miRNAspecificForwardPrimer
600740	miR-219 miRNAspecificForwardPrimer
600741	miR-331 miRNAspecificForwardPrimer
600742	miR-125a-5p miRNAspecificForwardPrimer
600743	miR-125b miRNAspecificForwardPrimer
600744	miR-142-3p miRNAspecificForwardPrimer
600745	miR-142-5p miRNAspecificForwardPrimer
600746	miR-181a miRNAspecificForwardPrimer
600747	miR-181b miRNAspecificForwardPrimer
600748	miR-210 miRNAspecificForwardPrimer
600749	miR-361-5p miRNAspecificForwardPrimer