

以下の製品は2023年9月30日をもって終売とさせていただきます。

型番	品名
200049	G418 Sterile Solution,10 ml,50 mg/ml
200160	Able Electro-Comp Cells 5 x 0.1 ml
200162	Able K Electro-Comp Cells 5 x 0.1 ml
200230	XL1-Blue MRF SuperComp Cells 5 x 0.2 ml
200232	AG1 Competent Cells 5 x 0.2 ml
200248	XL1-Blue MRF Kan SuperComp Cell 5x0.2ml
200317	XL10-Gold Kanr Comp Cells 10 x 0.1 ml
200339	RNAMaxx High Yield Transcription Kit
200388	Transfection MBS Mammalian Kit
200410	Klenow Fill-In Kit
200435	10x Taq DNA Polymerase Buffer, 10 ml
200473	RAT Anti-DYKDDDDK Tag Antibody
200474	RAT Anti-DYKDDDDK Tag Antibody
200488	ViraPack Transfection Reagent
200532	10x Cloned Pfu DNA Polymerase Buffr,10ml
201115	Kb DN+M3:M175A Ladder, 250 ug
201221	QuikHyb Rapid Hybridization Sol,1 L
204131	GeneJammer Transfection Reagent,4x1 ml
211173	pCMV-Tag3 Vector
211174	pCMV-Tag4 Vector
211521	pET-3 a-d Vectors
211523	pET-11a-d Vectors
211623	pET-11 a-d Expression System
213201	pMC1Neo and pMC1Neo Poly A Vectors
214301	pCAL-c Vector
214302	pCAL-n Vector
214310	pCAL-n-EK Vector
214311	pCAL-n-FLAG Vector
216201	pSG5 Expression Vector
217450	LacSwitch II Mammalian Expression Sys
217452	pESC-LEU Vector
217461	pEGSH Vector
217564	pFB-ERV Vector
230196	ArcticExpress RP Competent Cells
230265	BL21-CodonPlus (DE3)-RIL-X Comp Cells
234211	Lambda gt11/EcoR I/CIAP Treated Vctr Kit

239212	ZAP Express BamH I/CIAP-Treated Vector
240099	InterPlay TAP Purification Buffer Kit
240105	Streptavidin Resin 2.5 ml
240107	InterPlay Mammalian TAP Purification Kit
247211	Lambda DASH II/BamH I Vector Kit
248211	Lambda FIX II/Xho I Partial Fill-in Vect
251301	SuperCos I Vector Kit
300385	Prime-It II Random Primer Labeling Kit
300392	Prime-It RmT Random Primer Labeling Kit
300600	Alien QRT-PCR Inhibitor Alert
300602	Alien QPCR Reference Standard (FAM) Kit
300604	Alien QPCR Reference Standard ( VIC) Kit
302108	MycoSensor PCR Assay Kit (100 rxn)
302109	MycoSensor PCR Assay Kit (50 rxn)
400720	RNase-It Ribonuclease Cocktail
400761	StrataPrep Plasmid Miniprep Kit, 50 rxn
400763	StrataPrep Plasmid Miniprep Kit,250 rxn
400766	StrataPrep DNA Gel Extraction,50 rxn
400768	StrataPrep DNA Gel Extraction,250 rxn
400771	StrataPrep PCR Purification Kit, 50 rxn
400775	StrataPrep 96 PCR Purification, 2 plates
400814	Absolutely RNA miRNA Kit
400916	SideStep II Cell Lysis Analysis Kit
401193	Thermal Paper for Sony Printer
401342	THERMO-FAST96,PCRDctnPlate,BAYER250/BX
410098	40-Well Working Rack
410188	384 well plates, Pack of 50
540017	Total RNA, Liver, Human, 25 ug
540065	Total,Ileum,Nn Hdgkns Lymphma,Hmn,25ug
540077	Ttl,Colon,MtdSt(Nrml+Adncrcnma),Hmn,25ug
540079	Total RNA,Colon,Adenocarcinoma,Hmn,25 ug
540081	Ttl RNA,Kidney,Renl Cell Crnma,Hmn,25ug
540085	Total RNA,Bladder, Carcinoma,Hmn, 25 ug
540099	Total RNA,K-562 Cells/PMA Treated,25 ug
540103	Total RNA, K-562 Cells, 25 ug
540143	Total RNA,Brain Parietal Cortx,Hmn,25 ug
540169	Total RNA, Fetal Kidney, Human, 25 ug
600011	T4 DNA Ligase, 300U
600036	miRNA 1stStrand Synthss(50rxnat20uL/rxn)

600089	AccuScript Reverse Transcriptase,50 rxns
600090	AccuScript Reverse Transcriptase,200rxns
600123	T7 RNA Polymerase, 5,000U
600124	T7 RNA Polymerase, 25,000U
600135	Pfu DNA Polymerase (native), 100U
600136	Pfu DNA Polymerase (native), 500U
600151	SP6 RNA Polymerase, 3,000U
600163	Exo- Pfu DNA Polymerase, 250U
600180	AccuScript High Fidelity RT-PCR Kit
600182	EASY A 1-TUBE RT-PCR KIT
600184	ACCUSCRIPT PFU ULTRA II RT-PCR
600195	Taq2000 DNA Polymerase, 100U
600196	Taq2000 DNA Polymerase, 500U
600197	Taq2000 DNA Polymerase, 1000U
600201	Single-Stranded DNA BindingProtein,100ug
600210	TaqPlus Precision PCR System, 100U
600212	TaqPlus Precision PCR System, 1000U
600255	Pfu Turbo Alternative Detergent, 100U
600257	PfuTurbo Alternative Detergent, 500U
600260	Herculase Enhanced DNA Polymerase, 100U
600262	Herculase Enhanced DNA Polymerase, 500U
600264	Herculase Enhanced DNA Polymerase,1000U
600282	SureStart Taq DNA Polymerase, 500U
600284	SureStart Taq DNA Polymerase, 1000U
600310	Herculase Hotstart DNA Polymerase 100U
600312	Herculase Hotstart DNA Polymerase,500U
600314	Herculase Hotstart DNA Polymerase 1000U
600320	PfuTurbo Hotstart DNA Polymerase, 100U
600353	Cloned Pfu DNA Polymerase AD, 100U
600390	PfuUltra Hotstart DNA Polymerase, 100U
600396	PfuUltra Hotstart DNA Poly AD, 1000 U
600400	Easy-A High-FidelityPCR CloningEnzy,100U
600402	Easy-A High-FidelityPCR CloningEnzy,500U
600404	Easy-A High-FidelityPCRCloningEnzy,1000U
600420	PicoMaxx High Fidelity PCR System, 100U
600583	miRNA QPCR Master Mix
600630	PfuUltra Hotstart PCR Master Mix,100 rxn
600650	PicoMaxx High Fidelity PCR MstrMx,100rxn
600680	Paq5000TM DNA Polymerase, 500 U

600682	Paq5000TM DNA Polymerase, 1000 U
600835	Brilliant II QRT-PCR SYBRw/LowROX MstrMx
600837	Brilliant II QRT-PCR LwROXMstrMx,1-Step
600870	PAQ5000 2X PCR MASTER MIX 100
600872	PAQ5000 2X PCR MASTER MIX 400
600920	SureDirect Blood PCR Kit, 100 Rxns
736009	Total RNA, Liver, Mouse, BALB/c, Male
736509	Total RNA, Liver, C57BL, Male
737005	Total RNA,Heart,Sprague-Dawley Rat,Male
737013	Total,Skltl Mscle,Sprague DawleyRat,Male
929568	Amplified S7U SYBR Control Plate
930004	Custom electrocomp LC1062F'
930457	Custom cDNA Conversion Module, 16rxn
931331	BL21 Gold Competent Cells, 96-well plate
931708	Custom SSL Full Length QXT Blockers
5185-5989	HSA Standard, 20mg/mL, 1mL
5188-2750	Peptide Cleanup C18 Spin Tubes, pk/50
5188-5218	Mult Aff Rem Column 4.6 x100mm MS-3
5188-5253	Needles, PTFE, Luer-Lock, pk/10
5190-7717	SureGuide Cas9 Nuclease, 100rx
5500-0001	Agilent PCR-RFLP Fish Species ID Kit
G1990-85001	Trypsin Digest Methylated BSA Standard
G7515A	SureVector E. coli N-term Expansion Kit
G7515B	SureVector E. coli C-term Expansion Kit
G7517A	SureVector N-terminal Yeast
G7517B	SureVector C-terminal Yeast
G7518A	SureVector E. Coli Selection Kit
G7518B	SureVector E. coli N-term Promoter Kit
G7518C	SureVector E. coli C-term Promoter Kit
G7518D	SureVector E. coli N-terminal Tag Kit
G7518E	SureVector E. coli C-terminal Tag Kit

お問い合わせ先：

アジレント・テクノロジー株式会社

カスタマコンタクトセンタ

TEL 0120-477-111

<http://agilent.com/chem/genomics:jp>